CAMPYLOBACTER growth Supplement

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Rev.0 / 06.04.2005

CAMPYLOBACTER growth Supplement

Growth supplement for the isolation of *Campylobacter spp*

DESCRIPTION

CAMPYLOBACTER growth Supplement is a growth supplement for the isolation of *Campylobacter spp*, comprising a freeze-dried mixture of Sodium pyruvate, Sodium metabisulphite and ferrous sulphate heptahydrate.

KIT CONTENTS

Each kit contains:

- 10 bottles of CAMPYLOBACTER growth Supplement freeze-dried
- 1 Instruction sheet

PRINCIPLE OF THE METHOD

Micro-organisms in the genus *Campylobacter* are microaerophiles and are inhibited by normal atmospheric levels of oxygen: the optimum level of oxygen is 6%. These micro-organisms are extremely sensitive to the toxic forms of oxygen (superoxide anions, peroxides etc) that may appear in mediums incubated in aerobiosis, but the addition of ferrous sulphate, sodium metabisulphite and sodium pyruvate increases the aerotolerance of *Campylobacter jejuni/coli*. These compounds bind to the toxic oxygen molecules, inactivating them and permitting quick and efficacious isolation of *Campylobacter* species.

COMPOSITION

CAMPYLOBACTER growth Supplement		
	Content / bottle	Content / l of medium
Sodium Pyruvate	125.0 mg	250.0 mg
Sodium Metabisulphite	125.0 mg	250.0 mg
Ferrous Sulphate Heptahydrate	125.0 mg	250.0 mg

PROCEDURE FOR USE

1. Reconstitute the contents of a bottle of CAMPYLOBACTER growth Supplement aseptically with 5 ml of sterile distilled water. Shake until completely dissolved, avoiding foam formation.
2. Add the entire contents of a bottle (5 ml) aseptically to 500 ml of *Campylobacter* Agar Base medium code 611007-621007 autoclaved, cooled to 45-50 °C with a 5-10% addition of lysate of horse or ram blood, when required, and the contents of a bottle of the selective supplement chosen for the isolation of *Campylobacter*.
3. Mix with care.
4. Distribute onto Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the culture mediums chosen for the isolation of *Campylobacter*.

QUALITY CONTROL

1. Control of the appearance: greenish freeze-dried product.
2. Microbiological control.

Prepare the plates using as base CAMPYLOBACTER AGAR BASE code 611007 or 621007 with a 5% addition of ram blood, and enriched with CAMPYLOBACTER growth Supplement (1 bottle in 500 ml of medium).

The plates are seeded with the strain indicated in the microbiological control table.

Incubation conditions: 48 h at 42°C, in a microaerophilic atmosphere for *Campylobacter*.

Microbiological control:

Control strain	Growth
<i>Campylobacter fetus subsp. jejuni</i> ATCC 33291	good

PRECAUTIONS

The product CAMPYLOBACTER growth Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. CAMPYLOBACTER growth Supplement is a growth supplement to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store CAMPYLOBACTER growth Supplement at 2-8°C in its original packaging. In such conditions CAMPYLOBACTER growth Supplement maintains its validity until the expiry date indicated on the label. Do not use after that date. Eliminate without using if there are signs of deterioration.

REFERENCES

- Kiggins, E.M, and Plastringe, W.N. 1956. J. Bacteriol. **72**:397-400.
- George, H.A., Hoffman, P.S., and Krieg, N.R. 1978. J. Clin. Microb. **8**:36-41.
- Association of Official Analytical Chemists. 1995. Bacteriological analytical manual, 8th Ed.

PRESENTATION

product	REF	Σ
CAMPYLOBACTER growth Supplement	81050	10 bottles

One bottle is used to prepare 500 ml of medium

TABLE OF SYMBOLS

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code

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