VANCOMYCIN Supplement

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Rev.1 / 29.08.2012

VANCOMYCIN Supplement

Selective supplement for the isolation of Vancomycin-resistant Enterococci (VRE) from clinical samples

DESCRIPTION

VANCOMYCIN Supplement is a freeze-dried selective supplement constituted by Vancomycin, to use as a supplement of the culture medium VRE AGAR BASE, ref. 610501 for Vancomycin-Resistant Enterococci identification.

KIT CONTENTS

Each kit contains:

- 10 vials of VANCOMYCIN Supplement freeze-dried
- 1 instructions sheet

PRINCIPLE OF THE METHOD

Vancomycin is a bactericide antibiotic which inhibits the synthesis of the cell-walls. It is usually active against Gram-positive cocci and bacilli, included (with rare exceptions) *Staphylococcus aureus* and coagulase-negative staphylococci, which are resistant to cephalosporin and penicillin. Vancomycin is bacteriostatic against enterococci even if many strains of *E. faecium* are resistant. All Gram-negative microorganisms are resistant to vancomycin.

COMPOSITION

<u>VANCOMYCIN Supplement</u>		
	Content/vial	Content/liter of medium
Vancomycin	3.0 mg	6.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of VANCOMYCIN Supplement with 2 mL of sterile distilled water.
2. Mix to complete dissolution and add to 500 ml of VRE AGAR BASE, ref. 610501, autoclaved and cooled to 45-50°C.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for VRE AGAR BASE, ref. 610501.

QUALITY CONTROL

1. Control of the appearance: freeze-dried white in colour.
2. Microbiological control.

Prepare the plates using as base VRE AGAR BASE, ref. 610501 added with the content of 1 vial of VANCOMYCIN Supplement (1 vial in 500 ml of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1 °C for 18-48 hours in atmosphere. Microbiological control:

	Control strains	Growth
<i>Escherichia coli</i>	ATCC 25922	Inhibited
<i>Enterococcus faecalis</i>	ATCC 29212	Inhibited
<i>Enterococcus faecalis</i>	ATCC 51299	Good

PRECAUTIONS

The product VANCOMYCIN Supplement does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. VANCOMYCIN Supplement is intended for in vitro diagnostic use only and must be used by properly trained operators with approved methods of aseptic and safety against pathogens.

STORAGE AND TRANSPORT CONDITIONS

Store VANCOMYCIN Supplement at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- King, W. K. (1996) Bug Bytes Vol. 2 No. 19.
- CDC Preventing the spread of vancomycin resistance: a report from the Hospital Infection Control Practices Advisory Committee (1994). Fed Regist. May 17.
- Gold, H. S. & Moellering R.C. Jr. (1996) N. Engl. J. Med.; **335**(19): 1445-53.
- Weinbren, M. J., Johnson, A.P. & Woodford, N. (2000) J. Antimicrobial Chemotherapy; **45**:404-405.

PRESENTATION

Product	REF	
VANCOMYCIN Supplement	81062	10 vials

One vial is sufficient to prepare 500 ml of medium

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> diagnostic medical device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse