m- CP Supplement A

Phenolphthalein diphosphate

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m- CP Supplement A

Phenolphthalein diphosphate

Supplement for the preparation of the medium m- CP AGAR BASE
for the isolation of *Clostridium perfringens*.

ENGLISH

DESCRIPTION

m- CP Supplement A is a supplement constituted by phenolphthalein phosphate, to use as a supplement of the culture medium m- CP AGAR BASE ref. 610199 or 620199 for *Clostridium perfringens* isolation.

KIT CONTENTS

Each kit contains:

- 10 vials of **m- CP Supplement A**
- 1 instructions sheet

PRINCIPLE OF THE METHOD

m- CP Supplement A allows the presumptive identification of *Clostridium perfringens* colonies which can be tested for acid phosphatase activity by the exposure to ammonium hydroxide vapour for 20-30 seconds.

COMPOSITION

| m- CP Supplement A | | |
|---------------------------|----------------|-----------------------|
| | Content / vial | Content / L of medium |
| Phenolphthalein phosphate | 50.0 mg | 100.0 mg |

PROCEDURE FOR USE

1. Aseptically dissolve the content of one vial of **m- CP Supplement A** with 10 mL of sterile distilled water.
2. Mix to complete dissolution and add to 500 mL of m- CP AGAR BASE ref. 610199 or 620199, autoclaved, cooled to 45-50°C and added with the content of one vial of m- CP Supplement (ref. 81034), one vial of m- CP Supplement B (ref. 81071) and vial of m- CP Supplement C (ref. 81072).
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for m- CP AGAR BASE ref. 610199 or 620199.

QUALITY CONTROL

1. Microbiological control.

Prepare the plates using as base m- CP AGAR BASE ref. 610199 or 620199 added with the content of 1 vial of **m- CP Supplement A** (1 vial in 500 mL of medium) and with the content of one vial of m- CP Supplement (ref. 81034), one vial of m- CP Supplement B (ref. 81071) and vial of m- CP Supplement C (ref. 81072). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 44 ± 1°C for 21 ± 3 hours in anaerobic conditions. Microbiological control:

| Control strains | Growth/ Color of colonies | Color of colonies of NH ₄ OH |
|---|---------------------------|---|
| <i>Clostridium perfringens</i> ATCC 13124 | Good/ Yellow | Pink/ Red |
| <i>Clostridium sporogenes</i> ATCC 19404 | Good/ Purple | |
| <i>Escherichia coli</i> ATCC 25922 | Inhibited | |

PRECAUTIONS

m- CP Supplement A does not contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. **m- CP Supplement A** is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store **m- CP Supplement A** at 2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

REFERENCES

- Bisson, J. W. and Cabelli, V. J. (1979) *Applied and Environmental Microbiology*, Vol. 37, No. 1, pp 55-88.
- E.U. (1998) 98/83/EC of Council of 3rd of November 1998 on the quality of water intended for human consumption. *Off. J. Eur. Commun.*, L330, 32-54.

PRESENTATION

| Product | REF | |
|--------------------|-------|----------|
| m- CP Supplement A | 81070 | 10 vials |

One vial is sufficient to prepare 500 mL of medium

TABLE OF SYMBOLS

| | | | | |
|------------------|---------------------------|--------------|---|------------------------|
| Batch code | Do not reuse | Manufacturer | Contains sufficient for <n> tests | Temperature limitation |
| Catalogue number | Fragile, handle with care | Use by | Caution, consult accompanying documents | |

LIOFILCHEM s.r.l.

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

Rev.0/ 18.11.2008