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MEROPENEM Supplement

ENGLISH

Selective supplement for the isolation of Vancomycin-resistant Enterococci (VRE) from clinical samples

DESCRIPTION

MEROPENEM Supplement is a lyophilized supplement to add to VRE Broth Base and VRE Agar Base (ref. 610501) for the isolation of Vancomycin-resistant Enterococci from clinical samples.

PACKAGING

Each package contains:

MEROPENEM Supplement	Instructions Sheet
----------------------	--------------------

METHOD PRINCIPLE

The proliferation of Enterococci resistant to many commonly used antimicrobials is on the increase. Meropenem is a β -lactam ultra-broad spectrum antibiotic to add to VRE Broth Base and VRE Agar Base for the isolation of Vancomycin Resistant Enterococci (VRE) from clinical samples. Meropenem suppresses contaminating flora, particularly Gram-negatives and *Enterococcus gallinarum*.

COMPOSITION	Content/vial	Concentration in medium
Meropenem	1.0 mg	1.0 mg/l in VRE Agar Base 2.0 mg/l in VRE Broth Base

INSTRUCTIONS FOR USE

Resistant Enterococci can be isolated either directly by inoculation of the clinical sample onto supplemented VRE Agar, or indirectly isolated with a selective enrichment through VRE Broth followed by inoculation onto supplemented VRE Agar. Meropenem is used at 2 mg/l in VRE Broth Base and 1 mg/l in VRE Agar Base (ref. 610501).

TECHNIQUE AND RESULTS INTERPRETATION

VRE Agar Base contains an indicator system to detect the growth of aesculin-hydrolysing organisms. Enterococci produce black zones around the colonies from the formation of black iron phenolic compounds derived from aesculin-hydrolysis products and ferrous iron.

QUALITY CONTROL

- Appearance control: white lyophilized.
- Microbiological control:

Control strain		Growth
<i>Enterococcus faecalis</i>	NCTC 12201	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited
<i>Enterococcus faecalis</i>	ATCC 33186	Inhibited

PRECAUTIONS

The product MEROPENEM Supplement does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. MEROPENEM Supplement is intended for in vitro diagnostic use only and must be used by properly trained operators with approved methods of aseptic and safety against pathogens.

STORAGE AND TRANSPORT CONDITIONS

Store MEROPENEM Supplement at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

PRESENTATION

Product	REF	
MEROPENEM Supplement	81083	10 vials

One vial sufficient for 1 litre of VRE Agar Base (ref. 610501) or 500 ml of VRE Broth Base.

REFERENCES

- King, W.K. (1996) Bug Bytes Vol. 2 No. 19.
- CDC Preventing the spread of vancomycin in resistance: a report from the Hospital Infection Control Practice Advisory Committee (1994) Fed. Regist. May 17.
- Gold, H.S. And Moellering R.C. Jr. (1996) N. Engl. J. Med.; 335(19):1445-53.
- Weinbren, M.J., Johanson, A.P. And Woodford, N. (2000) J. Antimicrobial Chemotherapy; 45:404-405.

TABLE OF SYMBOLS

LOT Batch code	IVD In vitro diagnostic medical device	Manufacturer	Use by	Fragile, handle with care
REF Catalogue number	Temperature limitation	Contains Sufficient for <n> tests	Caution, consult accompanying documents	Do not reuse

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LIOFILCHEM® S.r.l.