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Selective supplement for the enrichment of *Escherichia coli* O157 from food and faecal samples.

DESCRIPTION

VCC MOD SELECTIVE Supplement is a selective supplement for the detection of *Escherichia coli* O157 from food and faecal samples, used for enriching media such as MODIFIED TRYPTONE SOYA BROTH (mTSB) ref. 610352 or 620352, or BUFFERED PEPTONE WATER ref. 611014 or 621014, or E.E. BROTH ref. 610017 or 620017.

PACKAGE CONTENT

Each package contains:

- 10 vials of VCC MOD SELECTIVE Supplement freeze-dried
- 1 instructions sheet

PRINCIPLE OF METHOD

Vancomycin inhibits synthesis of cell wall of Gram positive microorganisms. Cefixime has a good activity against Enterobacteria, some streptococci, included Group A Streptococci and pneumococci, and a poor activity against staphylococci. Cefsulodin is a wide broad antibiotic.

TYPICAL FORMULA

VCC MOD SELECTIVE Supplement		
	Content (mg)/ vial	Content (mg) / l of medium
Vancomycin	1.8	8.0
Cefixime	0.0028	0.0125
Cefsulodin	2.25	10.0

PROCEDURE

1. Aseptically reconstitute the content of one vial of VCC MOD SELECTIVE Supplement with 2.5 ml of sterile distilled water. Mix until completely dissolved, avoiding foam formation.
2. Aseptically add the content of one vial (2.5 ml) to 225 ml of MODIFIED TRYPTONE SOYA BROTH (mTSB) ref. 610352 or 620352, or BUFFERED PEPTONE WATER ref. 611014 or 621014, boiled, autoclaved at 121°C for 15 minutes and cooled to 45-50 °C. Or aseptically add the content of one vial (2.5 ml) to 225 ml of E.E. BROTH ref. 610017 or 620017, boiled and cooled to 45-50 °C.
3. Mix well.
4. Distribute into sterile tubes.

TECHNIQUE AND INTERPRETATION OF RESULTS

Refer to the technical sheet of MODIFIED TRYPTONE SOYA BROTH (mTSB) ref. 610352 or 620352, or BUFFERED PEPTONE WATER ref. 611014 or 621014, or E.E. BROTH ref. 610017 or 620017.

QUALITY CONTROL

1. Appearance control: white freeze-dried.
2. Microbiological control.

Prepare the tubes using as base the medium MODIFIED TRYPTONE SOYA BROTH (mTSB) ref. 610352 or 620352 enriched with VCC MOD SELECTIVE Supplement. The tubes are inoculated with the strains indicated in the table of microbiological control.

Incubation conditions: 24 ± 3 h at 37 ± 1 °C.

Microbiological control:

Control strains		Growth
<i>Escherichia coli</i> O157	NCTC 12900	Good
<i>Proteus mirabilis</i>	ATCC® 25933	Inhibited
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	Inhibited
<i>Staphylococcus aureus</i>	ATCC® 25923	Inhibited

PRECAUTIONS

The product VCC MOD SELECTIVE Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use.

VCC MOD SELECTIVE Supplement is a selective supplement intended for *in vitro* diagnostic use only. It must be used in a professional environment by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store VCC MOD SELECTIVE Supplement at 2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions VCC MOD SELECTIVE Supplement maintains its validity until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

REFERENCES

- Doyle, M. P. and J. L. Schoeni (1987) *Applied Environmental Microbiology* **53**: 2394-2396.
- Szabo R. E. et al. (1990) *Applied Environmental Microbiology* **56**: 3546-3549.
- Desmarchelier, P. M. and Grau, F. H. (1997) *Escherichia coli*. In: *Foodborne Microorganisms of Public Health Significance*. 5th Edition. pp.231-264. A. D. Hocking (Ed.). AIFST (NSW Branch) Food Microbiology Group, Australia.

PRESENTATION

Product	REF	Σ
VCC MOD SELECTIVE Supplement	81086	10 vials

One vial is sufficient to prepare 225 ml of medium.

TABLE OF SYMBOLS

 LOT	Batch code		Do not reuse
 REF	Catalogue number		Fragile, handle with care
	Manufacturer		Content sufficient for <n> tests
	Use by		Warning, see instructions for use
	Temperature limits		In Vitro Diagnostic Medical Device