Legionella BCYE Growth Supplement w/o L-Cysteine

Selective supplement omitting L-cysteine for presumptive identification of *Legionella* spp.

DESCRIPTION

Legionella BCYE Growth Supplement w/o L-Cysteine is a lyophilized supplement used for the preparation of Legionella CYE Agar Base (ref. 610125, 620125). The complete medium is used for detection and enumeration of *Legionella* spp from clinical specimens and environmental samples, according to ISO 11731 (all parts).

KIT CONTENTS

Each kit contains:

- 10 vials of freeze-dried Legionella BCYE Growth Supplement w/o L-Cysteine;
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

Legionella BCYE Growth Supplement w/o L-Cysteine is used for confirmation of presumed *Legionella* colonies by a comparative examination demonstrating the organism requirements for L-cysteine.

COMPOSITION

	Content / vial	Content / liter of medium
Alpha-Ketoglutarate	0.1 g	1.0 g
ACES Buffer (N-2-acetamido-2-aminoethanesulphonic acid)	1.0 g	10.0 g
Potassium Hydroxide	0.28 g	2.8 g
Iron Pyrophosphate	0.025 g	0.25 g

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of BCYE Growth Supplement w/o L-Cysteine with 10 ml sterile distilled water.
2. Mix to complete dissolution and add to 90 ml Legionella CYE Agar Base, autoclaved and cooled to 45-50°C.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the medium being prepared.

QUALITY CONTROL

1. Control of the appearance: freeze-dried pink to brown in colour.
2. Microbiological control.
Prepare the medium per label directions. Inoculate the plates with the strains indicated below and incubate at 36 ± 2°C for 2-5 days in humidified atmosphere.

Control strains

Legionella pneumophila WDCM 00107
Escherichia coli ATCC® 25922

Growth

Inhibited
Good

PRECAUTIONS

Legionella BCYE Growth Supplement w/o L-Cysteine is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- ISO 11731-2:2004. Water quality – Detection and enumeration of Legionella. Part 2: Direct membrane filtration method for waters with low bacterial count.
- ISO 11731:1998. Water quality – Detection and enumeration of Legionella.
- Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3):298.

PRESENTATION

Product	Ref.	Contents
Legionella BCYE Growth Supplement w/o L-Cysteine	81091	10 vials

One vial is sufficient to prepare 100 ml of medium.

TABLE OF SYMBOLS

 Batch code	 In Vitro Diagnostic Medical Device	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

Rev.1 / 12.01.2017

Legionella BCYE Growth Supplement w/o L-Cysteine

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