D-Cycloserine 4-MUP Supplement

ENGLISH

Selective supplement for detection of *Clostridium perfringens* from clinical specimens and other materials.

DESCRIPTION

D-Cycloserine 4-MUP Supplement is a selective lyophilized supplement used for the preparation of Clostridium perfringens Agar Base (ref. 610207, 620207). The complete medium is used for the identification and enumeration of *Clostridium perfringens* from clinical specimens, food, water and environmental samples.

KIT CONTENTS

Each kit contains:

- 10 vials of lyophilized D-Cycloserine 4-MUP Supplement
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

D-Cycloserine inhibits the accompanying bacterial flora and causes the colonies which develop to remain smaller. It also reduces a diffuse and thus disturbing blackening around the *Clostridium perfringens* colonies. 4-Methylumbelliferylphosphate (MUP) is a fluorogenic substrate for the acid phosphatase which is a high specific indicator for *C. perfringens*.

COMPOSITION

	Content / vial	Content / liter of medium
D-Cycloserine	200 mg	0.4 g
4-Methyl-umbelliferyl-phosphate Disodium Salt	50 mg	0.1 g

PROCEDURE FOR USE

1. Reconstitute aseptically the content of one vial of D-Cycloserine 4-MUP Supplement with 5 ml of sterile distilled water.
2. Mix to complete dissolution and add aseptically to 500 ml of Clostridium perfringens Agar Base autoclaved and cooled to 45-50°C (*).
3. Mix with care and pour into Petri dishes.

(*) If desired, 50 ml Egg Yolk Emulsion (ref. 80219) can be added to the medium to detect lecithinase activity.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical sheet of the medium being prepared.

QUALITY CONTROL

1. Visual inspection: whitish button, limpid colourless solution once reconstituted.
2. Microbiological control.

Prepare the medium per label directions. Inoculate the plates with the microbial strains indicated below and incubate at 44°C for 18-24 h under anaerobic conditions.

Control strains

Growth

Colony Appearance

<i>Clostridium perfringens</i>	WDCM 00007	Good	Black colonies; light blue fluorescence under UV lamp (with a halo if egg yolk emulsion has been added)
<i>Escherichia coli</i>	WDCM 00012	Inhibited	---

PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* Diagnostic Use only and must be used by properly trained operators.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

REFERENCES

- EN ISO 7937: 2004. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of *Clostridium perfringens* – Colony count technique.
- Hauschild, A.H.W., Hilsheimer, R., and Griffith, D.W. (1974) Enumeration of faecal *Clostridium perfringens* spores in egg yolk-free Tryptose – Sulfit – Cycloserine Agar. Appl. Microb., 27:527-530.

PRESENTATION

Product	Ref.	Content
D-Cycloserine 4-MUP Supplement	81098	10 vials

One vial is sufficient to prepare 500 ml of medium.

TABLE OF SYMBOLS

LOT	Batch code	IVD	<i>In vitro Diagnostic Medical Device</i>		Manufacturer		Contains sufficient for <n> tests		Temperature limitation
REF	Catalogue number		Fragile, handle with care		Use by		Caution, consult accompanying documents		Do not reuse

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LIOFILCHEM® S.r.l.