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E-Mail: liofilchem@liofilchem.net

Rev.0 / 07.06.2006

Biological indicators of irradiation sterilization processes containing *Bacillus pumilus* (ATCC 27142) spores in suspension

DESCRIPTION

USP (United States Pharmacopoeia), EP (European Pharmacopoeia) and DAB (Deutsches Arzneibuch) standards recommend to use bioindicators during irradiation sterilization processes.

Biological indicators **SUSPENSION SPORE PUM E8/E7/E6** are used for validation, re-validation and process monitoring of irradiation sterilizers.

Biological indicators **SUSPENSION SPORE PUM E8/E7/E6** are produced under strictly controlled conditions in order to satisfy the requirements indicated in the USP current edition and in accordance with ISO 11138 and EN 866 standards.

PRINCIPLE

SUSPENSION SPORE PUM E8/E7/E6 contain *Bacillus pumilus* (ATCC 7953) spores in distilled water in predefined concentrations: E8=1-5x10⁸CFU/ml; E7=1-5x10⁷CFU/ml; E6=1-5x10⁶CFU/ml.

Spores are completely killed off after the irradiation sterilization, if the process has been efficient and therefore during the following incubation in a suitable medium they are not able to grow and to modify the medium's aspect. On the contrary, if the sterilization process has not been efficient, spores partially survive and therefore during the incubation are able to grow and to modify the medium's aspect.

INSTRUCTIONS FOR USE

- Take a suitable volume of the suspension and put it on the selected material. Protect appropriately the inoculated material in order to avoid spore dispersion.
- Put the inoculated and protected on the critical points of material to sterilize.
- After the sterilization cycle transfer aseptically the material in a container with a culture medium suitable for microorganism growth as for example Steri-Test Medium (code 20199).
- Incubate the containers with the medium and the material at 30-35 °C (86-95 °F) for 7 days or for a shorter time validated by user.
- Incubate, at the same conditions of time and temperature, an unit of material not submitted to the sterilization cycle, inoculated with a spore suspension belonging to the same batch, as spore growth control (positive control).
- Examine medium's aspect in the container and interpret results as per **EVALUATION TABLE**: a change of medium's aspect indicates microbial growth and therefore an unsuccessful sterilization. On the contrary, the persistence of the medium's initial aspect indicates absence of microbial growth and therefore a successful sterilization.

EVALUATION

Bacillus pumilus spores are killed off if the sterilization cycle has been efficient: in this case the medium remains unchanged even after incubation at 30-35 °C (86-95 °F) for the selected time.

If the sterilization cycle has not been efficient, spores partially survive and the medium's aspect change after incubation at 30-35 °C (86-95 °F) for the selected time.

The medium in the container inoculated with the material not submitted to the sterilization cycle and used as spore growth control, has to change aspect after incubation. On the contrary, the test must be repeated after having investigated causes of the negative result.

EVALUATION TABLE

MEDIUM'S ASPECT	SPORES	STERILIZATION
Unchanged	Killed off	Successful
Changed	Vital	Unsuccessful

MATERIAL TREATMENT AFTER USE

After use, sterilize the positive containers (changed medium's aspect) in autoclave at 121 °C for at least 30 minutes and eliminate them in accordance with the procedures of the laboratory.

STORAGE

Store the product at 2-8 °C: in these conditions it maintains its validity until the expiry date indicated on the label.

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TABLE OF SYMBOLS

 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Caution, consult accompanying documents
 Use by	 Batch code	 Do not reuse

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745

Fax +390858930330

Website: www.liofilchem.net

E-Mail: liofilchem@liofilchem.net

Rev.0 / 07.06.2006