Citrate Test

Rapid test for the identification of Enterobacteriaceae and non-fermenting gram-negative bacteria.

INTENDED PURPOSE

Test for presumptive identification of Enterobacteriaceae and non-fermenting gram-negative bacteria. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Citrate Test is used in the identification of Enterobacteriaceae and non-fermenting gram-negative bacteria on the basis of citrate utilization.

KIT CONTENT

- 30 Tubes of Citrate Test

METHOD PRINCIPLE

Organisms able to utilize ammonium dihydrogen phosphate and sodium citrate as the sole sources of nitrogen and carbon respectively will grow on this medium and produce an alkaline reaction as evidenced by a change in the color of the bromthymol blue indicator from green (neutral) to blue (alkaline).

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, culture media, physiological solution (0.85% saline), 0.5 McFarland turbidity standard, quality control organisms.

REAGENTS

- A tube of Citrate Test contains ammonium salt, citrate, bromthymol blue.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Prepare a bacterial suspension (at least McFarland 0.5) of the strain to be tested in 0.25 ml physiological solution.
2. Add the entire microbial suspension to a Citrate Test tube.
3. Incubate at $35\pm 2^{\circ}\text{C}$ for 18-24 hours (positive reactions can sometimes be observed after 4-6 hours incubation).
4. Watch for the color changing of the suspension.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

Citrate utilization (positive reaction) is indicated by the suspension turning blue.

No citrate utilization (negative reaction) is indicated by the green color of the suspension.

STORAGE

Store at $2-8^{\circ}\text{C}$ away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance: an orange lyophilized button per tube.

Control strain		Incubation	Characteristic reactions
<i>Klebsiella pneumoniae</i>	ATCC® 13883	18-24 h/ 35 ± 2°C	The medium turns blue (citrate utilization)
<i>Proteus mirabilis</i>	ATCC® 25933		
<i>Salmonella Typhimurium</i>	ATCC® 14028		
<i>Enterobacter cloacae</i>	ATCC® 23355		The medium remains green (no citrate utilization)
<i>Escherichia coli</i>	ATCC® 25922		

PERFORMANCE CHARACTERISTICS

Performance testing of Citrate Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Luxuriant growth without an accompanying color change may indicate a positive test. However, if the agar does not turn blue on further incubation, the test should be repeated with less inoculum.

Tests with equivocal results should be repeated.

To avoid false-positive reactions, use a light inoculum to prevent the carryover of substances from the previous medium.

The reaction of this medium alone is not sufficient for identification at the species level.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Citrate Test	30 tests	88042

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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