Optochine Test

ITALIANO

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Optochine Test

ENGLISH

Diagnostic discs for pneumococci identification.

DESCRIPTION

Optochine Test is constituted by paper discs, each one containing 5 µg of Optochin (Ethylhydrocupreine hydrochloride), used for differentiating *Streptococcus pneumoniae* from the other alpha-haemolytic streptococci.

CONTENT OF THE PACKAGES

Each package contains:

- 2 cartridges with 50 discs each, packaged in a heat-sealed container.
- Dryer.

PRINCIPLE OF THE METHOD

Optochin is an agent specifically active against *Streptococcus pneumoniae*, the other alpha-haemolytic streptococci are resistant. The disc is placed onto the surface of a culture medium that is suitable for the growth of streptococci, inoculated with a pure liquid culture of the microorganism under examination. After the incubation all the plates are examined for the presence or absence of an inhibition halo around the disc of Optochin.

COMPOSITION

Each disc contains 5 µg of Optochin.

PREPARATION OF THE SPECIMEN

1. Mixed cultures or clinical specimens must not be used to determine susceptibility to Optochin.
2. Inoculate a tube of Brain Heart Infusion Broth (ref. 20104) with pure colonies of the microorganism under examination.
3. Incubate at $36 \pm 1^\circ\text{C}$ overnight.

TEST PROCEDURE

1. Take the cartridges container from the refrigerator and leave it on the test bench until it reaches room temperature (about 30 minutes). This will prevent humidity being deposited on the discs when the package is opened, which could prejudice their long-term stability.
2. Using a sterile swab, evenly inoculate the surface of a plate of blood agar such as Tryptic Soy Blood Agar (ref. 11037), Columbia Blood Agar (ref. 11025) or other blood medium, with the suspension of the streptococcus under examination.
3. Using sterile tools, gently press one disc of Optochin on the inoculated surface.
4. Turn the plate upside down and incubate at $36 \pm 1^\circ\text{C}$ for 18-24 hours in atmosphere containing 5% of CO_2 .
5. Check for presence or absence of an inhibition halo around the disc of Optochin.

INTERPRETATION OF THE RESULTS

The test organism is considered sensitive to Optochin and presumptively *Streptococcus pneumoniae* if the inhibition zone is ≥ 14 mm diameter. The presumptive identification must be confirmed by serological tests.

QUALITY CONTROL

Each batch of Optochine Test is tested for susceptibility to Optochin by using *Streptococcus pneumoniae* ATCC® 6305 for positive control, and *Streptococcus pyogenes* ATCC® 19615 for negative control, inoculated on Columbia Blood Agar with 5% of defibrinated sheep blood.

PRECAUTIONS

Optochine Test cannot be classified as being hazardous according to the current legislation. Optochine Test is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store Optochine Test at $-20^\circ\text{C}/+8^\circ\text{C}$ in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, Optochine Test and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

- https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/394193/TP_25i3.pdf
- J. Lab. Clin. Med., 49: 641, 1957.
- J. Clin. Path., 8: 58, 1955.
- Serological Studies on Pneumococci, Munksgaard, Copenhagen, Oxford University Press, London, 1943.
- J. Exp. Med. 22: 269, 1915.

PRESENTATION

Product	Ref.	Test
Optochine Test	9501	100

TABLE OF SYMBOLS

IVD	In Vitro Diagnostic Medical Device	Do not reuse
REF	Catalogue number	Fragile, handle with care
Manufacturer	Manufacturer	Contains sufficient for <n> tests
Use by	Use by	Caution, consult accompanying documents
Temperature limitation	Temperature limitation	LOT Batch code

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LIOFILCHEM® S.r.l.