STAPH LATEX KIT

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Rev. 0/ 10.12.2009

STAPH LATEX KIT

Rapid agglutination test for direct identification of *Staphylococcus aureus*.

DESCRIPTION

STAPH LATEX KIT is a rapid latex agglutination test on card for direct identification of *Staphylococcus aureus*.

STAPH LATEX KIT can be used both for the identification of colonies grown on plate and for the confirmation of *Staphylococcus aureus* directly from 3-STA well of the PATHOGENIC SYSTEM *double* (ref. 71678).

PACKAGE CONTENT

- **1 LATEX REAGENT (2.5 mL):** dropper bottle containing latex particles coated with fibrinogen (to which coagulase binds) and IgG (which binds with Protein A).
- **1 POSITIVE CONTROL (1 mL):** dropper bottle containing positive control composed of inactivated preparation of *S. aureus*.
- 10 x 6 wells disposable cards.
- 2 x 25 disposable mixing sticks.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

- Disposable pipettes
- Disposable loops
- Physiological solution

METHOD PRINCIPLE

Latex particles in **STAPH LATEX KIT** are individually coated with fibrinogen (to which coagulase binds) and IgG (which binds with Protein A). Mixed with a suspension of *S. aureus*, latex particles aggregate to give visible agglutination.

TEST PROCEDURE

- Allow the reagents and positive control to reach room temperature.
- Dispense one drop of physiological solution on an agglutination card. Using a disposable loop remove one well-isolated colony from an agar plate and emulsify it in the drop of physiological solution so to obtain a heavy and homogeneous suspension. Suspensions must only be made from colonies with morphologies resembling *Staphylococcus aureus*.
- In order to confirm the isolation on **PATHOGENIC SYSTEM double** take, using a disposable pipette, one drop from the suspension in the well 3- STA and drop it on an agglutination card.
- Dispense one drop of **STAPH LATEX REAGENT** on the agglutination card next to the bacterial suspension, and one drop on the well with the positive control. **Do not allow the dropper to touch bacterial suspension.**
- Mix the latex reagent and the bacterial suspension with a clean mixing stick and rock the slide gently two or three times. Excessive rocking of the slide is not necessary. Examine for agglutination within a maximum of 2 minutes.

INTERPRETATION OF RESULTS

Agglutination within 2 minutes is a positive result and indicates the presence of *Staphylococcus aureus* in the sample. Absence of agglutination indicates that coagulase/ Protein A positive *Staphylococcus aureus* is not present in the test culture. Reactions that are "curdy" or "stringy" may not be true agglutination.

LIMITATION OF USE

1. Results should be interpreted in the context of all available clinical and laboratory information.
2. Test only pure, single colonies since mixed colonies may give erroneous results.
3. Cultures older than 30 hours may autoagglutinate.
4. Rough strains of *Staphylococcus* may cause false positive reactions. These strains are rare and distinguishable from smooth strains by their colonial morphology. If suspected, these can be confirmed by emulsifying in a drop of saline and examining carefully for a smooth suspension.
6. Stringy reactions on the slide may not be true positive reactions and further biochemical tests are required.
7. All coagulase positive strains of *Staphylococcus* will react with **STAPH LATEX REAGENT** and *S. aureus* will therefore not be distinguishable from *S. intermedius* and *S. hyicus*. However, the latter two strains are infrequently isolated from human sources and are more commonly found in animals or as saprophytes.

PRECAUTION

STAPH LATEX KIT is designed for *in vitro* diagnostic use and must be used only by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store **STAPH LATEX KIT** at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions **STAPH LATEX KIT** will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use **STAPH LATEX KIT** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	REF	Σ
STAPH LATEX KIT	96017	50 tests

TABLE OF SYMBOLS

	<i>In vitro</i> diagnostic medical device		Batch code
	Catalogue number		Fragile, handle with care
	Manufacturer		Content sufficient for <n> tests
	Use by		Warning, see instructions for use
	Temperature limits		Do not reuse

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net

Rev. 0/ 10.12.2009