Rapid agglutination test for direct identification of *Campylobacter jejuni*.

DESCRIPTION

CAMPYLOBACTER LATEX KIT is a rapid agglutination test on card for direct identification of *Campylobacter jejuni*. The test can be used for both the identification of colonies grown on agar plates and the confirmation of *Campylobacter jejuni* directly from the relevant well in various Liofilchem identification panels.

PACKAGING

- 1 CAMPYLOBACTER LATEX REAGENT (2.5 mL): dropper bottle containing latex particles coated with rabbit antibodies raised against *Campylobacter jejuni*.
- 1 POSITIVE CONTROL (1 mL): dropper bottle containing positive control composed of polyvalent inactivated antigens of *Campylobacter jejuni*.
- 10 disposable cards each with 6 wells.
- 2 racks each with 25 disposable mixing sticks.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

Disposable pipettes and loops, Physiological solution, Timer.

METHOD PRINCIPLE

Latex particles in CAMPYLOBACTER LATEX REAGENT are individually coated with polyvalent antisera against a wide range of *Campylobacter* antigens. Mixed with a suspension of *Campylobacter* cells, latex particles aggregate to give visible agglutination.

TEST PROCEDURE

1. Bring reagents and positive control to room temperature.
2. From colonies on agar plate: dispense one drop of physiological solution onto an agglutination card. Using a disposable loop remove one well-isolated colony from an agar plate and emulsify it in the drop of physiological solution so to obtain a heavy and homogeneous suspension. Suspensions must only be made from colonies with morphologies resembling *Campylobacter jejuni*.
From Liofilchem Systems: using a disposable pipette take one drop of culture from a relevant well and dispense onto an agglutination card.
3. Dispense one drop of POSITIVE CONTROL into a well of the agglutination card next to the well containing the microbial suspension.
4. Dispense one drop of CAMPYLOBACTER LATEX REAGENT onto the well containing the test suspension and one drop onto the well containing the positive control. Do not allow the dropper to touch the bacterial suspension.
5. Mix the latex reagent both with the bacterial suspension and the positive control using the sticks and rock the card gently two or three times. Excessive rocking of the card is not necessary. Examine for agglutination within a maximum of 2 minutes. If there is agglutination it will be clearly visible to the naked eye.

INTERPRETATION OF RESULTS

Agglutination within 2 minutes is a positive result and indicates the presence of *Campylobacter jejuni* in the sample. Absence of agglutination indicates that *Campylobacter jejuni* is not present in the test culture. Reactions that are "curdy" or "stringy" may not be true agglutination.

LIMITATION OF USE

- Results should be interpreted in the context of all available clinical and laboratory information.
- Non-specific agglutination of the Control Latex reagent does not preclude the presence of *Campylobacter* but the result has to be reported as inconclusive. The sample should be tested by an alternative method.
- Very low numbers of *Campylobacter* may result in a negative test result. 48-hour cultures should be used to maximise growth of the bacteria.

PRECAUTION

CAMPYLOBACTER LATEX KIT is designed for in vitro diagnostic use and must be used only by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store CAMPYLOBACTER LATEX KIT at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions CAMPYLOBACTER LATEX KIT will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After the use CAMPYLOBACTER LATEX KIT, and the material that comes into contact with the sample must be decontaminated and disposed according to the techniques used in the laboratory for decontamination and disposal of potentially infectious materials.

PRESENTATION

Product	REF	Σ
CAMPYLOBACTER LATEX KIT	96143	50 tests

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handel with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM® S.r.l.