LISTERIA LATEX KIT

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Rev. 0.1 / 10.12.2009

LISTERIA LATEX KIT

Rapid agglutination test for direct identification of *Listeria spp.*

DESCRIPTION

LISTERIA LATEX KIT is a rapid latex agglutination test on card for direct identification of *Listeria spp.*, both from culture media in plate, and from enrichment broths.

PACKAGE CONTENT

- **1 LATEX REAGENT (2.5 mL):** dropper bottle containing latex particles coated with antibodies raised against *Listeria spp.*
- **1 POSITIVE CONTROL (1 mL):** dropper bottle containing positive control composed of polyvalent inactivated antigenes of *Listeria spp.*
- 10 x 6 wells disposable cards.
- 2 x 25 disposable mixing sticks.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

- Disposable pipettes
- Disposable loops
- Physiological solution

METHOD PRINCIPLE

Latex particles in **LISTERIA LATEX KIT** are individually coated with polyvalent antisera against a wide range of *Listeria* antigens. Mixed with a suspension of *Listeriae*, latex particles aggregate to give visible agglutination.

TEST PROCEDURE

- Allow the reagents and positive control to reach room temperature.
- **METHOD FROM CULTURE MEDIUM IN PLATE:** dispense one drop of physiological solution on an agglutination card. Using a disposable loop remove one well-isolated colony from an agar plate and emulsify it in the drop of physiological solution so to obtain a heavy and homogeneous suspension. Suspensions must only be made from colonies with morphologies resembling *Listeria spp.*
- **METHOD FROM ENRICHMENT BROTHS:** using a disposable pipette, take one drop from the suspension in enrichment broth and drop it on an agglutination card.
- Dispense one drop of **POSITIVE CONTROL** on a well of the agglutination card.
- Dispense one drop of **LISTERIA LATEX REAGENT** on the agglutination card next to the bacterial suspension, and one drop on the well with the positive control. **Do not allow the dropper to touch bacterial suspension.**
- Mix the latex reagent and the bacterial suspension with a clean mixing stick and rock the slide gently two or three times. Excessive rocking of the slide is not necessary. Examine for agglutination within a maximum of 2 minutes.

INTERPRETATION OF RESULTS

Agglutination within 2 minutes is a positive result and indicates the presence of *Listeria spp.* in the sample. Absence of agglutination indicates that *Listeria spp.* is not present in the test culture. Reactions that are "curdy" or "stringy" may not be true agglutination.

LIMITATION OF USE

1. Results should be interpreted in the context of all available clinical and laboratory information.
2. Identification with **LISTERIA LATEX KIT** is presumptive and all positive results should be confirmed by further identification tests and serotyping of pure cultures.

PRECAUTION

LISTERIA LATEX KIT is designed for use in laboratory by properly trained operators only, and with approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store **LISTERIA LATEX KIT** at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions **LISTERIA LATEX KIT** will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use **LISTERIA LATEX KIT** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	REF	
LISTERIA LATEX KIT	96146	50 tests

TABLE OF SYMBOLS

	Batch code		Do not reuse
	Catalogue number		Fragile, handle with care
	Manufacturer		Content sufficient for <n> tests
	Use by		Warning, see instructions for use
	Temperature limits		

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net

Rev. 0.1 / 10.12.2009