## Testseren Anti-Yersinia enterocolitica

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## Test Sera Anti-Yersinia enterocolitica

### Instructions for professional use

#### Intended use

The test sera are used to identify the serovar of *Yersinia enterocolitica* strains isolated from test material of human and animal origin by means of slide agglutination.

#### Principle of the test

If the *Yersinia enterocolitica* strain contains an antigen covered by the test serum's range of specificity, this antigen will be bound with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

#### Composition

The test sera are absorbed sera from immunised rabbits.  
Preservative: Sodium azide (NaN<sub>3</sub>) 0.9 mg/ml

#### Available specificities

	REF
Anti-Yersinia enterocolitica O 3	TS 1701
Anti-Yersinia enterocolitica O 5	TS 1704
Anti-Yersinia enterocolitica O 8	TS 1705
Anti-Yersinia enterocolitica O 9	TS 1703
Anti-Yersinia enterocolitica O 27	TS 1706

#### Shelf life and storage

The test sera are lyophilized. They must be dissolved in 1 ml distilled water before use and properly sealed using the enclosed pipette.

When stored at 2–8 °C, the unopened test sera remain stable up to the expiration date specified on the label. Once opened and dissolved, they remain stable for at least 18 months when stored at 2–8 °C. Do not use after the expiration date printed on the label.

Occasionally, the test sera may become turbid in the absence of microbial contamination. This has no effect on efficacy; such turbidity can be resolved by centrifugation or filtration.

#### Warnings and precautions

The test sera are biological materials that should be treated as potentially infectious materials and handled accordingly.

As they contain sodium azide, contact with the skin and mucous membranes must be avoided. If contact occurs, rinse the affected area with plenty of water.

Since the slide agglutination test involves working with native pathogenic agents (risk of infection!), all applicable work protection procedures must be followed.

#### Materials and equipment not supplied

Glass slides, stirring rods, distilled water, physiological saline (9 g/L NaCl), disposal containers for infectious material, pipettes, pipette tips.

#### Test material and methodology

Transfer a small amount of bacterial mass from a single colony onto a slide (the slide should be placed on a dark surface). Mix well with a drop of test serum (ca. 25 µl) until a homogenous, slightly milky suspension is obtained. In rare cases, selective culture media may impair the agglutinability of the bacteria. This factor can be eliminated by removing the bacteria from the nutrient or blood agar or from the Kligler medium.

Read the reaction with the naked eye by holding the slide in front of a light source against a black background while tilting it back and forth.


To test for spontaneous agglutination, a negative control (NC) must be run with the samples using physiological saline instead of test serum.

#### Interpretation of the results of the analysis


The test results are valid only if the negative control (NC) remains milky-opaque.

**Positive:** Visible agglutination after the sample has been swirled back and forth 1 to 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is added. In

a weakly positive result, agglutination only appears after the slide has been swirled 10 to 20 times.

**Negative:** The result is negative (NR) if the suspension remains milky-opaque or if the reaction does not start before the slide has been swirled back and forth more than 20 times.

#### Limits of the procedure

The test sera react with *Yersinia enterocolitica* strains containing antigens of the specificity stated in the declaration.

Anti-Yersinia enterocolitica O 3, O 5, O 8 and O 9 test sera permit identification of *Yersinia enterocolitica* isolates belonging to group O. Anti-Yersinia enterocolitica O 27 test serum makes it possible to test O:5-positive strains for the presence or absence of O:27 antigen.

Serotyping alone does not permit an assessment of the pathogenicity of the examined strains in humans. For instance, not all O:5,27 or O:8 serovars are necessarily pathogenic. The results of biochemical testing (biotype identification) and/or phenotypic or genotypic detection of the virulence plasmid must be considered.

Because the antigens may be similar or related, cross-reactions with other *Yersinia* spp. or with other genera of the Enterobacteriaceae family may occur in isolated cases.

#### Explanation of the symbols used

IVD	In Vitro Diagnostic Medical Device
LOT	Batch code (Lot)
REF	Catalogue number
	Temperature limitation
	Use by YYYY-MM (MM = end of month)

Date of revision: July 2005