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Rapid agglutination test for direct identification of *Salmonella* spp.

DESCRIPTION

SALMONELLA LATEX KIT is a rapid agglutination test on card for direct identification of *Salmonella* spp. The test can be used both for the identification of colonies grown on agar plates and for the confirmation of *Salmonella* spp. directly from the relevant well in various Liofilchem identification panels.

PACKAGING

- 1 SALMONELLA LATEX REAGENT (2.5 mL): dropper bottle containing latex particles coated with rabbit antibodies raised against *Salmonella* spp.
- 1 POSITIVE CONTROL (1 mL): dropper bottle containing positive control composed of polyvalent inactivated antigens of *Salmonella*.
- 10 disposable cards each with 6 wells.
- 2 racks each with 25 disposable mixing sticks.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

Disposable pipettes and loops, Physiological solution, Timer.

METHOD PRINCIPLE

Latex particles in SALMONELLA LATEX KIT are individually coated with polyvalent antisera against a wide range of *Salmonella* antigens. Mixed with a suspension of salmonellae, latex particles aggregate to give visible agglutination.

Note, latex test for salmonella contains polyvalent anti-flagella/surface antigen antisera. The non-motile *S. gallinarum* and *S. pullorum* may be detected as they retain the protein basal bodies which are the attachment points of flagella.

TEST PROCEDURE

- Bring reagents and positive control to room temperature.
- From colonies on agar plate: Dispense one drop of physiological solution onto an agglutination card. Using a disposable loop remove one well-isolated colony from an agar plate and emulsify it in the drop of physiological solution so to obtain a heavy and homogeneous suspension. Suspensions must only be made from colonies with morphologies resembling *Salmonella* spp.
- From Liofilchem Systems: using a disposable pipette take one drop of culture from a relevant well and dispense onto an agglutination card.
- Dispense one drop of POSITIVE CONTROL into a well of the agglutination card next to the well containing the microbial suspension.
- Dispense one drop of SALMONELLA LATEX REAGENT onto the well containing the test suspension and one drop onto the well containing the positive control. Do not allow the dropper to touch the bacterial suspension.
- Mix the latex reagent both with the bacterial suspension and the positive control using the sticks and rock the card gently two or three times. Excessive rocking of the card is not necessary. Examine for agglutination within a maximum of 2 minutes. If there is agglutination it will be clearly visible to the naked eye.

INTERPRETATION OF RESULTS

Agglutination within 2 minutes is a positive result and indicates the presence of *Salmonella* spp. in the sample. Absence of agglutination indicates that *Salmonella* spp. is not present in the test culture. Reactions that are "curdy" or "stringy" may not be true agglutination.

LIMITATION OF USE

- The clinician should be interpreted the results in the context of all available clinical and laboratory informations.
- Some rough strains of *Salmonella* cause a non-specific agglutination in saline alone, hence cannot be tested.
- It is possible SALMONELLA LATEX KIT cannot detect some host-specific strains. Some oxidase-positive organisms may give false positive reactions.
- Old stock cultures of Enterobacteriaceae on nutrient agar slopes may cause non-specific agglutination whereas old stocks of *Salmonella* may give false negative results. Fresh sub-cultures should be prepared for testing.
- Identification with SALMONELLA LATEX KIT is presumptive and all positive results should be confirmed by further identification tests and serotyping of pure cultures.
- *S. bergen* and *S. dublin* result negative with SALMONELLA LATEX KIT.

PRECAUTION

SALMONELLA LATEX KIT is designed for in vitro diagnostic use and must be used only by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store SALMONELLA LATEX KIT at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions SALMONELLA LATEX KIT will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After the use SALMONELLA LATEX KIT, and the material that comes into contact with the sample must be decontaminated and disposed according to the techniques used in the laboratory for decontamination and disposal of potentially infectious materials. The reagent contains sodium azide must be handled with care. Sodium azide may react with copper and lead plumbing giving rise to explosive azides. Dispose of rinse with plenty of water to prevent azide build-up.

PRESENTATION

Product	REF	Σ
SALMONELLA LATEX KIT	96151	50 tests

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM® S.r.l.