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### DESCRIPTION

**STREPTO B LATEX KIT** is a rapid latex agglutination slide test for Group B Streptococci identification according to Lancefield classification. **STREPTO B LATEX KIT** can be used either for identification of colonies from culture plates or for confirmation of *Streptococcus agalactiae* (Group B) directly from **23-STG** well of the system A.F. GENITAL SYSTEM (ref.74156).

### CONTENT OF THE PACKAGE

- **1 LATEX REAGENT (2.5 mL)**: bottle with dropper containing latex particle sensitised with rabbit antibodies specific to antigens of Group B Streptococci.
- **1 EXTRACTION ENZYME (10 mL)**: bottle containing freeze dried extraction enzyme to reconstitute by adding 10 mL of distilled water.
- **1 POSITIVE CONTROL (1 mL)**: bottle with dropper containing positive control containing inactivated polyvalent antigenic extracts to Group B Streptococci.
- 10 disposable cards with 6 wells.
- 2 rack with disposable sticks for mixing

### NECESSARY ITEMS NON INCLUDED

Disposable pipettes, tubes and loops.  
Physiological solution. Distilled sterile water.

### PRINCIPLE OF THE METHOD

Latex particle are individually sensitised with rabbit antibodies specific to Streptococcal antigens of groups B. Streptococci colonies or bacterial suspension to test are incubated in an enzymatic suspension to extract antigen. The extract enzyme preparation is tested on a reaction card. In the presence of homologous antigen, particles will aggregate to deliver visible agglutination.

### TEST PROCEDURE

- Allow the **LATEX REAGENT** and **POSITIVE CONTROL** to reach room temperature.
- Reconstitute the bottle **EXTRACTION ENZYME** with 10 mL of distilled water and gently mix to ensure complete reconstitution.
- Dispense 0.2 mL of **EXTRACTION ENZYME** into a test tube.
- Add streptococci colonies picked from an agar surface or a drop of bacterial suspension from the **23-STG** of **A.F. GENITAL SYSTEM** (ref.74156) and emulsify gently.
- Incubate the tube for 10-15 minutes at 37°C.
- Shake the tube after the first 5 minutes of incubation and shake vigorously prior to testing to obtain even suspension of antigen.
- Dispense one drop of **LATEX REAGENT** separately into one of three circles on a reaction card.
- Transfer one drop of well mixed extract (or of **POSITIVE CONTROL**) next to the drop of **LATEX REAGENT**.
- Mix the content using mixing stick and spread the liquid to cover the area of the circle.
- Slowly and gently, rock and rotate the reaction card to mix the reagents for a maximum of one minute.
- Watch the card for agglutination. If present, agglutination should be clearly visible with the naked eye.

### INTERPRETATION OF RESULTS

#### Reactions from colonies grown on plate

In case of positive result the latex particles aggregate into visible clumps. If extract contains high quantity of antigens, the agglutination may be very rapid giving large clumps.

With weaker extracts agglutination may take longer to appear and give smaller clumps.

In case of negative result the latex particles retain their original milky appearance without any significant aggregation.

### Reactions from well 23- STG A.F. GENITAL SYSTEM (ref. 74156)

WELL COLOUR 23-STG	LATEX AGGLUTINATION REACTION	INTERPRETATION
GREEN	POSITIVE AGGLUTINATION	Presence of <i>Streptococcus agalactiae</i> (Group B)
GREEN	NEGATIVE AGGLUTINATION	Absence of <i>Streptococcus agalactiae</i> (Group B)
YELLOW	-	Absence of <i>Streptococcus agalactiae</i> (Group B)

### LIMITATION OF USE

Results must be evaluated in the light of other available clinical and laboratory information. Accurate results depend on testing an appropriate amount of growth.

### STORAGE

Store all reagents at 2-8°C. Do not freeze. Under these conditions the reagents will be usable until the expiry date.

Extraction Enzyme is stable for 3 months after reconstitution if stored at 2-8°C. To prolong the life of the enzyme, it may be dispensed into suitable test tubes in 0,4 mL volumes and stored frozen, at -20°C or below when it will be stable for 6 months.

Enzyme should not be frozen and thawed more than once.

### WARNINGS AND PRECAUTIONS

**STREPTO B LATEX** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

### DISPOSAL OF USED MATERIAL

After use, all materials that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

The reagent in this kit contain 0,099% sodium azide as preservative which should be handled with care. Sodium azide can react with lead and copper plumbing to form explosive azides. Upon disposal of reagents, flush with copious quantities of water to prevent azide build-up.

### PRESENTATION

Product	Ref.	Package
STREPTO B LATEX KIT	96153	50 test

### REFERENCES

- Facklam, R.R. (1985). Serologic identification of *streptococci* : how useful is serologic grouping? *Clin Microbiol Newsletter* 7, 91-94.
- Facklam, R.R. and Smith, P.B. (1976). The Gram positive cocci. *Human Pathology* 7, 187-194.
- Parker, M.T. and Ball, L.C. (1976). Streptococci and Aerococci associated with infection in man. *J. Med Microbiol* 9, 275-302.

### TABLE OF SYMBOLS

	In Vitro Diagnostic Medical Device		Fragile, handle with care
	Batch code		Do not reuse
	Catalogue number		Temperature limitation
	Manufacturer		Contains sufficient for <n> tests
	Use by		Caution, consult accompanying documents

**Liofilchem s.r.l.**

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39. 0858930745 - Fax +39 0858930330  
Web site: <http://www.liofilchem.net> E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

Rev. 0/ 15.11.2009