STREPTO A LATEX KIT

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Rev. 0/17.02.2010

STREPTO A LATEX KIT

Rapid agglutination test for direct identification of Group A Streptococcus.

DESCRIPTION

STREPTO A LATEX KIT is a rapid latex agglutination test on card for direct identification of Group A Streptococcus, both from culture media in plate, and from enrichment broths.

PACKAGE CONTENT

- **1 LATEX REAGENT (2.5 mL):** dropper bottle containing latex particles coated with antibodies raised against Group A Streptococcus.
- **1 EXTRACTION ENZYME (10 mL):** bottle containing freeze dried extraction enzyme to reconstitute by adding 10 mL of distilled water.
- **1 POSITIVE CONTROL (1 mL):** dropper bottle containing positive control composed of polyvalent inactivated antigens of Group A Streptococcus.
- 10 x 6 wells disposable cards.
- 2 x 25 disposable mixing sticks.

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

- Disposable pipettes.
- Disposable loops.
- Physiological solution.

METHOD PRINCIPLE

Latex particles in **STREPTO A LATEX KIT** are individually coated with antisera against a wide range of Group A Streptococcus antigens. Mixed with a suspension of these streptococci, latex particles aggregate to give visible agglutination.

TEST PROCEDURE

- Allow the reagents and positive control to reach room temperature.
- Reconstitute the bottle **EXTRACTION ENZYME** with 10 mL of distilled water and gently mix to ensure complete reconstitution.
- Dispense 0.2 mL of **EXTRACTION ENZYME** into a test tube.
- Add streptococci colonies picked from an agar surface and emulsify gently.
- Incubate the tube for 10-15 minutes at 37°C.
- Shake the tube after the first 5 minutes of incubation and shake vigorously prior to testing to obtain even suspension of antigen.
- Dispense one drop of **LATEX REAGENT** separately into one of three circles on a reaction card.
- Transfer one drop of well mixed extract (or of **POSITIVE CONTROL**) next to the drop of **LATEX REAGENT**.
- Mix the content using mixing stick and spread the liquid to cover the area of the circle.
- Slowly and gently, rock and rotate the reaction card to mix the reagents for a maximum of one minute.
- Watch the card for agglutination. If present, agglutination should be clearly visible with the naked eye.

INTERPRETATION OF RESULTS

Agglutination within 2 minutes is a positive result and indicates the presence of Group A Streptococcus in the sample. Absence of agglutination indicates that Group A Streptococcus is not present in the test culture. Reactions that are "curdy" or "stringy" may not be true agglutination.

LIMITATION OF USE

1. Results should be interpreted in the context of all available clinical and laboratory information.
2. Identification with **STREPTO A LATEX KIT** is presumptive and all positive results should be confirmed by further identification tests and serotyping of pure cultures.

PRECAUTION

STREPTO A LATEX KIT is designed for *in vitro* diagnostic use and must be used only by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store **STREPTO A LATEX KIT** at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions **STREPTO A LATEX KIT** will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use **STREPTO A LATEX KIT** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	REF	Σ
STREPTO A LATEX KIT	96154	50 tests

TABLE OF SYMBOLS

IVD	<i>In vitro</i> diagnostic medical device	LOT	Batch code
REF	Catalogue number		Fragile, handle with care
	Manufacturer	Σ	Content sufficient for <n> tests
	Use by		Warning, see instructions for use
	Temperature limits		Do not reuse

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net

Rev. 0/17.02.2010