## CRYSTAL VIOLET AGAR (Sheep Blood 5%)

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## CRYSTAL VIOLET AGAR (Sheep Blood 5%)

Selective medium for streptococci isolation.

TYPICAL FORMULA	(g/l)
Casein Peptone Plus	4.6
Soy Peptone	3.0
Meat Peptone	7.4
Sodium Chloride	3.0
Corn Starch	0.6
Agar	10.0
Sodium Azide	0.2
Crystal Violet	4.0 ml
Sheep Defibrinated Blood	50.0 ml
Final pH 7.3 ± 0.2	

### DESCRIPTION

CRYSTAL VIOLET AGAR (*Sheep Blood 5%*) is a selective medium suitable for isolation and growth of streptococci from clinical samples.

### PRINCIPLE

Casein peptone plus, soy peptone e meat peptone provide for proteinic requirements of microorganisms. Sodium chloride maintains the osmotic balance of the medium. Starch is a source of energy while agar is the solidifying agent. Sodium azide is inhibitory against Gram-negative microorganisms. Crystal violet inhibits Gram-positive microorganisms except streptococci. Sheep defibrinated blood supplies further nutritional factors for microorganisms growth.

### TECHNIQUE

Heavily inoculate by streaking the specimen onto the surface of the medium with a sterile loop. Incubate at 36±1°C for 36-48 hours, in 10% CO<sub>2</sub> atmosphere.

### INTERPRETATION OF RESULTS

At the end of the incubation time the colored colonies are easily recognized; longer incubation or storage in the refrigerator tend to cause all the colonies to become colored. Strains producing violet or golden growths are likely to be pathogenic while about 95% of those producing white or pale violet growths are not pathogenic.

### STORAGE AND TRANSPORT CONDITIONS

6-12°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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CRYSTAL VIOLET AGAR (Sheep Blood 5%)

### PRESENTATION

Ready to use plates (90 mm) containing 22+/-1 ml of medium

### STORAGE

6-12°C

### PACKAGING

Ref.	Content	Packaging
10042	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>
10042*	100 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 10 piles (10 x 10 plates) in cardboard box</li> </ul>

### pH OF THE MEDIUM

7.3 ± 0.2

### USE

CRYSTAL VIOLET AGAR (*Sheep Blood 5%*) is a selective medium suitable for isolation and growth of streptococci

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Purple, clear to slightly opalescent

### SHELF LIFE

45 days

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at 22 ± 1°C, in aerobiosis  
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control  
Inoculum for productivity: 10-100 CFU/ml  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> CFU/ml  
Inoculum for specificity: ≤10<sup>4</sup> CFU/ml  
Incubation Conditions: 36-48 hours at 36 ± 1°C, in 10% CO<sub>2</sub> atmosphere

Microorganism		Growth
<i>Streptococcus pyogenes</i>	ATCC 19615	good
<i>Enterococcus faecalis</i>	ATCC 19433	good
<i>Escherichia coli</i>	ATCC 25922	inhibited
<i>Staphylococcus aureus</i>	ATCC 25923	inhibited

### TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

**LIOFILCHEM® S.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY  
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net