## Yersinia Selective Agar

Differential and selective medium for the isolation of *Y. enterocolitica* from clinical and nonclinical specimens, according to ISO 10273.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Gelatin	17.0
Enzymatic Digest of Casein and Animal Tissues	3.0
Yeast Extract	2.0
Sodium Chloride	1.0
Sodium Pyruvate	2.0
Magnesium Sulfate	0.01
Mannitol	20.0
Sodium Deoxycholate	0.5
Crystal Violet	0.001
Neutral Red	0.03
Agar	14.0
Cefsulodin	0.015
Irgasan	0.004
Novobiocin	0.0025
Final pH 7.4 ± 0.2 at 25°C	

### DESCRIPTION

Yersinia Selective Agar is a medium used for the selective isolation and differentiation of *Yersinia enterocolitica*.

This medium (CIN agar) is recommended by ISO 10273 for the examination of food and animal feed stuffs as well as environmental samples in the area of food production and food handling.

### PRINCIPLE

Enzymatic digest of gelatin and enzymatic digest of casein and animal tissues provide amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Sodium pyruvate and magnesium sulfate stimulate organisms growth. Mannitol is the carbohydrate which allows to differentiate between mannitol fermenting and non-fermenting bacteria. Sodium deoxycholate and crystal violet inhibit Gram-positive bacteria. Neutral red is the pH indicator. Agar is the solidifying agent. Cefsulodin, irgasan (triclosan) and novobiocin are the selective agents inhibiting the growth of most Gram-negative enteric bacteria.

### TECHNIQUE

Inoculate the specimen onto the medium by either direct plating or pour plating (\*). Incubate aerobically at 30 ± 1°C for 18-24 h.

(\*) The ISO method for the detection of presumptive pathogenic *Yersinia enterocolitica* recommends to first perform enrichment in Peptone, Sorbitol and Bile Salts (PSB) Broth for 48-72 hours at 22-25°C with agitation, or 5 days without agitation.

### INTERPRETATION OF RESULTS

Organisms capable of fermenting mannitol cause a localized pH reduction, forming colonies with red centre surrounded by a transparent border (characteristic "bull's-eye" colony). Organisms that do not ferment mannitol form colorless, translucent colonies. Some strains of *Serratia*, *Citrobacter* and *Enterobacter* may give a colonial morphology resembling *Yersinia enterocolitica*. Final identification should be confirmed by standard biochemical tests.

### STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

### REFERENCES

- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 10273:2003. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of presumptive pathogenic *Yersinia enterocolitica*.
- Schieman, D.A. (1979) Synthesis of a selective agar medium for *Yersinia enterocolitica*. Can. J. Microbiol. 25:1298-1304.
- Schieman, D.A. (1980) *Yersinia enterocolitica*: Observation on some growth characteristics and response to selective agents. Can. J. Microbiol. 43:14-27.
- Devenish, J.A., and D.A. Schieman (1981) An abbreviated scheme for identification of *Yersinia enterocolitica* isolated from food enrichments on CIN (cefsulodin-irgasan-novobiocin) agar. Can. J. Microbiol. 27:937-941.

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Yersinia Selective Agar

### PRESENTATION

Ready to use plates (90 mm) containing  $22 \pm 1$  ml of medium

### STORAGE

10-25°C

### PACKAGING

Ref.	Content	Packaging
10052	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>
10052*	100 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 10 x 10 plates in cardboard box</li> </ul>

### pH OF THE MEDIUM

$7.4 \pm 0.2$

### USE

Yersinia Selective Agar is a differential and selective medium used for the isolation of *Yersinia enterocolitica* from clinical specimens and other types of samples, according to ISO 10273

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Slightly opalescent, reddish-orange

### SHELF LIFE

180 days

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at  $22 \pm 2^\circ\text{C}$ , in aerobiosis  
7 days at  $35 \pm 2^\circ\text{C}$ , in aerobiosis
- Microbiological control  
Inoculum for productivity: 50-100 CFU  
Inoculum for selectivity:  $10^4$ - $10^6$  CFU  
Incubation Conditions: 18-24 h at  $30 \pm 1^\circ\text{C}$ , in aerobiosis

Microorganism	WDCM	Growth	Colony Appearance
<i>Yersinia enterocolitica</i>	WDCM 00038	Good	Colonies with red center
<i>Escherichia coli</i>	WDCM 00012	Partially to totally inhibited	---
<i>Staphylococcus aureus</i>	WDCM 00034	Inhibited	---

### TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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1. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 10273:2003. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of presumptive pathogenic *Yersinia enterocolitica*.
3. Schieman, D.A. (1979) Synthesis of a selective agar medium for *Yersinia enterocolitica*. Can. J. Microbiol. 25:1298-1304.
4. Schieman, D.A. (1980) *Yersinia enterocolitica*: Observation on some growth characteristics and response to selective agents. Can. J. Microbiol. 43:14-27.
5. Devenish, J.A., and D.A. Schieman (1981) An abbreviated scheme for identification of *Yersinia enterocolitica* isolated from food enrichments on CIN (cefusulodin-irgasan-novobiocin) agar. Can. J. Microbiol. 27:937-941.

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**LIOFILCHEM® S.r.l.**