Schaedler KKV Agar (Sheep Blood 5%)

Selective medium for the isolation and cultivation of Gram-negative anaerobic bacteria from clinical specimens.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	8.2
Peptic Digest of Animal Tissue	2.5
Papaic Digest of Soybean Meal	1.0
Glucose	5.8
Yeast Extract	5.0
Sodium Chloride	1.7
L-Cystine	0.4
Dipotassium Phosphate	0.8
Hemin	0.01
Vitamin K1	0.01
Kanamycin Sulphate	0.1
Vancomycin	0.0075
Defibrinated Sheep Blood	50.0 ml
Agar	13.5
Final pH 7.6 ± 0.2	

DESCRIPTION

Schaedler KKV Agar (Sheep Blood 5%) is a selective medium used for the isolation and cultivation of Gram-negative anaerobic bacteria from clinical specimens.

PRINCIPLE

Pancreatic digest of casein, peptic digest of animal tissue and papaic digest of soybean meal provide amino acids, nitrogen, carbon, vitamins and minerals. Glucose is an energy source. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. L-cystine is a reducing agent. Dipotassium phosphate is the buffer. Sheep blood and hemin stimulate the growth of fastidious organisms. Vitamin K1 enhances the recovery of certain anaerobes. Kanamycin inhibits protein synthesis in susceptible microorganisms, while vancomycin inhibits Gram-positive bacteria by interfering with cell wall synthesis. Agar is the solidifying agent

TECHNIQUE

Inoculate the plates by streaking the specimen directly onto the agar surface. Incubate at 37°C for 48-72 hours in anaerobic atmosphere.

INTERPRETATION OF RESULTS

Observe for the growth of colonies.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. MacFaddin, J.D. (1985) Media for isolation-cultivation-identification-identification-maintenance of medical bacteria, p. 695-699, vol. 1. Williams & Wilkins, Baltimore, MD.
2. Schaedler, R.W., R. Dubos, and R. Costello (1965) The development of bacterial flora in the gastrointestinal tract of mice. J. Exp. Med. 122: 59.

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Schaedler KKV Agar (Sheep Blood 5%)

PRESENTATION

Ready to use plates (90 mm) containing 22±1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10065	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
10065*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 x 10 plates in cardboard box

pH OF THE MEDIUM

7.6 ± 0.2

USE

Schaedler KKV Agar (Sheep Blood 5%) is a selective medium used for the isolation and cultivation of Gram-negative anaerobic bacteria from clinical specimens

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Cherry red medium, opaque

SHELF LIFE

60 days

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
 - 7 days at 22 ± 1°C, in aerobiosis
 - 7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
 - Inoculum for productivity: 10-100 CFU/ml
 - Inoculum for selectivity: 10⁴-10⁵ CFU/ml
 - Incubation Conditions: 48-72 hours at 35 ± 2°C, in anaerobiosis

Microorganism		Growth
<i>Bacteroides fragilis</i>	ATCC® 25285	Good
<i>Bacteroides thetaiotaomicron</i>	ATCC® 29741	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Inhibited
<i>Escherichia coli</i>	ATCC® 25922	Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.