## Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L $\beta$ -NAD)

Medium recommended by EUCAST for antimicrobial susceptibility testing of *Streptococcus* spp, *Haemophilus* spp and several other fastidious organisms.

TYPICAL FORMULA	(g/L)
Meat Extract	2.0
Casamino Acids, Technical	17.5
Starch	1.5
$\beta$ -NAD	0.02
Agar	15.0
Horse Defibrinated Blood	50 mL
Final pH 7.3 $\pm$ 0.1	

### DESCRIPTION

Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L  $\beta$ -NAD) is recommended by EUCAST for antimicrobial susceptibility testing of *Streptococcus* spp (including *S. pneumoniae*), *Haemophilus* spp and several other fastidious organisms.

### PRINCIPLE

Meat extract and casamino acids supply amino acids and other nitrogenous substances, minerals, vitamins, carbon and other nutrients to support the growth of microorganisms. Starch acts as a protective colloid against toxic substances that may be present in the medium. Beside hydrolysis of starch during autoclaving provides a small amount of glucose which is a source of energy. Agar is the solidifying agent. Horse blood and  $\beta$ -NAD are a further growth supplementation for fastidious microorganisms.

### TECHNIQUE

- Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 mL of a suitable suspension medium.
- Compare the turbidity to the appropriate McFarland standard.
- Introduce a sterile swab into the broth and inoculate a plate of Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L  $\beta$ -NAD) by streaking the swab 2 or 3 times over the entire agar surface.
- Apply the antibiotic disc or the MIC Test Strip on the agar surface.
- Incubate plates in an inverted position under conditions appropriate for the microorganism.
- Measure the diameter of the inhibition halo using a pair of compasses for the antibiotic disc or read the MIC value where the edge of the inhibition ellipse intersects the strip.

### INTERPRETATION OF RESULTS

For antibiotic disc zone diameter measured around discs should be compared with those in the EUCAST breakpoint tables.

Results obtained with specific organisms may then be reported as resistant, intermediate or susceptible.

It is important to measure the diameter of the zone of growth inhibition, not of hemolysis inhibition.

MIC break points for defining susceptibility categories as provided by the EUCAST could be used for interpreting MIC values.

Results obtained with specific organisms may then be reported as resistant, intermediate or susceptible.

### STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### REFERENCES

- Current EUCAST Breakpoints and QC Tables.
- Media preparation for EUCAST disk diffusion testing and for determination of MIC values by the broth microdilution method. Version 4, June 2014.
- Barry, Garcia and Thrupp. (1970). Am. J. Clin. Pathol. 53: 149.
- Bauer et al. (1966) J. Clin. Pathol. 45:493-496.
- Mueller, J.H., and Hinton (1941) Proc. Soc. Exp. Biol. Med. 48: 330-333.

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Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L  $\beta$ -NAD)

### PRESENTATION

Ready-to-use plates with  $4 \pm 0.5$  mm agar depth

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
10132	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>
10132*	100 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 10 x 10 plates in cardboard box</li> </ul>
11132	10 plates (140 mm)	<ul style="list-style-type: none"> <li>• 5 plates in thermally soldered film</li> <li>• 2 x 5 plates in cardboard box</li> </ul>
12033	5 square plates (120 mm)	<ul style="list-style-type: none"> <li>• 5 plates in thermally soldered film and in cardboard box</li> </ul>

### pH OF THE MEDIUM

$7.3 \pm 0.1$

### USE

Mueller Hinton Fastidious Agar (Horse blood 5% + 20 mg/L  $\beta$ -NAD) is recommended by EUCAST for antimicrobial susceptibility testing of *Streptococcus* spp (including *S. pneumoniae*), *Haemophilus* spp and several other fastidious organisms

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Cherry red, opaque

### SHELF LIFE

60 days

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at  $22 \pm 2^\circ\text{C}$ , in aerobiosis  
7 days at  $35 \pm 2^\circ\text{C}$ , in aerobiosis
- Microbiological control  
Incubation conditions: set according to EUCAST and shown on the quality control certificate of each batch

#### Microorganism

#### Growth

<i>Streptococcus pneumoniae</i>	ATCC® 49619	Good
<i>Haemophilus influenzae</i>	NCTC 8468	Good
<i>Haemophilus influenzae</i>	ATCC® 49766	Good
<i>Haemophilus influenzae</i>	ATCC® 49247	Good
<i>Campylobacter jejunii</i>	ATCC® 33560	Good

### TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro Diagnostic Medical Device</i>		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for $<n>$ tests		Caution, consult instruction for use		Do not reuse

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**LIOFILCHEM® S.r.l.**