Blood Agar (Sheep Blood 7%)

Medium for isolation and cultivation of fastidious organisms, especially of pathogenic species, and detection of haemolytic reactions.

TYPICAL FORMULA	(g/l)
Proteose Peptone	15.0
Liver Digest	2.5
Yeast Extract	5.0
Sodium Chloride	5.0
Agar	15.0
Defibrinated Sheep Blood	70.0 ml
Final pH 7.2 ± 0.2	

DESCRIPTION

Blood Agar (Sheep Blood 7%) is a medium used for the cultivation of fastidious pathogenic organisms and detection of their haemolytic activity.

This medium is recommended in ISO 7932 and ISO 11290-1 for confirmation of *Bacillus cereus* and *Listeria monocytogenes*, respectively.

PRINCIPLE

Proteose peptone and liver digest provide amino acid, carbon, nitrogen, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Blood is included to provide additional growth factors for fastidious microorganisms and allows to demonstrate the haemolytic reactions.

TECHNIQUE

Inoculate the plates by the streak/spread method. Incubate under optimal conditions, usually 35 ± 2°C for 24-48 hours.

To determine the haemolytic activity, after obtaining a pure culture, streak, stab or spot the selected colonies onto the medium in a manner which allows good interpretation of the haemolytic reactions.

INTERPRETATION OF RESULTS

Examine the plates for growth and haemolytic reactions. Four types of haemolysis on blood agar media can be described:

1. α-haemolysis is the reduction of hemoglobin to methemoglobin in the medium surrounding the colony, causing a greenish decolorization of the medium;
2. β-haemolysis is the lysis of red blood cells, producing a clear zone surrounding the colony;
3. γ-haemolysis indicates no destruction of red blood cells and no change in the color of the medium;
4. δ-haemolysis indicates a partial lysis.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. ISO 7932:2003. Microbiology of food and animal feeding stuffs - Horizontal Methods for the enumeration of presumptive *Bacillus cereus*. Colony count technique at 30°C.
2. ISO 11290-1:1996. Microbiology of food and animal feeding stuffs - Horizontal Methods for the detection and *Listeria monocytogenes*. - Part 1 Detection method.
3. U.S. Food and Drug Administration. (1995). Bacteriological analytical manual, 8th ed., AOAC International, Gaithersburg, Md.

PRODUCT SPECIFICATIONS

NAME

Blood Agar (Sheep Blood 7%)

PRESENTATION

Ready to use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10142	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
10142*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 piles (10 x 10 plates) in cardboard box

pH OF THE MEDIUM

7.2 ± 0.2

USE

Blood Agar (Sheep Blood 7%) is a medium used for the cultivation of fastidious pathogenic organisms and detection of their haemolytic activity, according to ISO 7932 and ISO 11290-1

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Cherry red, opaque

SHELF LIFE

75 days

QUALITY CONTROL


- Control of general characteristics, label and print
- Sterility control
7 days at $22 \pm 2^\circ\text{C}$, in aerobiosis
7 days at $35 \pm 2^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Incubation Conditions: 18-24 hours at $35 \pm 2^\circ\text{C}$

Microorganism		Growth	Haemolysis
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	β -haemolysis
<i>Streptococcus pyogenes</i> *	ATCC® 19615	Good	β -haemolysis
<i>Streptococcus pneumoniae</i> *	ATCC® 6303	Good	α -haemolysis
<i>Neisseria meningitidis</i> *	ATCC® 13090	Good	γ -haemolysis
<i>Listeria monocytogenes</i>	ATCC® 19118	Good	β -haemolysis
<i>Bacillus cereus</i>	ATCC® 11778	Good	β -haemolysis
<i>Escherichia coli</i>	ATCC® 25922	Good	none
<i>Clostridium perfringens</i> **	ATCC® 13124	Good	β -haemolysis

*10% CO₂

**Anaerobic

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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1. ISO 7932:2003. Microbiology of food and animal feeding stuffs - Horizontal Methods for the enumeration of presumptive *Bacillus cereus*. Colony count technique at 30°C.
2. ISO 11290-1:1996. Microbiology of food and animal feeding stuffs - Horizontal Methods for the detection and *Listeria monocytogenes*. - Part 1 Detection method.
3. U.S. Food and Drug Administration. (1995). Bacteriological analytical manual, 8th ed., AOAC International, Gaithersburg, Md

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LIOFILCHEM® S.r.l.