Campylobacter Karmali Agar

Selective medium for detection of thermotolerant *Campylobacter* spp from clinical specimens and food.

TYPICAL FORMULA	(g/l)
Peptospecial	23.0
Corn Starch	1.0
Sodium Chloride	5.0
Activated Charcoal	4.0
Haemin	0.032
Agar	10.0
Vancomycin	0.02
Sodium Pyruvate	0.1
Cefoperazone	0.032
Cycloheximide	0.1
Final pH 7.4 ± 0.2	

DESCRIPTION

Campylobacter Karmali Agar is a selective medium used for the isolation and cultivation of *Campylobacter jejuni* and *Campylobacter coli* from clinical specimens and food.

PRINCIPLE

Peptospecial provides amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Corn starch and activated charcoal absorb and neutralize any toxic metabolites. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. The selectivity is due to the presence of the antimicrobial agents, Vancomycin inhibits Gram-positive bacteria, Cefoperazone suppress Gram-negative organisms other than *Campylobacter* spp and Cycloheximide inhibits fungal growth. Sodium pyruvate and haemin represent further growth supplementation for *Campylobacter* spp.

TECHNIQUE

Inoculate the medium by streaking the sample onto the agar surface. Incubate at 42 ± 1°C for 40-48 hours in microaerobic atmosphere.

INTERPRETATION OF RESULTS

Examine for growth of typical colonies of *Campylobacter* spp which appear greyish, often with a metallic sheen, are flat and moist, with a tendency to spread. Other forms of colonies may occur. Suspected colonies have to be submitted to biochemical and serological confirmatory tests.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Karmali, Simor, Roscoe, Fleming, Smith and Lane (1986) J. Clin. Microbiol. 23: 456.
2. Blaser, M.J., V. Berkowitz, F.M. Laforce (1979) Campylobacter enteritidis: Clinical and Epidemiologic Features. Ann. Intern. Med. 91:179-185.
3. Association of Official Analytical Chemists (1995) Bacteriological analytical manual, 8th Ed.

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Campylobacter Karmali Agar

PRESENTATION

Ready to use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10145	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.4 ± 0.2

USE

Campylobacter Karmali Agar is a selective medium used for the isolation and cultivation of *Campylobacter jejuni* and *Campylobacter coli* from clinical specimens and food

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Opaque, black

SHELF LIFE

4 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at $22 \pm 2^\circ\text{C}$, in aerobiosis
7 days at $35 \pm 2^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10^4 - 10^6 CFU
Incubation Conditions: 40-48 hours at $42 \pm 1^\circ\text{C}$ in atmosphere enriched with 5-10% carbon dioxide

Microorganism		Growth
<i>Campylobacter jejuni</i>	ATCC® 33291	Good
<i>Escherichia coli</i>	ATCC® 25922	Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	In vitro Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.