## Campylobacter Preston Agar

Selective medium for the isolation of *Campylobacter* spp, from clinical and non clinical samples.

TYPICAL FORMULA	(g/l)
Peptone	7.5
Enzymatic Digest of Casein	7.5
Enzymatic Digest of Liver	2.5
Yeast Extract	5.0
Sodium Chloride	5.0
Sodium Pyruvate	0.25
Iron(II) Sulfate	0.25
Sodium Metabisulfite	0.25
Agar	15.0
Lysed Horse Blood	50 ml
Polymyxin B Sulfate	5,000 IU
Trimethoprim	0.01
Rifampicin	0.01
Cycloheximide	0.1
Final pH 7.5 ± 0.2	

### DESCRIPTION

Campylobacter Preston Agar is a medium used for the selective isolation of *Campylobacter* spp from human, animal and environmental specimens.

### PRINCIPLE

Peptone, enzymatic digest of casein and enzymatic digest liver provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Sodium pyruvate, ferrous sulfate and sodium metabisulfite act as oxygen quencher. Agar is the solidifying agent. Lysed horse blood provides additional nutrients for fastidious bacteria. Antimicrobics suppress normal enteric flora and enhance the growth of *Campylobacter* spp.

### TECHNIQUE

Inoculate the plates by using the enrichment culture obtained in a suitable medium such as Bolton Broth (ref. 470340) or streak the sample directly over the agar surface to achieve single colonies. Incubate microaerobically at 37°C or 41.5 ± 1°C for 40-48 hours.

### INTERPRETATION OF RESULTS

Observe typical colonies of *Campylobacter* species which appear whitish, flat and moist, often with a metallic sheen, and with a tendency to spread. Other forms of colonies may occur.

For confirmation, subculture a suspected colony to a Columbia Blood Agar plate (ref. 11025) and examine the pure culture for morphology, motility, microaerobic growth at 25°C, aerobic growth at 41.5°C and the presence of oxidase.

### STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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Campylobacter Preston Agar

### PRESENTATION

Ready to use plates (90 mm) containing  $22 \pm 1$  ml of medium

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
10146	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>

### pH OF THE MEDIUM

$7.5 \pm 0.2$

### USE

Campylobacter Preston Agar is a medium used for the selective isolation of *Campylobacter* spp from all types of specimens

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Slightly opalescent, red-brown

### SHELF LIFE

4 months

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at  $22 \pm 2^\circ\text{C}$ , in aerobiosis  
7 days at  $35 \pm 2^\circ\text{C}$ , in aerobiosis
- Microbiological control  
Inoculum for productivity: 50-100 CFU  
Inoculum for selectivity:  $10^4$ - $10^6$  CFU  
Incubation Conditions: 40-48 hours at  $41.5 \pm 1^\circ\text{C}$ , in microaerobic atmosphere

#### Microorganism

*Campylobacter jejuni* ATCC® 29428  
*Escherichia coli* ATCC® 25922  
*Candida albicans* ATCC® 10231

#### Growth

Good  
Inhibited  
Inhibited

### TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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**LIOFILCHEM® S.r.l.**