CAMPYLOBACTER AGAR (Sheep blood 10%)

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1. Skirrow, M.D. 1977 *Campylobacter* enteritidis: A New Disease. Br. Med. J. 2:9-11.
2. Blaser, M.J., V. Berkowitz, F.M. Laforce, 1979. *Campylobacter* enteritidis: Clinical and Epidemiologic Features. Ann. Intern. Med. 91:179-185.
3. Association of Official Analytical Chemists. 1995. Bacteriological analytical manual, 8th Ed.

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CAMPYLOBACTER AGAR (Sheep blood 10%)

Selective medium for the isolation of *Campylobacter* spp.

TYPICAL FORMULA	(g/l)
Peptone	7.5
Tryptone	7.5
Enzymatic Digest of Liver	2.5
Potassium Phosphate Monobasic	1.0
Yeast Extract	5.0
Sodium Chloride	5.0
Agar	15.0
Defibrinated Sheep Blood	100 ml
Polymyxin B	2500 IU
Trimethoprim	5.0 mg
Amphotericin B	2.0 mg
Cefalotin	15.0 mg
Vancomycin	10.0 mg
Final pH 7.4 ± 0.2	

DESCRIPTION

CAMPYLOBACTER AGAR (Sheep blood 10%) is a selective medium for *Campylobacter* spp isolation.

PRINCIPLE

Peptone, tryptone and enzymatic digest of liver are a source of amino acids and proteins. Yeast extract is a source of group B vitamins. Sodium chloride maintains the osmotic balance of the medium. Potassium Phosphate is the buffer system that maintains constant the pH of the medium. Agar is the solidifying agent. The blood constitutes an additional source of nutrients for fastidious microorganisms. The addition of antimicrobial agents suppresses the growth of the normal microbial flora in fecal specimens, thereby facilitating isolation of *C. jejuni*.

TECHNIQUE

Inoculate the plates streaking the sample onto the surface of the medium, using a sterile loop. Incubate at 42°C for 24-48 hours in microaerobic atmosphere.

INTERPRETATION OF RESULTS

C. jejuni produces two types of colonies:

- colonies small, brown-gray, smooth, glossy, translucent edges.
- mucoid appearance colonies, translucent, grayish, with irregular borders.

STORAGE

6-12°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Skirrow, M.D. 1977 *Campylobacter* enteritidis: A New Disease. *Br. Med. J.* 2:9-11.
2. Blaser, M.J., V. Berkowitz, F.M. Laforce, 1979. *Campylobacter* enteritidis: Clinical and Epidemiologic Features. *Ann. Intern. Med.* 91:179-185.
3. Association of Official Analytical Chemists. 1995. *Bacteriological analytical manual*, 8th Ed

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

CAMPYLOBACTER AGAR (Sheep blood 10%)

PRESENTATION

90 mm ready plates containing 22+/-1 ml of medium

STORAGE

6-12°C

PACKAGING

Ref.	Content	Packaging
10148	20 plates	<ul style="list-style-type: none"> • 5 plates in thermally soldered film • 4 x 5 plates in cardboard boxes

pH OF THE MEDIUM

7.4 ± 0.2

USE

CAMPYLOBACTER AGAR (Sheep blood 10%) is a selective medium for *Campylobacter* spp isolation

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Cherry red medium, opaque

SHELF LIFE

75 days

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 CFU/ml
Inoculum for selectivity: 10⁴-10⁵ CFU/ml
Inoculum for specificity: ≤10⁴ CFU/ml
Incubation Conditions: 48 hours at 42°C, in microaerobic atmosphere

Microorganism

Escherichia coli ATCC 25922
Campylobacter jejuni ATCC 33291

Growth

Inhibited
Good

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@lioilchem.net