Vancomycin Screen Agar

Medium for detection of vancomycin-resistant enterococci (VRE) and *Staphylococcus aureus* (VRSA) from pure cultures.

TYPICAL FORMULA	(g/l)
Brain Heart Infusion Solids	17.5
Peptone	10.0
Glucose	2.0
Sodium Chloride	5.0
Disodium Phosphate	2.5
Vancomycin	0.006
Agar	15.0
Final pH 7.4 ± 0.2	

DESCRIPTION

Vancomycin Screen Agar is a medium used for testing vancomycin susceptibility in enterococci and *Staphylococcus aureus*.

This medium, Brain Heart Infusion (BHI) Agar with 6 mg/l of vancomycin, is recommended by the Clinical and Laboratory Standards Institute (CLSI document M07).

PRINCIPLE

The meat infusion and peptone provide amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Glucose is the fermentable carbohydrate. Disodium phosphate acts as buffer. Vancomycin is the selective agent included to detect resistance. Agar is the solidify agent.

TECHNIQUE

1. Prepare a standardized suspension (McFarland 0.5) of the test organism by the direct colony suspension method.
2. Dip a sterile cotton swab into the adjusted suspension.
3. Inoculate the plate by streaking the swab over the entire agar surface.
4. Incubate aerobically at 35 ± 2°C for a full 24 hours.

INTERPRETATION OF RESULTS

Greater than one colony or a film growth suggest reduced susceptibility (*S. aureus*) or resistance (enterococci) to vancomycin. Organisms not growing on this medium are considered sensitive.

Notes:

- There are insufficient data to recommend using this agar screen test for coagulase-negative staphylococci (CoNS).
- All staphylococci isolates that grow on this medium should be inspected for purity and be tested using broth microdilution (BMD) or other MIC method for confirmation.
- This screening test for vancomycin-resistant enterococci has been shown to be comparable in reliability to standard methods in detecting clinically significant resistance, and additional confirmatory tests are unnecessary.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. CLSI (2017) Performance Standards for Antimicrobial Susceptibility Testing; 27th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute.
2. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 7.0, 2017. <http://www.eucast.org>.
3. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 7.0, 2017. <http://www.eucast.org>.
4. CLSI M07-A10 (2015) Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically; Approved Standard - Tenth Edition. Wayne, PA: Clinical and Laboratory Standards Institute.
5. https://www.cdc.gov/hai/pdfs/vrsa-investigation-guide-05_12_2015.pdf

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Vancomycin Screen Agar

PRESENTATION

Ready-to-use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10407	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.4 ± 0.2

USE

Vancomycin Screen Agar is a medium used for detection of vancomycin resistance

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Slightly opalescent, amber

SHELF LIFE

3 months

QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control
 - 7 days at 22 ± 2°C, in aerobiosis
 - 7 days at 35 ± 2°C, in aerobiosis
3. Microbiological control
 - Inoculum: McFarland 0.5
 - Incubation Conditions: 24 h at 35 ± 2°C, in aerobiosis

Microorganism

Enterococcus faecalis ATCC® 51299
Enterococcus faecalis ATCC® 29212

Growth

Good
Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.