## Campylobacter CCDA Agar

Selective medium for detection and enumeration of *Campylobacter* spp, according to ISO 10272 and ISO 17995.

TYPICAL FORMULA	(g/l)
Meat Extract	10.0
Enzymatic Digest of Animal Tissues	10.0
Enzymatic Digest of Casein	3.0
Sodium Chloride	5.0
Charcoal	4.0
Sodium Deoxycholate	1.0
Iron(II) Sulfate	0.25
Sodium Pyruvate	0.25
Agar	12.0
Cefoperazone	0.032
Amphotericin B	0.01
Final pH 7.4 ± 0.2	

### DESCRIPTION

Campylobacter CCDA Agar is a medium used for the selective isolation and colony counting of *Campylobacter* spp from food, water and environmental samples.

This medium, also known as Modified Charcoal Cefoperazone Deoxycholate Agar (mCCD agar), complies with the recommendations of ISO 10272 and ISO 17995 for the microbial examination of food and water, respectively.

### PRINCIPLE

Meat extract, enzymatic digest of animal tissues and enzymatic digest of casein provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Charcoal, ferrous sulfate and sodium pyruvate act as oxygen quencher. Sodium deoxycholate suppresses the growth of Gram-positive bacteria. Cefoperazone is a cephalosporin antibiotic which is effective against Gram-negative enteric bacilli and some Gram-positive organism. Amphotericin B is inhibitory for fungi. Agar is the solidifying agent.

### TECHNIQUE

Inoculate the plates by using the enrichment culture obtained in a suitable medium such as Bolton Broth (ref. 470340) or streak the sample directly over the agar surface to achieve single colonies.

Incubate microaerobically at 37°C or 41.5 ± 1°C for 40-48 hours.

### INTERPRETATION OF RESULTS

Observe typical colonies of *Campylobacter* species which appear greyish, flat and moist, often with a metallic sheen, and with a tendency to spread. Other forms of colonies may occur.

For confirmation, subculture a suspected colony to a Columbia Blood Agar plate (ref. 11025) and examine the pure culture for morphology, motility, microaerobic growth at 25°C, aerobic growth at 41.5°C and the presence of oxidase.

### STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

**For professional use only.** Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

### REFERENCES

1. Public Health England (2018). Detection and enumeration of *Campylobacter* species National Infection Service. Food, Water & Environmental Microbiology Standard Method FNES15 (F21) Version 4.
2. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
3. ISO 17995:2019. Water quality -- Detection and enumeration of thermotolerant *Campylobacter* spp.
4. ISO 10272-1:2017. Microbiology of the food chain – Horizontal method for detection and enumeration of *Campylobacter* spp – Part 1: Detection method.
5. ISO 10272-2:2017. Microbiology of the food chain – Horizontal method for detection and enumeration of *Campylobacter* spp – Part 2: Colony count technique.
6. Bolton, F. J., D. N. Hutchinson, and G. Parker (1988) Eur. J. Clin. Microbiol. Infect Dis. 7:155-160.
7. Bolton, F. J., D. N. Hutchinson, and D. Coates (1984) J. Clin. Microbiol. 19:169-171.

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Campylobacter CCDA Agar

### PRESENTATION

Ready to use plates (90 mm) containing  $22 \pm 1$  ml of medium

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
10409	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>

### pH OF THE MEDIUM

7.4  $\pm$  0.2

### USE

Campylobacter CCDA Agar is a medium used for the selective isolation and colony counting of *Campylobacter* spp from food, water and environmental samples, according to ISO 10272 and ISO 17995

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Opaque, black

### SHELF LIFE

3 months

### QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control  
7 days at  $22 \pm 2^\circ\text{C}$ , in aerobiosis  
7 days at  $35 \pm 2^\circ\text{C}$ , in aerobiosis
3. Microbiological control  
Inoculum for productivity: 50-100 CFU  
Inoculum for selectivity:  $10^4$ - $10^6$  CFU  
Incubation Conditions: 40-48 hours at  $41.5 \pm 1^\circ\text{C}$ , in microaerobic atmosphere

#### Microorganism

*Campylobacter jejuni* WDCM 00156  
*Escherichia coli* WDCM 00013  
*Staphylococcus aureus* WDCM 00034

#### Growth

Good  
Partially or totally Inhibited  
Inhibited

### TABLE OF SYMBOLS

 LOT	Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
 REF	Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instructions for use	

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- Public Health England (2018). Detection and enumeration of *Campylobacter* species National Infection Service. Food, Water & Environmental Microbiology Standard Method FNES15 (F21) Version 4.
- EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 17995:2019. Water quality -- Detection and enumeration of thermotolerant *Campylobacter* spp.
- ISO 10272-1:2017. Microbiology of the food chain – Horizontal method for detection and enumeration of *Campylobacter* spp – Part 1: Detection method.
- ISO 10272-2:2017. Microbiology of the food chain – Horizontal method for detection and enumeration of *Campylobacter* spp – Part 2: Colony count technique.
- Bolton, F. J., D. N. Hutchinson, and G. Parker (1988) Eur. J. Clin. Microbiol. Infect Dis. 7:155-160.
- Bolton, F. J., D. N. Hutchinson, and D. Coates (1984) J. Clin. Microbiol. 19:169-171.

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**LIOFILCHEM® S.r.l.**