Legionella BCYE + AB Agar

Selective medium for detection and enumeration of *Legionella* spp, according to ISO 11731.

TYPICAL FORMULA	(g/l)
Yeast Extract	10.0
Activated Charcoal	2.0
ACES Buffer	10.0
Potassium Hydroxide	2.8
α -Ketoglutarate	1.0
L-Cystine	0.4
Ferric Pyrophosphate	0.25
Agar	12.0
Polymyxin B Sulfate	80 000 IU
Sodium Cefazolin	0.009
Pimaricin	0.07
Final pH 6.8 \pm 0.2	

DESCRIPTION

Legionella BCYE + AB Agar is a medium used for the selective isolation and numeration of *Legionella* spp from water samples and water related matrices.

This medium complies with ISO 11731:2017.

PRINCIPLE

Yeast extract provides amino acids, nitrogen, carbon, vitamins and minerals. Activated charcoal decomposes hydrogen peroxide, a metabolic product toxic to *Legionella* spp, and may also collect carbon dioxide and modify surface tension. ACES buffer (N-2-acetamido-2-aminoethane sulfonic acid) and potassium hydroxide maintain the proper pH for optimal growth. Alpha-ketoglutarate, cysteine and ferric pyrophosphate are incorporated to satisfy the specific nutritional requirements of *Legionella* species. Agar is the solidifying agent. Polymyxin B and cefazolin are selective agents against contaminating bacteria while pimaricin (natamycin) is included as antifungal agent.

TECHNIQUE

Membrane filtration is recommended as concentration technique when the number of legionellae in the water sample is unknown or not greater than 10⁴ CFU per litre. Direct plating of the unconcentrated sample is carried out if a higher concentration of *Legionella* spp is expected. Highly contaminated samples should be diluted and direct plating used before and after the pre-treatment with heat or acid solution.

Legionella BCYE + AB Agar is used in conjunction with Legionella BCYE Agar (ref. 10051), which has the same formulation without antibiotics. **Note:** Both these two media are inoculated with samples with a low concentration of contaminating organisms; samples containing a high or extremely high concentration of interfering microorganisms should be tested by using highly selective culture media such as Legionella GVPC Agar (ref. 10128) and Legionella MWY Agar (ref. 10127).

Incubate the inoculated plates at 36 \pm 2°C for 7 to 10 days in humidified atmosphere (air with 2.5% CO₂ can be beneficial for the growth of some Legionella but is not essential).

INTERPRETATION OF RESULTS

Examine the plates either on day 2, 3, 4 or 5 for monitoring overgrowth of contaminant flora. At the end of the incubation period observe growth and fluorescence under long-wave UV light.

Colonies of *Legionella* are often white-grey but can also appear in other colours. They are smooth with an entire edge and exhibit a characteristic ground-glass appearance. Under UV light colonies usually autofluoresce brilliant white.

For confirmation, subculture from the plate(s) showing the highest counts of presumptive colonies of *Legionella* per water volume onto plates of Legionella BCYE Agar and Legionella BCYE Agar w/o cysteine (ref. 10412). Regard as *Legionella* those colonies which grow on Legionella BCYE Agar but fail to grow on Legionella BCYE Agar w/o Cysteine.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. ISO 11731:2017. Water quality – Enumeration of *Legionella*.
2. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
3. Clesceri L.S., A.E. Greenberg and A.D. Eaton (1998) Standard methods for the examination of water and wastewater, 20th ed. American Public Health Association (APHA), Washington, D.C.
4. Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3): 298.

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Legionella BCYE + AB Agar

PRESENTATION

Ready to use plates (90 mm) containing 22±1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
10448	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
10448*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 x 10 plates in cardboard box

pH OF THE MEDIUM

6.8 ± 0.2

USE

Legionella BCYE + AB Agar is a medium formulated according to ISO 11731 and used for the isolation of *Legionella* species from clinical specimens and environmental samples

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Black, opaque

SHELF LIFE

6 months

QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control
7 days at 22 ± 2°C, in aerobiosis
7 days at 35 ± 2°C, in aerobiosis
3. Microbiological control
Inoculum for productivity: 50-100 CFU
Incubation Conditions: 2-5 days at 36 ± 2°C

Microorganism

Legionella pneumophila

WDCM 00107

Growth

Good

Specification

White-grey-blue-purple colonies with entire edge and ground-glass appearance

Escherichia coli

WDCM 00013

Totally or partially inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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- ISO 11731:2017. Water quality – Enumeration of *Legionella*.
- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- Clesceri L.S., A.E. Greenberg and A.D. Eaton (1998) Standard methods for the examination of water and wastewater, 20th ed. American Public Health Association (APHA), Washington, D.C.
- Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3): 298.

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LIOFILCHEM® S.r.l.