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OXACILLIN RESISTANCE STAPHYLOCOCCUS AGAR

Medium for oxacillin-resistant staphylococci determination.

TYPICAL FORMULA	(g/l)
Meat Extract	2.0
Casamino Acids Technical	19.0
Soluble Potatoes Starch	1.5
Sodium Chloride	40.0
Oxacillin (874 µg/mg)	6.5 mg
Agar	17.0
Final pH 7.4 ± 0.2	

DESCRIPTION

OXACILLIN RESISTANCE STAPHYLOCOCCUS AGAR is a medium for oxacillin-resistant staphylococci determination.

PRINCIPLE

Meat extract and casamino acids technical are a source of nitrogen, carbon, sulphur and other essential growth factors. Starch is a source of energy. Sodium chloride carries out a selective action in the medium. Oxacillin is an antibiotic belonging to isoazolipenicillins group, semisynthetic penicillins particularly indicated in the therapy of infections caused by penicillase-producing staphylococci.

TECHNIQUE

1. Prepare the inoculum by suspending several well-isolated colonies of *S. aureus* to test, isolated from a 18-24 hours plate culture, into a tube of suitable broth medium (for example Tryptic Soy Broth).
2. Inoculate 10 µl of above-mentioned suspension using a micropipette.
3. Alternatively saturate a cotton swab with the test suspension, then gently press out excess fluid against the inner wall of the tube and so streak it on the plate.
4. Inoculate also a plate of TSA with 5% of sheep blood as non-selective growth control.
5. Incubate at 30-35°C for 24 hours. DO NOT EXCEED OVER 35°C.

INTERPRETATION OF RESULTS

After incubation observe for the formation of colonies. Any growth, even one colony, indicates that the isolated colonies are oxacillin-resistant. No growth indicates that the microorganism is susceptible to PRPs (methicillin, nafcillin, oxacillin). The microorganism which grows on OXACILLIN RESISTANCE STAPHYLOCOCCUS AGAR should be reported as resistant all β-lactamic antimicrobial agents.

STORAGE AND TRANSPORT CONDITIONS

6-12°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Chain, Florey and Jennings. 1949. In Florey, Chain, Heatley, Jennings, Sanders and Abraham (ed.), Antibiotics, vol. II. Oxford University Press, London.
2. Barrett, McGehee and Finland. 1968. N. Engl. J. Med. 279: 444.

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PRODUCT SPECIFICATIONS

NAME

OXACILLIN RESISTANCE STAPHYLOCOCCUS AGAR

PRESENTATION

Ready to use plates (90 mm) containing 22+/-1 ml of medium

STORAGE

6-12°C

PACKAGING

Ref.	Content	Packaging
10600	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box

pH OF THE MEDIUM

7.4 ± 0.2

USE

OXACILLIN RESISTANCE STAPHYLOCOCCUS AGAR is a medium for oxacillin-resistant staphylococci determination

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Opaque, light amber medium

SHELF LIFE

4 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 CFU/ml
Inoculum for selectivity: 10⁴-10⁵ CFU/ml
Inoculum for specificity: ≤10⁴ CFU/ml
Incubation Conditions: 18-24 hours at 34 ± 1°C, in aerobiosis

Microorganism		Growth
<i>Escherichia coli</i>	ATCC 25922	Inhibited
<i>Staphylococcus aureus</i>	ATCC 33592	Good
<i>Staphylococcus aureus</i>	ATCC 25923	Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 REF	Catalogue number	 Temperature limitation		 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net