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## HELICOBACTER PYLORI EGG YOLK EMULSION AGAR

Medium for the cultivation of *Helicobacter pylori*.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	10.0
Meat Peptic Digest	5.0
Heart Pancreatic Digest	3.0
Yeast Extract	5.0
Corn Starch	1.0
Sodium Chloride	5.0
Agar	14.0
Vitalex	10.0 ml
Egg Yolk Emulsion	100.0 ml
Triphenyltetrazolium Chloride (TTC)	0.04
Final pH 7.3 ± 0.2	

### DESCRIPTION

HELICOBACTER PYLORI EGG YOLK EMULSION AGAR is a medium containing Columbia Agar as base and used for the cultivation of *Helicobacter pylori*.

### PRINCIPLE

Peptones provide nitrogen, carbon, amino acids and vitamins. Yeast extract and corn starch are included to supply a carbon energy source in the form of B-complex vitamins. Sodium Chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Supplementation with vitalex and egg yolk emulsion provides additional growth factors for fastidious microorganisms. TTC is a colorless compound that is reduced by *H. pylori* to formazan, a red compound

### TECHNIQUE

Inoculate the plate by streaking the microorganism onto the agar surface. Incubate at 35 ± 2°C for 3-5 days in microaerobic atmosphere.

Note: *H. pylori* is extremely sensitive to desiccation and exposure to oxygen. The recovery of the organism from gastric biopsies is improved by direct cultivation, immediately after collection, onto a selective medium such as Helicobacter pylori Agar (ref. 10082). If a delay is expected, transport media such as Gesa Medium (ref. 27001) is recommended.

### INTERPRETATION OF RESULTS

*Helicobacter pylori* cultivates with translucent and non-coalescent colonies. The red color induced by the TTC makes the colonies easy to spot on this medium, even early in the incubation. Contamination can easily be identified by the lack of the typical golden sheen of *H. pylori*. Also the microorganism grows at 35°C and does not grow at 25°C, it is urease, catalase and oxidase positive and hippurate negative. Final identification should be done performing appropriate biochemical tests.

### STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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## PRODUCT SPECIFICATIONS

### NAME

HELICOBACTER PYLORI EGG YOLK EMULSION AGAR

### PRESENTATION

Ready-to-use plates (90 mm) containing 22 ± 1 ml of medium

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
10605	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>

### pH OF THE MEDIUM

7.3 ± 0.2

### USE

HELICOBACTER PYLORI EGG YOLK EMULSION AGAR is a medium containing Columbia Agar as base and used for the cultivation of *Helicobacter pylori*

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Light amber, slightly opalescent medium

### SHELF LIFE

120 days

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
7 days at 22 ± 1°C, in aerobiosis  
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control  
Inoculum for productivity: 10-100 CFU/ml  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> CFU/ml  
Incubation Conditions: 24-48 hours at 35 ± 2°C, in microaerobic atmosphere

Microorganism		Growth	Specification
<i>Helicobacter pylori</i>	ATCC® 43504	Good	Red colonies with golden sheen
<i>Escherichia coli</i>	ATCC® 25922	Good	Red colonies

### TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

**LIOFILCHEM® S.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY  
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

## HELICOBACTER PYLORI EGG YOLK EMULSION AGAR

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