## Columbia Agar (Sheep blood 5%)

Medium for the cultivation of fastidious and non-fastidious microorganisms from clinical specimens.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	12.0
Starch	1.0
Peptic Digest of Animal Tissue	5.0
Sodium Chloride	5.0
Yeast Extract	3.0
Beef Extract	3.0
Sheep Blood Defibrinated	50.0 ml
Agar	14.0
Final pH 7.3 ± 0.2	

### DESCRIPTION

Columbia Agar (Sheep blood 5%) is a medium used for the growth of fastidious and non-fastidious organisms and determination of haemolytic reactions

### PRINCIPLE

Pancreatic digest of casein, beef extract and peptic digest of animal tissue and yeast extract provide nitrogen, carbon, sulfur and other essential growth factors. Sodium chloride maintains the osmotic balance of the medium. Sheep blood defibrinated supplies additional growth factors for fastidious microorganisms and allows to evidence the haemolytic reactions. Agar is the solidifying agent.

### TECHNIQUE

Inoculate the plates by streaking the sample to examine onto the surface of the medium using a sterile loop in order to isolate single colonies. Then stab the agar several times to deposit beta-haemolytic streptococci below the agar surface. Incubate plates at 36 ± 1°C for 18-48 hours, in aerobic, anaerobic or microaerophilic atmosphere.

### INTERPRETATION OF RESULTS

Observe for growth and haemolytic reactions. Four different kinds of haemolysis can be distinguished:

1. alfa haemolysis: haemoglobin is reduced to metahaemoglobin in the medium that surrounds the colony with a consequent greenish decolouring of the medium;
2. beta haemolysis: it is the lysis of erythrocytes which is evident in a light zone around the colony;
3. gamma haemolysis: it does not occur any destruction of erythrocytes or any change in the medium;
4. alfa-prime haemolysis: it's evident a small zone of complete haemolysis surrounded by an area of partial lysis

### STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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Columbia Agar (Sheep blood 5%)

### PRESENTATION

Ready to use plates (90 mm) containing 22 ± 1 ml of medium

### STORAGE

2-8°C

### PACKAGING

Ref.	Content	Packaging
11025	20 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 2 x 10 plates in cardboard box</li> </ul>
11025*	100 plates	<ul style="list-style-type: none"> <li>• 10 plates in thermally soldered film</li> <li>• 10 piles (10 x 10 plates) in cardboard box</li> </ul>

### pH OF THE MEDIUM

7.3 ± 0.2

### USE

Columbia Agar (Sheep blood 5%) is a medium used for the cultivation of fastidious microorganisms from clinical specimens

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

Cherry red, opaque

### SHELF LIFE

90 days

### QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control  
7 days at 22 ± 2°C, in aerobiosis  
7 days at 35 ± 2°C, in aerobiosis
3. Microbiological control  
Inoculum for productivity: 50-100 CFU  
Incubation Conditions: 18-24 hours at 36 ± 1°C

Microorganism		Growth	Haemolysis
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	β-haemolysis
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Good	α-haemolysis
<i>Streptococcus pyogenes</i>	ATCC® 19615	Good	β-haemolysis
<i>Escherichia coli</i>	ATCC® 25922	Good	None

### TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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**LIOFILCHEM® S.r.l.**