Thayer Martin Agar

Selective medium for detection of *Neisseria* spp and other fastidious organisms from clinical specimens.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Casein	7.5
Enzymatic Digest of Animal Tissues	7.5
Corn Starch	1.0
Potassium Phospahte Bibasic	4.0
Potassium Phospahte Monobasic	1.0
Sodium Chloride	5.0
Agar	17.0
Defibrinated Horse Blood	50.0 ml
Cyanocobalamine	0.05 mg
Glutamine	5.5 mg
Thiamine HCl	1.0 ml
L-Cysteine HCl	0.12 ml
Adenine	5.0 mg
p-Aminobenzoic Acid	0.065 mg
Coccarboxylase	0.5 mg
Guanine	0.15 mg
NAD	1.25 mg
Ferric Nitrate	0.1 mg
Vancomycin	1.0 mg
Colistin Sulphate	3.75 mg
Amphotericin B	0.5 mg
Trimethoprim	1.5 mg
Final pH 7.2 ± 0.2	

DESCRIPTION

Thayer Martin Agar is a selective medium used for the isolation and cultivation of pathogenic *Neisseria* spp, from clinical specimens containing mixed flora of bacteria and fungi.

PRINCIPLE

Enzymatic digest of casein and enzymatic digest of animal tissues provide amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Corn starch absorbs any toxic metabolites produced. The phosphates are buffering agents. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Blood, cofactors, vitamins and NAD are included to enhance recovery of *Heemophilus* and *Neisseria* spp. Antimicrobial supplements inhibit contaminating organisms improving selectivity.

TECHNIQUE

Inoculate by rolling a culture swab across a segment of the plate or in a large "Z" pattern so that an adequate area of the plate is inoculated. Streaking of the plate is carried out with a sterile loop to ensure an appropriate dispersion of the organisms. Incubate at 35 ± 2°C for 18-48 hours, in a sealed jar with at least 70% humidity and 5-10% carbon dioxide.

INTERPRETATION OF RESULTS

Examine the plates after 24 hours incubation and, if negative, reincubate for further 24 hours. Presumptive gonococcus colonies are identified by the Gram stain, oxidase and sugar fermentation reactions.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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Thayer Martin Agar

PRESENTATION

Ready to use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
11040	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
11040*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 piles (10 x 10 plates) in cardboard box

pH OF THE MEDIUM

7.2 ± 0.2

USE

Thayer Martin Agar is a selective medium used for the isolation and cultivation of pathogenic *Neisseria* spp, from clinical specimens containing mixed flora of bacteria and fungi

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Opaque, chocolate brown

SHELF LIFE



4 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at $22 \pm 2^\circ\text{C}$, in aerobiosis
7 days at $35 \pm 2^\circ\text{C}$, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10^4 - 10^6 CFU
Incubation Conditions: 18-24 hours at $35 \pm 2^\circ\text{C}$, under 5-10% CO_2 enriched atmosphere.

Microorganism		Growth
<i>Neisseria gonorrhoeae</i>	ATCC® 43069	Good
<i>Neisseria meningitidis</i>	ATCC® 13090	Good
<i>Escherichia coli</i>	ATCC® 25922	Inhibited
<i>Candida albicans</i>	ATCC® 10231	Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.