Mycoplasma Agar

Selective medium for the detection of *Mycoplasma* spp from clinical specimens and nonclinical samples.

TYPICAL FORMULA	(g/l)
Tryptone	15.7
Soy Peptone	2.7
Glucose	2.5
Sodium Chloride	4.5
Potassium Phosphate Bibasic	1.8
Agar	11.8
Citric Acid	0.2
Manganous Sulphate Monohydrate	0.2
Penicillin G Potassium Salt	0.6
Thallium Acetate	0.05
Yeast Extract 25%	20 ml
L-Cysteine 2%	5.2 ml
Vitalax	10.0 ml
Urea 40%	2.5 ml
Horse Serum	20.0 ml
Final pH 6.4 ± 0.2	

DESCRIPTION

Mycoplasma Agar is a medium used for the selective isolation of urogenital mycoplasmas.

This medium can also be used for examination of samples from animals or the environment, such as for detection of avian mycoplasmas like *M. gallisepticum*, *M. synoviae*, *M. iowe* and *M. meleagridis*, which are of major concern to the poultry industry.

PRINCIPLE

Bacteriological research of mycoplasma, organisms which are devoid of cell wall, has to be carried out on culture media which allow to satisfy their metabolic requirements. Tryptone, soy peptone and yeast extract provide nitrogen, vitamins, and minerals which favour bacterial growth, whilst glucose is a source of carbohydrates. Sodium chloride and potassium phosphate maintain osmotic balance of the medium. Agar is the solidifying agent. Citric acid is the forerunner of important metabolic enzymes and is an anticoagulant agent. Manganous sulphate monohydrate is used as an indicator and allows the identification of *Ureaplasma urealyticum*. The presence of penicillin prevents the development of bacterial contaminants. Thallium acetate is the other selective agent incorporated in the medium to enhance the recovery of mycoplasmas over faster-growing organisms. Vitalax is a growth supplement which supports the growth of fastidious microorganisms. Urea 40% is a supplement for the determination of urease activity of microorganisms: urea provides a source of nitrogen and supports a vigorous growth for microorganisms that produce urease. Horse serum provides essential lipoproteins.

TECHNIQUE

Inoculate the medium by streaking the sample on the agar surface immediately after the collection. Incubate at 35 ± 2°C for up to 10 days in an anaerobic or microaerophilic atmosphere.

INTERPRETATION OF RESULTS

Observe colonies morphology by direct transmitted illumination under low magnification.

Mycoplasma hominis is recognised by typical "fried-egg" colonies or finely granular colonies with a berry-like appearance that penetrate the agar surface. *Ureaplasma* show typical sea-urchin-like colonies.

Refer to specific guidances for more details.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for In vitro diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Shepard, M.C., Lunceford, C.D. (1976) J. Clin. Microbiol. 3:613.
2. Adler HE, Fabricant J, Yamamoto R, Berg J.(1958) Isolation and identification of pleuropneumonia-like organisms of avian origin. In: Symposium on chronic respiratory diseases of poultry. American Journal of Veterinary Research; 19:440-7.

PRODUCT SPECIFICATIONS

NAME

Mycoplasma Agar

PRESENTATION

Ready-to-use plates (90 mm) containing 22±1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
11205	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box

pH OF THE MEDIUM

6.4 ± 0.2

USE

Mycoplasma Agar is a selective medium for the detection of *Mycoplasma* spp

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Light amber, very slightly opalescent

SHELF LIFE


4 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 2°C, in aerobiosis
7 days at 35 ± 2°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Inoculum for specificity: 10³-10⁴ CFU
Incubation conditions: 18-24 h at 35 ± 2°C, in anaerobic or microaerophilic atmosphere

Microorganism		Growth
<i>Candida albicans</i>	ATCC® 10231	good
<i>Ureaplasma urealyticum</i>	ATCC® 27618	good
<i>Mycoplasma hominis</i>	ATCC® 23114	good
<i>Proteus mirabilis</i>	ATCC® 25933	inhibited
<i>Staphylococcus aureus</i>	ATCC® 25923	inhibited
<i>Escherichia coli</i>	ATCC® 25922	inhibited

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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LIOFILCHEM® S.r.l.