MUELLER HINTON II AGAR (Sheep Blood 5%)

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MUELLER HINTON II AGAR (Sheep Blood 5%)

Medium recommended for antimicrobial disc diffusion susceptibility testing of *Streptococcus pneumoniae*.

TYPICAL FORMULA	(g/l)
Beef Extract	2.0
Acid Digest of Casein	17.5
Starch	1.5
Agar	17.0
Sheep Defibrinated Blood	50 ml
Final pH 7.3 ± 0.2	

DESCRIPTION

MUELLER HINTON II AGAR (Sheep Blood 5%) is recommended for antimicrobial disc diffusion susceptibility testing of *Streptococcus pneumoniae*.

PRINCIPLE

Beef extract and acid digest of casein supply amino acids and other nitrogenous substances, minerals, vitamins, carbon and other nutrients to support the growth of microorganisms. Starch acts as a protective colloid against toxic substances that may be present in the medium. Beside hydrolysis of starch during autoclaving provides a small amount of glucose which is a source of energy. Agar is the solidifying agent. Sheep blood is a further growth supplementation for fastidious microorganisms.

TECHNIQUE

Transfer 4-5 colonies into a liquid culture medium and incubate at 36±1°C until obtaining a 0.5 Mac Farland opacity. Introduce a sterile swab into the broth and inoculate a plate of MUELLER HINTON II AGAR (Sheep Blood 5%) by streaking the swab 2 or 3 times onto the entire surface. Press the disc containing the antimicrobial agent on the agar surface. Incubate at 36±1°C for 18-24 hours. Measure the diameter of the inhibition halo using a pair of compasses.

INTERPRETATION OF RESULTS

Zone diameter measured around discs should be compared with those in the NCCLS Document M100 (M2). Results obtained with specific organisms may then be reported as resistant, intermediate or susceptible. It is important to measure the diameter of the zone of growth inhibition, not of hemolysis inhibition.

STORAGE AND TRANSPORT CONDITIONS

6-12°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Bauer et al. (1966). J. Clin. Pathol. 45:493-496.
2. Mueller, J.H., and Hinton. 1941. Proc. Soc. Exp. Biol. Med. 48: 330-333.
3. CLSI. Performance standards for susceptibility testing; Twelve Informational Supplement. CLSI Document M100-S22, January 2012.
4. Barry, Garcia and Thrupp. (1970). Am. J. Clin. Pathol. 53: 149.

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

MUELLER HINTON II AGAR (Sheep Blood 5%)

PRESENTATION

Ready to use plates

STORAGE

6-12°C

PACKAGING

Ref.	Content	Packaging
10131	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
10131*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 piles (10 x 10 plates) in cardboard box
11231	10 plates (140 mm)	<ul style="list-style-type: none"> • 5 plates in thermally soldered film • 2 x 5 plates in cardboard box
12032	5 plates (120 x 120 mm)	<ul style="list-style-type: none"> • 5 plates in thermally soldered film, in cardboard box

pH OF THE MEDIUM

7.3 ± 0.2

USE

MUELLER HINTON II AGAR (Sheep Blood 5%) is recommended for antimicrobial disc diffusion susceptibility testing of *Streptococcus pneumoniae*

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Cherry red medium, opaque

SHELF LIFE

75 days

QUALITY CONTROL

1. Control of general characteristics, label and print
2. Sterility control
7 days at 22 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
3. Microbiological control
Incubation Conditions: 18-24 hours at 36 ± 1°C

Microorganism		Growth
<i>Streptococcus pyogenes</i>	ATCC 19615	Good
<i>Enterococcus faecalis</i>	ATCC 19433	Good
<i>Escherichia coli</i>	ATCC 25922	Good
<i>Staphylococcus aureus</i>	ATCC 25923	Good

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net