Tinsdale Agar

Selective medium for the isolation and differentiation of *Corynebacterium diphtheriae*.

TYPICAL FORMULA	(g/l)
Peptone	10.0
Tryptone	10.0
Sodium Chloride	5.0
Sodium Thiosulphate Anhydrous	0.5
L-Cystine	0.24
Agar	15.0
Horse Serum	150 ml
Potassium Tellurite	0.3
Final pH 8.0 ± 0.2	

DESCRIPTION

Tinsdale Agar is a medium used for the selective isolation and differentiation of *Corynebacterium diphtheriae* from clinical specimens.

PRINCIPLE

Peptone and tryptone provide nitrogen, vitamins, carbon and amino acids. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Horse serum provides essential growth factors. Sodium thiosulphate and L-cystine are H₂S indicators. Potassium tellurite is a selective agent. The formation of black to brown halos surrounding the colony results from the reduction of potassium tellurite to metallic tellurite.

TECHNIQUE

Inoculate the medium to obtain well separated colonies. Stab deep into the agar in several parts of the plate to stimulate browning. Incubate at 35 ± 2°C for 24-48 hours.

INTERPRETATION OF RESULTS

Brown-black colonies surrounded by a brown-black halo presumptively indicates the presence of *C. diphtheriae*.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Tinsdale (1947). J. Pathol. Bacteriol. 59: 461.
2. Billings (1956) An investigation of Tinsdale Tellurite Medium: its usefulness and mechanisms of halo-formation. M.S. Thesis. University of Michigan, Ann Arbor, Mich.

PRODUCT SPECIFICATIONS

NAME

Tinsdale Agar

PRESENTATION

Ready to use plates (90 mm) containing 22 ± 1 ml of medium

STORAGE

2-8°C

PACKAGING

Ref.	Content	Packaging
11250	20 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 2 x 10 plates in cardboard box
11250*	100 plates	<ul style="list-style-type: none"> • 10 plates in thermally soldered film • 10 x 10 plates in cardboard box

pH OF THE MEDIUM

8.0 ± 0.2

USE

Tinsdale Agar is a medium used for the selective isolation and differentiation of *Corynebacterium diphtheriae* from clinical specimens

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Light beige medium, opalescent

SHELF LIFE

4 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 2°C, in aerobiosis
7 days at 35 ± 2°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Inoculum for specificity: 10³-10⁴ CFU
Incubation Conditions: 18-48 hours at 35 ± 2°C, in aerobiosis

Microorganism

Corynebacterium diphtheriae ATCC® 27010
Klebsiella pneumoniae ATCC® 13883
Streptococcus pyogenes ATCC® 19615

Growth

Good
Markedly inhibited
Poor

Colony Appearance

Brown with halos

Brown to black without halos

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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LIOFILCHEM® S.r.l.