R.P.M.I. AGAR

Medium for antifungal susceptibility testing.

TYPICAL FORMULA	(g/l)
RPMI1640 + MOPS 0.165M + L-glutamine	46.19
Glucose	20.0
Agar	15.0
Final pH 7.0 ± 0.2	

DESCRIPTION

R.P.M.I. AGAR is a medium for antifungal susceptibility testing and is recommended for tests performed using the Antifungal MIC TEST STRIP.

TECHNIQUE

Pick 4-5 colonies from 24-hours old cultures and suspend in saline solution or sterile water in order to obtain a 0.5 McFarland cell density. Introduce a sterile swab into the culture and inoculate R.P.M.I. AGAR by streaking the swab two or three times onto the entire surface. Allow the plates to dry, then press the discs or strips containing the antifungal onto the agar surface. Incubate at 35 °C for 24-48 hours.

INTERPRETATION OF RESULTS

Measure the inhibition zone around the discs with a compass. In case of usage of MIC TEST STRIP, read the M.I.C. value where the edge of the inhibition ellipse intersects the strip (intersection between two scale segments should be round up to the higher value).

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING and PRECAUTIONS

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of $\geq 1\%$.

The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

1. Szekely A. et al. Comparison of Etest and broth microdilution for antifungal drug susceptibility testing of moulds. JCM, vol. 37, no. 5, p. 1480-1483, 1999.
2. Johnson E.M. et al. Lack of correlation of in vitro amphotericin B susceptibility testing with outcome in a murine model of aspergillus infection. JAC; vol. 45, no. 1, p. 85-93, 2000.
3. Pfaller M.A. et al. In vitro susceptibility testing of filamentous fungi: comparison of Etest and reference microdilution methods for determining itraconazole MICs. JCM, vol. 38, no. 9, p. 3359-3361, 2000.
4. CLSI M38-A2: Reference Method for Broth Dilution Antifungal Susceptibility Testing of Conidium-Forming Filamentous Fungi; Proposed Standard (1998).
5. CLSI M27-A3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard - Third Edition.
6. CLSI M27-S3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Third Informational supplement.

PRODUCT SPECIFICATIONS**NAME**
R.P.M.I. AGAR**PRESENTATION**

Plates containing ready to use solid medium.

STORAGE

10-25°C

PACKAGING

Code	Content	Packaging
11509	20 plates	10 plates in thermally soldered film 2 x 10 plates in cardboard box
10233	10 plates (140mm)	10 plates in thermally soldered film

pH OF THE MEDIUM

7.0 ± 0.2

USE

R.P.M.I. AGAR is a medium for antifungal susceptibility testing and is recommended for tests performed using the Antifungal MIC TEST STRIP.

TECHNIQUE

Refer to technical sheet of the product.

APPEARANCE OF THE MEDIUM

Clear medium, light rose in colour.

SHELF LIFE

6 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 20 ± 2°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Incubation conditions: 18-24 h at 36 ± 1°C

Microorganism		Growth
<i>Candida albicans</i>	ATCC 90028	Good
<i>Candida parapsilosis</i>	ATCC 22019	Good
<i>Candida krusei</i>	ATCC 6258	Good
<i>Saccharomyces cerevisiae</i>	ATCC 9763	Good

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	LOT Batch code	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	

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BIBLIOGRAPHIE

1. Szekely A. et al. Comparison of Etest and broth microdilution for antifungal drug susceptibility testing of moulds. JCM, vol. 37, no. 5, p. 1480-1483, 1999.
2. Johnson E.M. et al. Lack of correlation of *in vitro* amphotericin B susceptibility testing with outcome in a murine model of aspergillus infection. JAC; vol. 45, no. 1, p. 85-93, 2000.
3. Pfaller M.A. et al. *In vitro* susceptibility testing of filamentous fungi: comparison of Etest and reference microdilution methods for determining itraconazole MICs. JCM, vol. 38, no. 9, p. 3359-3361, 2000.
4. CLSI M38-A2: Reference Method for Broth Dilution Antifungal Susceptibility Testing of Conidium-Forming Filamentous Fungi; Proposed Standard (1998).
5. CLSI M27-A3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard - Third Edition.
6. CLSI M27-S3. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Third Informational supplement.

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