Violet Red Bile Glucose Agar + Neutralizing

Selective medium for detection and enumeration of Enterobacteriaceae with inactivation of disinfectants.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Animal Tissues	7.0
Yeast Extract	3.0
Glucose	10.0
Sodium Chloride	5.0
Bile Salts	1.5
Neutral Red	0.03
Crystal Violet	0.002
Agar	14.0
Histidine	1.0
Lecithin	0.7
Polysorbate 80	5.0
Sodium Thiosulfate	0.5
Final pH 7.4 ± 0.2	

DESCRIPTION

Violet Red Bile Glucose Agar + Neutralizing is a selective culture medium in contact plates used for the detection of bile-tolerant Gram-negative bacteria in environmental and personnel hygiene monitoring. Neutralizing agents are included in the medium to inactivate residual disinfectants allowing also comparative results before and after cleaning.

The formulation of the medium (not considering the neutralizing components) complies with the recommendations of the harmonized USP/EP/JP method for the microbiological examination of non-sterile pharmaceutical products and with ISO 21528 for analysis of food and environmental samples in the area of food production and handling.

PRINCIPLE

Enzymatic digest of animal tissues provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Glucose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium. Bile salts and Crystal violet are selective agents effective against Gram-positive cocci. Neutral red is the pH indicator. Agar is the solidifying agent. Lecithin neutralizes quaternary ammonium compounds. Polysorbate 80 (Tween 80) is effective against phenolic compounds and mercurial derivatives. Sodium thiosulfate neutralizes halogen compounds.

TECHNIQUE

For active air monitoring, insert the plate without the lid in an air sampler and draw a volume of air from 100 to 1000 liters.

For surfaces and personnel hygiene monitoring, such as for sampling of clothing and face masks, firmly press the agar medium against the test area for about 10 sec.

Incubate the plates aerobically at 30-35°C for 18-24 hours or 37°C for 24 ± 2 hours, depending on the method used

INTERPRETATION OF RESULTS

Count characteristic pink to red colonies (with or without precipitation haloes).

Confirm by subculturing onto a non selective agar medium looking for oxidase reaction (ref. 88029) and glucose fermentation (ref. 88202). Colonies that are oxidase-negative and glucose-positive are confirmed as Enterobacteriaceae.

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used by properly trained operators only.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 1: Detection of *Enterobacteriaceae*.
2. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique.
3. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
5. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
6. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 www.liofilchem.com

PRODUCT SPECIFICATIONS

NAME

Violet Red Bile Glucose Agar + Neutralizing

PRESENTATION

55 mm contact plate containing 17 ± 1 ml of medium

STORAGE

10-25°C

PACKAGING

Ref.	Content	Packaging
15382	20 plates	<ul style="list-style-type: none"> individually packed in transparent blister of 2 pieces double-wrapped

pH OF THE MEDIUM

7.4 ± 0.2

USE

Violet Red Bile Glucose Agar + Neutralizing is a selective medium in contact plates used for the detection and enumeration of bile-tolerant Gram-negative bacteria in food and other sample materials with inactivation of residual disinfectants

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Slightly opalescent, reddish-purple

SHELF LIFE

9 months

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22.5 ± 2.5°C, in aerobiosis
7 days at 32.5 ± 2.5°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Incubation Conditions: 18-24 h at 30-35°C for *E. coli* and *P. aeruginosa* (Pharmacopoeia growth promotion);
24 ± 2 h at 37 ± 1°C for *E. coli*, *S. Typhimurium*, *S. Enteritidis* and *E. faecalis*.

Microorganism		Growth
<i>Escherichia coli</i>	WDCM 00012	Good
<i>Salmonella</i> Typhimurium	WDCM 00031	Good
<i>Salmonella</i> Enteritidis	WDCM 00030	Good
<i>Enterococcus faecalis</i>	WDCM 00009	Inhibited
<i>Pseudomonas aeruginosa</i>	ATCC® 9027	Good

TABLE OF SYMBOLS

 Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	

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- ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 1: Detection of *Enterobacteriaceae*.
- ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique.
- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
- United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
- Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.

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1. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 1: Detection of Enterobacteriaceae.
2. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique.
3. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. European Pharmacopoeia 6.5 (2009) 2.6.13. Microbiological examination of non-sterile products: Test for specified microorganisms.
5. United States Pharmacopoeia 32 NF 27 (2009) <62> Microbiological examination of non-sterile products: Test for specified microorganisms.
6. Japanese Pharmacopoeia 4.05 (2008) Microbiological examination of non-sterile products: Test for specified microorganisms.

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