Mycoplasma Selective Broth

Liquid medium for detection of urogenital mycoplasmas.

DESCRIPTION

Mycoplasma Selective Broth is a liquid medium used for the detection of *Ureaplasma* spp. (*Ureaplasma urealyticum* and *Ureaplasma parvum*) and *Mycoplasma hominis*.

PRINCIPLE OF THE METHOD

The medium is formulated to allow growth of mycoplasmas most commonly found in the urogenital tract while inhibiting contaminating microbial flora. *Ureaplasma* spp. and *M. hominis* metabolise urea and arginine, respectively, resulting in a colour change of the medium from yellow to red.

COMPOSITION (per litre of purified water)

Meat Peptone	8.0 g	Urea	1.0 g
Casein Peptone	8.0 g	Phenol Red	0.05 g
Yeast Extract	4.0 g	Vitox Supplement	10 ml
Sodium Chloride	3.5 g	Horse Serum	100 ml
Arginine Hydrochloride	5.0 g	Antibiotic Mixture	10 ml
Cystine Hydrochloride	0.1 g	Final pH 6.3 ± 0.2	

SPECIMENS

This product can be used with the following specimens:

- Vaginal and urethral swabs
- Seminal fluid

Specimen/swab should be placed in a tightly sealed transport container with a transport medium to prevent drying (e.g. Mycoplasma Transport Broth, ref. 20158). Samples should be taken directly to the laboratory and processed as soon as possible. The specimen should be collected before administering any antibiotic treatment.

PROCEDURE FOR USE

1. Allow the vial of Mycoplasma Selective Broth to come to room temperature.
2. Inoculate by placing the swab into the vial, or by adding 100 µl of a transport medium containing swab, or 200 µl of the liquid specimen.
3. Incubate at 36 ± 2°C for 24-48 hours.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Watch for the colour change in the broth after 24 hours and 48 hour of incubation by comparison with an inoculated vial.

Yellow (**negative**): absence of *Ureaplasma* spp. and/or *M. hominis* Red (**positive**): presence of *Ureaplasma* spp. and/or *M. hominis*

QUALITY CONTROL

The following microbial strain can be used: *Ureaplasma parvum* (*U. urealyticum* serovar 1) ATCC® 27813/7.

FACTORS THAT MAY INVALIDATE THE RESULTS

Condition, timing and volume of specimen collected for culture are significant variable in obtaining reliable culture results.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use only and must be used by properly trained operators.

STORAGE

2-8°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

REFERENCES

- Pereyre S, Sirand-Pugnet P, Beven L, Charron A, Renaudin H, et al. (2009) PLoS Genet 5(10): e1000677. doi:10.1371/journal.pgen.1000677.
- Lennette E.H, A. Balows, W.J. Hausler, H.J. Shadomy (1985) Manual of Clinical Microbiology 4th Ed.
- Shepard, M.C., Lunceford, C.D. (1976) J. Clin. Microbiol. 3:613.

PRESENTATION

Product	Ref.	Packaging
Mycoplasma Selective Broth	20162	20 x 2 ml vials

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Do not reuse

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LIOFILCHEM® S.r.l.