MALONATE BROTH

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MALONATE BROTH

Liquid medium for the differentiation of enterobacteria

TYPICAL FORMULA	(g/l)
Bromothymol Blue	0.024
Sodium Malonate	1.4
Sodium Chloride	1.9
Potassium Phosphate Monobasic	1.3
Potassium Phosphate Dibasic	1.2
Ammonium Sulphate	1.9
Final pH 6.7 ± 0.2	

DESCRIPTION

MALONATE BROTH is a liquid medium for the differentiation of enterobacteria.

PRINCIPLE

MALONATE BROTH contains ammonium sulphate as sole nitrogen source and sodium malonate as sole carbon source. The potassium phosphate is the buffer system of the medium. The use of malonate by enterobacteria causes alkalization of the medium and then the pH indicator bromothymol blue turns from green to blue.

PREPARATION

Verify that the content of the tube is homogeneous and clear, possibly repeatedly inverting the tube. MALONATE BROTH can be used as is or split, under aseptic conditions in smaller aliquots.

TECHNIQUE

1. Bring to room temperature or, preferably at $36\pm 1^{\circ}\text{C}$, the tubes of MALONATE BROTH needed.
2. Open the tubes only immediately before the inoculum (tubes open, closed and not used within a few hours should be discarded).
3. Inoculate introducing the sample into the tube.
4. Immediately close the tube.
5. Incubate at $36\pm 1^{\circ}\text{C}$ for 18-24 hours.

INTERPRETATION OF RESULTS

The differentiation of enterobacteria is based on the ability of these microorganisms to utilize malonate with resultant alkalinity and blue color in the medium.

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Association of Official Analytical Chemists (1995) Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD.
2. Vanderzant, C., and D.F. Splittstoesser (ed.) (1992) Compendium of methods for the microbiological examination of food, 3rd ed. American Public Health Association, Washington D.C.

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MALONATE BROTH

PRESENTATION

Glass tubes containing 10 ml of medium

STORAGE

10-25°C

PACKAGING

Ref.	Content	Packaging
24137	20 tubes x 10 ml	20 tubes in cardboard box

pH OF THE MEDIUM

6.7 ± 0.2

USE

MALONATE BROTH is a liquid medium for the differentiation of enterobacteria

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Clear green medium without precipitates

SHELF LIFE

1 year

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 10-100 UFC/ml
Incubation Conditions: 18-24 h at 35 ± 2°C, in aerobiosis

Microorganism		Growth	Features
<i>Escherichia coli</i>	ATCC 25922	Good	Green
<i>Enterobacter cloacae</i>	ATCC 13047	Good	Blue
<i>Klebsiella pneumoniae</i>	ATCC 13883	Good	Blue
<i>Salmonella typhimurium</i>	ATCC 14028	Good	Green

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net