SIM MEDIUM

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SIM MEDIUM

Differential medium based on hydrogen sulphide production, indole production and motility

TYPICAL FORMULA	(g/l)
Tryptone	8.0
Peptone	5.0
Meat Extract	4.0
Yeast Extract	4.0
Sodium Thiosulphate	0.2
Ferric Ammonium Sulphate	0.2
Agar	3.0
Final pH 7.3 ± 0.2	

DESCRIPTION

SIM MEDIUM is a medium used for the differentiation of microorganisms on the basis of hydrogen sulphide production, indole production and motility.

PRINCIPLE

Tryptone contains tryptophan, which is converted to indole. Indole is detected after incubation by addition of Kovac's Reagent (ref. 80271). Peptone and meat extract are sources of amino acids and proteins. Yeast extract supplies amino acids and group B vitamins. Sodium thiosulphate and ferric ammonium sulphate are substrates used to detect hydrogen sulphide producing black precipitates. Agar is the solidifying agent. SIM MEDIUM is a semi-solid, due to the low concentration of agar. The semi-solid nature of this medium allows for easy visual determination of motility which appears as growth extending outward from the original line of incubation.

TECHNIQUE

1. Using a wire, inoculate test organism two-thirds into the medium with stab motion.
2. Incubate with loose caps at 35 ± 2°C for 18-24 hours.
3. Examine tubes to check for motility and H₂S production.
4. Add 3-4 drops of Kovac's Reagent (ref. 80271) to each tube. Record as indole positive if a pink or red color appear, or as indole negative if there is no color change.

Note: add Kovac's Reagent after determining motility and hydrogen sulphide production.

INTERPRETATION OF RESULTS

Motility is indicated by turbidity of the medium or growth extending from inoculating stab line. H₂S production is shown by a blackening along the stab line. Indole production is seen as the production of a red color after the addition of Kovac's Reagents (ref. 80271).

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

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SIM MEDIUM

PRESENTATION

Slant medium containing 12 ± 1 ml of medium

STORAGE

10-25°C

PACKAGING

Ref.	Content	Packaging
30095	10 tubes x 12 ml	10 tubes in cardboard box

pH OF THE MEDIUM

7.3 ± 0.2

USE

SIM MEDIUM is a medium used for the differentiation of microorganisms on the basis of hydrogen sulphide production, indole production and motility

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Slightly opalescent amber medium

SHELF LIFE

1 year

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 22 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 UFC/ml
Inoculum for selectivity: 10⁴-10⁵ UFC/ml
Inoculum for specificity: ≤10⁴ UFC/ml
Incubation Conditions: 18-24 h at 36 ± 1°C, in aerobiosis

Microorganism		Growth	H ₂ S	Motility	Indole
<i>Escherichia coli</i>	ATCC 25922	Good	-	+	+
<i>Salmonella typhimurium</i>	ATCC 14028	Good	+	+	-
<i>Salmonella typhi</i>	ATCC 6539	Good	+	+	-
<i>Shigella flexneri</i>	ATCC 12022	Good	-	-	-

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@lioilchem.net