SIMMONS CITRATE AGAR

Differential medium for enterobacteria identification.

TYPICAL FORMULA (g/l)

| | |
|-------------------------------|------|
| Dipotassium Phosphate | 1.6 |
| Magnesium Sulphate | 0.2 |
| Ammonium Dihydrogen Phosphate | 1.0 |
| Sodium Citrate | 2.0 |
| Sodium Chloride | 5.0 |
| Brom Thymol Blue | 0.08 |
| Agar | 14.0 |
| Final pH 6.8 ± 0.2 | |

DESCRIPTION

SIMMONS CITRATE AGAR is recommended for the differentiation and identification of *Enterobacteriaceae* on the basis of citrate utilization.

PRINCIPLE

Organisms able to utilize ammonium dihydrogen phosphate and sodium citrate as the sole source of nitrogen and carbon respectively will grow on this medium and produce an alkaline reaction as evidenced by a change in the color of the bromthymol blue indicator from green (neutral) to blue (alkaline).

PREPARATION

Melt the content of one bottle in a boiling water-bath at 100°C (with cap partially unscrewed) until completely dissolved. After molten, verify the proper homogeneity of medium, turning upside down the bottle after screwing cap. If necessary, put the bottle in the water-bath to reach the total melting. Cool to 45-50°C, carefully mix avoiding bubbles formation. Aseptically distribute into Petri dishes. Carefully allow the plates to dry in order to avoid condensation formation that may cause medium contamination. Store the plates upside-down in tightly closed containers or polyethylene envelopes. Otherwise introduce the molten medium into sterile tubes and allow to solidify in a slant position.

TECHNIQUE

Inoculate the plates by streaking the sample to examine onto the entire surface using a sterile loop, or the tubes by stabbing the butt and streaking the slope. Incubate at 36±1°C for 24-48 hours.

INTERPRETATION OF RESULTS

Organisms able to utilize sodium citrate will grow with colonies.

STORAGE

10-25°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING and PRECAUTIONS

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of ≥1%. The product is designed for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

1. Ewig, W.H. and Edwards P.R. (1960). Bull. Bact. Nomen. and Taxon. 10: 1-12.
2. American Public Health Association (1981). Standards Methods for the Examination of Water and Wastewater, 15th ed. APHA Inc, Washington DC.
3. Matsen J.M., and Sherris J.C. (1969) Appl. Microbiol. 18: 452-454.

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME

SIMMONS CITRATE AGAR

PRESENTATION

Glass bottles containing 200 mL of medium.

STORAGE

10-25°C

PACKAGING

| Code | Content | Packaging |
|--------|-------------------|------------------------------|
| 412030 | 6 bottles x 200mL | 6 bottles in cardboard boxes |

pH OF THE MEDIUM

6.8 ± 0.2

USE

SIMMONS CITRATE AGAR is recommended for the differentiation and identification of *Enterobacteriaceae* on the basis of citrate utilization.

TECHNIQUE

Refer to technical sheet of the product.

APPEARANCE OF THE MEDIUM

Forest green medium, slightly opalescent, might have a slight precipitate.

SHELF LIFE

2 years

QUALITY CONTROL

Control of general characteristics, label and print

Sterility control

7 days at 25 +/- 1°C, in aerobiosis

7 days at 36 +/- 1°C, in aerobiosis

Microbiological control

Inoculum for productivity: 10-100 UFC/ml

Inoculum for selectivity: 10⁴-10⁵ UFC/ml

Inoculum for specificity: >10⁴ UFC/ml

Incubation conditions: 36 +/- 1°C for 24-48 hours

| Microorganisms | | Growth | Characteristics |
|-------------------------------|------------|-----------|-----------------|
| <i>Enterobacter aerogenes</i> | ATCC 13048 | Good | Blue colonies |
| <i>Escherichia coli</i> | ATCC 25922 | Inhibited | |
| <i>Salmonella typhimurium</i> | ATCC 14028 | Good | Blue colonies |
| <i>Salmonella typhi</i> | ATCC 19430 | Good | Green colonies |

TABLE OF SYMBOLS

| | | | |
|------------------------------------------------------|------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| IVD <i>In vitro</i> Diagnostic Medical Device | LOT Batch code |  Manufacturer |  Contains sufficient for <n> tests |
| REF Catalogue number |  Temperature limitation |  Use by |  Caution, consult accompanying documents |

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