CLOSTRIDIUM DIFFICILE AGAR

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CLOSTRIDIUM DIFFICILE AGAR

Medium for the isolation of *Clostridium difficile*.

TYPICAL FORMULA	(g/l)
Proteose Peptone	40.0
Fructose	6.0
Disodium Hydrogen Phosphate	5.0
Potassium Dihydrogen Phosphate	1.0
Magnesium Sulphate	0.1
Sodium Chloride	2.0
Agar	15.0
Final pH 7.4 ± 0.2 at 25°C	

DESCRIPTION

CLOSTRIDIUM DIFFICILE AGAR is a medium used for the isolation of *Clostridium difficile*.

PRINCIPLE

Proteose peptone provides nitrogen, vitamins, minerals and amino acid essential for growth. Fructose is the fermentable carbohydrate used to facilitate the recovery and growth of *C. difficile*. Potassium dihydrogen phosphate and disodium hydrogen phosphate act as buffering system. Magnesium sulphate is source of magnesium ions essential for several enzymatic reactions and DNA duplication. Sodium chloride maintains the osmotic balance of the medium. This basal medium needs to be added with horse blood that supplies essential growth factors and selective agents to inhibit the growth of most of the microorganisms present in fecal sample other than *C. difficile*.

PREPARATION

Suspend 34.6 g of powder in 500 ml of distilled water. Heat until completely dissolved. Autoclave at 121°C for 15 minutes. Cool to 45-50°C. Aseptically add 1 vial of CLOSTRIDIUM difficile *Supplement* (ref. 81007) previously reconstituted with 2 ml of sterile distilled water and 35 ml of HORSE BLOOD DEFIBRINATED (ref. 83395). Mix gently and dispense in Petri dishes.

TECHNIQUE

Direct inoculum of fecal samples

Lightly inoculate the medium with fecal sample spreading part of the original inoculum in order to obtain well separated colonies. Incubate the plates at 36±1°C for 18-24 hours anaerobically.

Treatment for alcohol shock

Mix equal part of absolute ethyl alcohol and the fecal specimen. Homogenize using a vortex mixer. Leave at room temperature of 1 hour. Inoculate the medium and incubate the plates at 36±1°C for 18-24 hours anaerobically.

INTERPRETATION OF RESULTS

C. difficile grows with whitish opaque colonies of 4-6mm diameter irregular.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

1. Levett (1985) J. Clin. Pathol. 38: 233-234.
2. Barlett, J.G. et al. (1978) N. Eng. J. Med., 298, 531.
3. Boriello, S.P. et al (1981) J. Antimicrob. Chemother. 7 Supp. A. 53-62.
4. George, R.H. et al (1976) J. Clin. Microbiol. 6, 214-219.

PRODUCT SPECIFICATIONS

NAME

CLOSTRIDIUM DIFFICILE AGAR

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGE

Ref.	Content	Packaging
610115	500 g	500 g of powder in plastic bottle
620115	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.4 ± 0.2

USE

CLOSTRIDIUM DIFFICILE AGAR is a medium used for the isolation of *Clostridium difficile*

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Dehydrated medium

Appearance: free-flowing, homogeneous

Colour: beige

Prepared medium

Appearance: opalescent

Colour: cherry red

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 10-100 CFU/ml
Inoculum for selectivity: 10⁴-10⁵ CFU/ml
Inoculum for specificity: ≤10⁴ CFU/ml
Incubation conditions: 18-24 h at 36 ± 1°C anaerobically

Microorganism		Growth
<i>Escherichia coli</i>	ATCC® 25922	Inhibited
<i>Clostridium difficile</i>	ATCC® 9869	Good

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> diagnostic medical device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instruction for use	 Do not reuse
 Keep away from heat sources				

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LIOFILCHEM® S.r.l.