VRE AGAR BASE

Medium for the isolation of vancomycin-resistant enterococci.

TYPICAL FORMULA	(g/l)
Tryptone	17.0
Peptone	3.0
Yeast Extract	5.0
Sodium Chloride	5.0
Ox Bile	10.0
Esculin	1.0
Ferric Ammonium Citrate	0.5
Sodium Citrate	1.0
Sodium Azide	0.25
Agar	15.0
Final pH 7.1 ± 0.2 at 25°C	

DESCRIPTION

VRE AGAR BASE is a medium used for the isolation of vancomycin-resistant enterococci.

PRINCIPLE

Tryptone and peptone are sources of carbon, nitrogen, amino acids and proteins. Yeast extract provides vitamins of B-group which promote bacterial growth. Sodium chloride maintains the osmotic balance of the medium. Ox bile inhibits gram-positive bacteria, except for enterococci. Esculin is a glycoside hydrolysed by enterococci to yield esculetin and glucose. Esculetin reacts with the ferric ions to produce a black colored complex that surrounds the colonies. Sodium azide inhibits gram-negative bacteria. Agar is the solidifying agent.

VRE AGAR BASE is supplemented with antibiotics, i.e. vancomycin and meropenem. Vancomycin is effective against gram-positive bacteria, other than vancomycin-resistant enterococci. Meropenem suppress the contaminant microbial flora, in particular gram-negative bacteria and *Enterococcus gallinarum*.

PREPARATION

Suspend 57 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Aseptically, add rehydrated content of 2 vials of Vancomycin Supplement (ref. 81062) and 1 vial of Meropenem Supplement (ref. 81083). Mix well and pour in Petri dishes.

TECHNIQUE

Inoculate the medium to obtain well isolated colonies. Incubate at 36 ± 1°C for 18-48 h.

INTERPRETATION OF RESULTS

Enterococcus spp. cultivate with colonies surrounded by a black or dark-brown halo.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain harmful substances in concentrations higher than the limits set by current legislation, therefore it is not classified as dangerous; however, for its correct use, it is recommended to consult the safety data sheet. The product is intended for *in vitro* diagnostic use only and must be used by qualified personnel.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- King, W. K. (1996) Bug Bytes Vol. 2 No. 19.
- CDC Preventing the spread of vancomycin resistance: a report from the Hospital Infection Control Practices Advisory Committee (1994). Fed Regist. May 17.
- Gold, H. S. & Moellering R.C. Jr. (1996) N. Engl. J. Med.; 335(19): 1445-53.
- Weinbren, M. J., Johnson, A.P. & Woodford, N. (2000) J. Antimicrobial Chemotherapy; 45 :404-405.

PRODUCT SPECIFICATIONS

NAME

VRE AGAR BASE

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
610501	500 g	500 g of powder in plastic bottle

pH OF THE MEDIUM

7.1 ± 0.2

USE

VRE AGAR BASE is a medium used with supplements for the selective isolation of vancomycin-resistant enterococci

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

Ready-to-use medium

Appearance: slightly opalescent

Colour: amber

SHELFLIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
 Inoculum for productivity: 10-100 CFU/ml
 Inoculum for selectivity: 10⁴-10⁵ CFU/ml
 Inoculum for specificity: ≤10⁴ CFU/ml
 Incubation Conditions: 18-48 h at 35 ± 2°C

Microorganism		Growth	Specification
<i>Escherichia coli</i>	ATCC 25922	Inhibited	---
<i>Enterococcus faecalis</i>	ATCC 51299	Good	Blackening of the medium
<i>Enterococcus faecalis</i>	ATCC 29212	Inhibited	---

TABLE OF SYMBOLS

 Batch code	 In Vitro Diagnostic Medical device	 Manufacturer	 Use by
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use

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LIOFILCHEM® S.r.l.