Sabouraud CAF Agar + Actidione

Selective medium for isolation and cultivation of pathogenic fungi from clinical specimens.

TYPICAL FORMULA	(g/l)
Papaic Digest of Soybean Meal	10.0
Glucose	10.0
Chloramphenicol	0.05
Actidione (Cycloheximide)	0.4
Agar	15.0
Final pH 7.0 ± 0.2 at 25°C	

DESCRIPTION

Sabouraud CAF Agar + Actidione is a highly selective medium used for the isolation and maintenance of fungi, including dermatophytes, from clinical specimens or materials with mixed flora.

PRINCIPLE

Papaic digest of soybean meal provides amino acids, carbon, nitrogen, vitamins and minerals. Glucose is an energy source. Chloramphenicol is a broad-spectrum antibiotic inhibitory to a wide range of Gram-negative and Gram-positive bacteria. Cycloheximide suppresses the growth of saprophytic fungi. Agar is the solidifying agent.

PREPARATION

Suspend 35.5 g of powder in 1 liter of distilled water. Mix well and dissolve by heating with frequent agitation. Sterilize in autoclave at 118°C for 15 minutes. Cool to 45-50°C before distributing into Petri dishes.

TECHNIQUE

Streak the specimen directly over the agar surface. Incubate aerobically at 25-30°C for 2-7 days.

Note: Streak the specimen as soon as possible after it is received in the laboratory

For isolation of fungi causing systemic mycoses, two sets of plates should be inoculated, one set incubated at 25-30°C and the other one at 35 ± 2°C.

To provide an indication of all pathogens present in the specimen, two other media might be inoculated as well, such as Sabouraud Dextrose Agar (ref. 10035) and Columbia Sheep Blood (ref. 11025).

INTERPRETATION OF RESULTS

Examine the plates for fungal colonies exhibiting typical color and morphology.

Biochemical tests and serological procedures should be performed to confirm findings.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store the prepared medium at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product contains hazardous substances. It is recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used by properly trained operators only.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

- Wehr and Frank (2004) Standard methods for the examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
- Larone (1995) Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
- Stein, R.J.(1973) Handbook of Phycological Methods, Cambridge:University Press.
- Booth, C. (1971) Methods in microbiology, Vol. 4, London: Academy press.
- Sabouraud. 1892. Ann. Dermatol. Syphil. 3:1061.

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10-30°C

PACKAGE

Ref.	Content	Packaging
610179	500 g	500 g of powder in plastic bottle
620179	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.0 ± 0.2

USE

Sabouraud CAF Agar + Actidione is a selective medium used for the isolation and cultivation of pathogenic fungi from clinical specimens or other materials containing mixed flora

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Dehydrated medium

Appearance: free-flowing, homogeneous

Colour: off-white

Prepared medium

Appearance: slightly opalescent

Colour: amber

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Incubation conditions: 2-7 days at 25-30°C under aerobic atmosphere

Microorganism		Growth
<i>Candida albicans</i>	ATCC® 10231	Good
<i>Trichophyton mentagrophytes</i>	ATCC® 9533	Good
<i>Aspergillus niger</i>	ATCC® 16404	Partially to completely inhibited
<i>Saccharomyces cerevisiae</i>	ATCC® 9763	Inhibited
<i>Escherichia coli</i>	ATCC® 8739	Inhibited

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instruction for use	 Keep away from heat sources

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LIOFILCHEM® S.r.l.