COPRO SYSTEM *Plus*

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COPRO SYSTEM Plus

System for detection and presumptive identification of intestinal pathogenic microorganisms

DESCRIPTION

COPRO SYSTEM Plus is a 18 wells system containing desiccated biochemical substrata for the detection, presumptive identification of microorganisms from fecal samples (fecal culture).
The system is inoculated with the suspension of the specimen and incubated at 36 ± 1 °C for 18-24 hours.
The tests for the detection and presumptive identification of microorganisms present in the specimen are interpreted by evaluating the color change in the wells and performing immunologic and microscopic confirmation tests.

CONTENT OF THE PACKAGE

The package contains:

20 Systems COPRO SYSTEM Plus	1 Instruction sheet
20 Physiological Solution vials (4.5 mL/vial)	1 TEST RESULTS FORM

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

<i>Salmonella</i> antiserum (ref.96149)	OXIDASE TEST DISC (ref. 88004)
<i>Yersinia enterocolitica</i> antiserum (ref.96147)	ENTEROSYSTEM 18R REAGENT (ref. 80252)
<i>Shigella</i> antiserum (ref.96148)	Sterile containers for fecal culture (ref. 96708)
E.COLI O157 LATEX (ref.96150)	Collection swabs
CAMPYLOBACTER LATEX (ref.96143)	Cover slides and slides
	Microscope

CONFIGURATION

The system has the configuration indicated in table no.1.

Table no.1

Well	PRESUMPTIVE IDENTIFICATION of <i>Salmonella</i> spp.
1- LDC	Lysine decarboxylase
2- ORN	Ornithine decarboxylase
3- H₂S	Dihydrogen sulphide production
4- SAL*	Immunoserological test for confirmation <i>Salmonella</i> spp.
Well	PRESUMPTIVE IDENTIFICATION of <i>Proteus</i> spp.
5- PRO	<i>Proteus</i> spp./ <i>Providencia</i> spp.
6- UR	Urea hydrolysis
Well	PRESUMPTIVE IDENTIFICATION of <i>Pseudomonas</i> spp.
7- PSE	<i>Pseudomonas</i> spp.
8- OX*	Cytochrome oxidase test for confirming <i>Pseudomonas</i> spp
Well	PRESUMPTIVE IDENTIFICATION of <i>Yersinia enterocolitica</i>
9- CIN	Cultural medium for growth of <i>Yersinia enterocolitica</i>
10- YER*	Immunoserological test for confirming <i>Yersinia enterocolitica</i>
Well	PRESUMPTIVE IDENTIFICATION of <i>Shigella</i> spp.
11- SHI*	Immunoserological test for confirming <i>Shigella</i> spp.
Well	PRESUMPTIVE IDENTIFICATION of <i>Campylobacter jejunii</i>
12- CAM*	Immunoserological test for confirming <i>Campylobacter jejunii</i>
Well	PRESUMPTIVE IDENTIFICATION of <i>Escherichia coli</i> and enteropathogenic <i>E.coli</i> O157
13- ESC*	<i>Escherichia coli</i> , enteropathogenic <i>E.coli</i> O157
14- IND*	Indole test for confirming <i>Escherichia coli</i>
Well	PRESUMPTIVE IDENTIFICATION of KES group microorganisms
15- KES	KES Group (<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i>)
16- VP*	Voges Proskauer test for confirming of KES group microorganisms
Well	PRESUMPTIVE IDENTIFICATION of <i>Candida</i> spp.
17- CAN	<i>Candida</i> spp.
18- OBS	Microscopical examination for confirmation of <i>Candida</i> spp.

: After inoculation, cover the wells with three drops of Vaseline Oil

* : After incubation, add the indicated reagent to performe the test

• : After incubation, performe the agglutination test

PRINCIPLE OF THE METHOD

COPRO SYSTEM Plus allows detection and presumptive identification of pathogenic and non-pathogenic microorganisms, that can be present in fecal specimens, such as: *Salmonella spp.*, *Proteus spp.*, *Pseudomonas spp.*, *Yersinia spp.*, *Shigella spp.*, *Campylobacter jejunii*, *Escherichia coli*, enteropathogenic *E.coli O157*, microorganisms of KES group (*Klebsiella*, *Enterobacter*, *Serratia*), *Candida spp.*

- The presence of ***Salmonella spp.*** is detected by the yellow to red color change of well **1-LDC**; by the yellow to red color change of well **2-ORN**; by the yellow to black color change of well **3-H₂S** and confirmed through the serum agglutination test performed with ***Salmonella antiserum*** (ref.96149) from microbial culture of the well **4-SAL**¹.
- The presence of ***Proteus spp./ Providencia spp.*** is detected by the yellow to brown-black color change of well **5-PRO** and confirming by the yellow to red-fuchsia of well **6-UR**².
- The presence of ***Pseudomonas spp.*** is detected by the yellow to turbid green color change of well **7-PSE** and confirmed through the biochemical test with **OXIDASE TEST DISCS** (ref. 88004) performed from microbial culture of the well **8-OX**⁴.
- The presence of ***Yersinia enterocolitica*** is detected by the pink to turbid red color change of well **9-CIN** and confirmed through serum agglutination test performed with ***Yersinia enterocolitica antiserum*** (ref.96147) from microbial culture of the well **10-YER**³.
- The presence of ***Shigella spp.*** is detected by the serum agglutination test performed with ***Shigella antiserum*** (ref.96143) from microbial culture of well **11-SHI**¹.
- The presence of ***Campylobacter jejunii*** is detected through serum agglutination test performed with **CAMPYLOBACTER LATEX** (ref.96143) from microbial culture of the well **12-CAM**.
- The presence of ***Escherichia coli*** and of enteropathogenic ***E. coli O157*** is detected by the red to blue color change of well **13-ESC** and by the indole test performed in the well **14-IND**. The presence of enteropathogenic *E.coli O157* is confirmed through serum agglutination with **E.COLI O157 LATEX** (code 96150) from microbial culture of the well **13-ESC**¹.
- The presence of **KES group microorganisms** (*Klebsiella spp.*, *Enterobacter spp.*, *Serratia spp.*) is detected by the purple to yellow color change of well **15-KES** and through Voges-Proskauer test performed in the well **16-VP**. The identification of various microorganisms has to be confirmed through isolation on selective media for enterobacteria and through biochemical tests².
- The presence of ***Candida spp.*** is detected by the green to yellow color change and by the formation of a sediment in the well **17-CAN**⁵. Watch for the presence of chlamydo spores and mycelial hyphae at the microscope taking a drop from the microbial culture of the well **18-OBS**.

COMPOSITION

Table n°2

Well	Content
1-LDC	Culture medium with substrate for lysine decarboxylation
2-ORN	Culture medium with substrate for ornithine decarboxylation
3-H₂S	Culture medium with substrate for dihydrogen sulphide production
4-SAL	Culture medium for <i>Salmonella spp.</i> immunoserological test
5-PRO	Culture medium for the growth of <i>Proteus spp./ Providencia spp.</i>
6-UR	Culture medium with substrate for urea hydrolysis
7-PSE	Culture medium for the growth of <i>Pseudomonas spp.</i>
8-OX	Culture medium for the performance of Cytochrome-oxidase test
9-CIN	Selective culture medium for the growth of <i>Yersinia enterocolitica</i>
10-YER	Culture medium for <i>Yersinia enterocolitica</i> immunoserological test
11-SHI	Culture medium for <i>Shigella spp.</i> immunoserological test
12-CAM	Culture medium for <i>Campylobacter jejunii</i> immunoserological test
13-ESC	Culture medium for the growth of <i>Escherichia coli</i> and enteropathogenic <i>E.coli O157</i>
14-IND	Culture medium for the performance of Indole Test
15-KES	Culture medium for the growth of KES group microorganisms (<i>Klebsiella spp.</i> , <i>Enterobacter spp.</i> , <i>Serratia spp.</i>)
16-VP	Culture medium for the performance of Voges Proskauer Test
17-CAN	Selective culture medium for the growth of <i>Candida spp.</i>
18-OBS	Culture medium for the growth of spores and mycelial hyphae of yeasts <i>Candida spp.</i>

Physiological Solution (g/L): Sodium Chloride 0,9 g; Distilled water 1000.0 mL; pH 7.0 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLE

Take the fecal sample with proper sterile containers for the collection.

In special cases the collection of the fecal specimen with a swab may be useful.

The samples to be tested must be collected according to standard microbiological methods suggested for fecal culture ⁷.

TEST PROCEDURE

PREPARATION OF THE CLINICAL SPECIMEN

A) FECES

- Homogenize with care a portion of feces of about 1g in a tube containing 5 mL of Physiological Solution.
- Transfer 0.2 mL of homogenized sample into a vial of Physiological Solution contained in the kit.
- Shake and wait for 5 minutes before inoculating the system.

B) RECTAL SWAB

- Perform the rectal swab as reported in clinical protocols and methods.
- Dip the swab into a vial of *Physiological Solution* contained in the kit.
- Leave the swab immersed for 5 minutes before inoculating the system.

INOCULATION OF THE SYSTEM

1. Take a system from its pouch and leave it at room temperature.
2. Note name of the patient, date of beginning of the examination and type of clinical material.
3. Transfer 0.5 mL (4 drops) of the resulting suspension into each well of the system.
4. Cover with three drops of vaseline oil for microbiological use wells: **1-LDC, 2-ORN, 3-H₂S, 6-UR**.
5. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.

INTERPRETATION OF RESULTS

Subsequently to the incubation:

1. Watch for the color change in the well **1-LDC, 2-ORN** and **3-H₂S** and interpret the result using the table no.3
2. A positive reaction of **1-LDC, 2-ORN, 3-H₂S** wells indicates a presumptive identification of *Salmonella spp.*
3. Perform the rapid serum agglutination test **Salmonella antiserum** (ref.96149) from the microbial culture of **4-SAL** well following the test procedure as in the pack insert to confirm the presence of *Salmonella spp.*
4. Watch for the color change in the well **5-PRO** and **6-UR** and interpret the result using the table no.3
5. A positive reaction of **5-PRO** and **6-UR** wells indicates a presumptive identification of *Proteus spp.*
6. Watch for the color change in the well **7-PSE** and interpret the result using the table no.3.
7. Perform the oxidase test by adding a disc of **OXIDASE TEST DISC** (ref.88004) into the well **8-OX** and watch for the formation of a blue-purple coloring within 2 minutes to confirm the presence of *Pseudomonas spp.*
8. Watch for the color change in the well **9-CIN** and interpret the result using the table no.3.
9. A positive reaction of the well **9-CIN** indicates a presumptive identification of *Yersinia enterocolitica*
10. Perform the rapid serum agglutination test with **Yersinia enterocolitica antiserum** (ref.96147) from microbial culture of **10-YER** well following the test procedure as in the pack insert to confirm the presence of *Yersinia enterocolitica*
11. Perform the rapid serum agglutination test with **Shigella antiserum** (ref.96148) from microbial culture of **11-SHI** well following the test procedure as in the pack insert to confirm the presence of *Shigella spp.*
12. Perform the rapid serum agglutination test with **CAMPYLOBACTER KIT LATEX** (ref.96143) from microbial culture of **12-CAM** well following the test procedure as in the pack insert to confirm the presence of *Campylobacter jejuni*.
13. Watch for the color change in the well **13-ESC** and interpret the result using the table no.3 .
14. Add 2 drops of Kovac's Reagent (ref.80270) to the well **14-IND** and wait for the appearance of a red color.
15. Perform the rapid serum agglutination test with **E.COLIO157 LATEX** (ref.96150) from microbial culture of **13-ESC** well following the test procedure as in the pack insert to confirm the presence of enteropathogenic *E.coli* O157.
16. Watch for the color change in the well **15-KES**, interpret the result using the table no.3 .
17. Add 2 drops of alpha-naphthol and 1 drop of NaOH 40% into the well **16-VP** (ref. 88035) and wait 15-20 minutes until a pink-red color appears.
18. Perform the microbial isolation test: take a loop of culture broth from the well **15-KES** and inoculate on the surface of a selective plate for enterobacteria. Incubate the plate at 36 ± 1 °C and carry out biochemical tests for the cultivated colonies to confirm the presence of *KES* group microorganisms
19. Watch for the color change in the well **17-CAN**, interpret the result using the table no.3.
20. Take a drop from the well **17-CAN**, deposit it on a glass slide, cover it, and watch for chlamydo spores and mycelial hyphae at the microscope (40x) to confirm the presence of *Candida spp.*
21. Note the results on the **TEST RESULTS FORM** (copy as many forms as necessary).

Table No. 3

Well	PRESUMPTIVE IDENTIFICATION of <i>Salmonella spp.</i>	Well color	
		Positive reaction	Negative reaction
1-LDC	Lysine decarboxylation	red	yellow
2-ORN	Ornithine decarboxylation	red	yellow
3-H ₂ S	Dihydrogen sulphide production	black	yellow
4-SAL	<i>Salmonella spp.</i> : immunoserological agglutination test	agglutination	no agglutination
Well	PRESUMPTIVE IDENTIFICATION of <i>Proteus spp.</i>	Well color	
		Positive reaction	Negative reaction
5-PRO	<i>Proteus spp./ Providencia spp.</i>	brown-black	yellow
6-UR	Urea hydrolysis	red-fuchsia	yellow
Well	PRESUMPTIVE IDENTIFICATION of <i>Pseudomonas spp.</i>	Well color	
		Positive reaction	Negative reaction
7-PSE	<i>Pseudomonas spp.</i>	turbid green	Yellow
8-OX	Cytochrome-oxydase test to confirm <i>Pseudomonas spp.</i>	blue	colorless
Well	PRESUMPTIVE IDENTIFICATION of <i>Yersinia enterocolitica</i>	Well color	
		Positive reaction	Negative reaction
9-CIN	<i>Yersinia enterocolitica</i>	turbid red	clear pink
10-YER	<i>Yersinia enterocolitica</i> : immunoserological agglutination test	agglutination	no agglutination
Well	PRESUMPTIVE IDENTIFICATION of <i>Shigella spp.</i>	Well color	
		Positive reaction	Negative reaction
11-SHI	<i>Shigella spp.</i> : immunoserological agglutination test	agglutination	no agglutination
Well	PRESUMPTIVE IDENTIFICATION of <i>Campylobacter jejunii</i>	Well color	
		Positive reaction	Negative reaction
12-CAM	<i>Campylobacter jejunii</i> : immunoserological agglutination test	agglutination	no agglutination
Well	PRESUMPTIVE IDENTIFICATION of <i>Escherichia coli, E.coli O157</i>	Well color	
		Positive reaction	Negative reaction
13-ESC	<i>Escherichia coli, E.coli O157</i>	blue	grey-red
	<i>E.coli O157</i> : immunoserological agglutination test	agglutination	no agglutination
14-IND	Indole test	pink-red ring	colorless
Well	PRESUMPTIVE IDENTIFICATION of KES group microorganisms (<i>Klebsiella spp., Enterobacter spp., Serratia spp.</i>)	Well color	
		Positive reaction	Negative reaction
15-KES	KES group (<i>Klebsiella spp., Enterobacter spp., Serratia spp.</i>)	yellow	purple
16-VP	Voges-Proskauer test	pink-red	colorless
Well	PRESUMPTIVE IDENTIFICATION of <i>Candida spp.</i>	Well color	
		Positive reaction	Negative reaction
17-CAN	<i>Candida spp.</i>	yellow	green
18-OBS	Microscopic examination (40x)	presence of chlamydospores and mycelial hyphae	absence of chlamydospores and mycelial hyphae

QUALITY CONTROL

Every batch of **COPRO SYSTEM Plus** is subjected to the quality control using the following reference microorganisms:

<i>Candida albicans</i>	ATCC 10231	<i>Shigella flexneri</i>	ATCC 12022
<i>Proteus mirabilis</i>	ATCC 25933	<i>Escherichia coli O157</i>	ATCC 35150
<i>Escherichia coli</i>	ATCC 25922	<i>Salmonella typhimurium</i>	ATCC 14028
<i>Klebsiella pneumoniae</i>	ATCC 13883	<i>Yersinia enterocolitica</i>	ATCC 9610
<i>Pseudomonas aeruginosa</i>	ATCC 27853	<i>Campylobacter jejunii</i>	ATCC 33291

FACTORS THAT MAY INVALIDATE THE RESULTS

- Poor standardization of the inoculum.
- Unsuitable clinical material.
- Use of expired systems and/or reagents.
- Non compliance with temperatures and times of incubation.

CAUTIONS

The product, **COPRO SYSTEM Plus**, is not classified as hazardous under current legislation, however refer to the safety data sheet for a correct use. **COPRO SYSTEM Plus** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **COPRO SYSTEM Plus** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Package
COPRO SYSTEM Plus	71675	20 tests

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

BIBLIOGRAFIA / BIBLIOGRAPHY

1. Cheryl A. Bopp, Frances W. Brenner, Joy G. Wells and Nancy A. Strockbine. *Escherichia, Shigella and Salmonella* . pg. **459-474** in Manual of Clinical Microbiology 7th Edition (1999).
2. JJ Farmer III. *Enterobacteriaceae: Introduction and Identification* pg. **442-458** in Manual of Clinical Microbiology 7th Edition (1999).
3. Stojanka Aleksic and Jochen Bockemuhl. *Yersinia and other Enterobacteriaceae*.pg. **483-496** in Manual of Clinical Microbiology 7th Edition (1999).
4. Deanna L:Kiska and Peter H.Gilligan. *Pseudomonas* pg **517-525** in Manual of Clinical Microbiology 7th Edition (1999).
5. Davise H. Larone. Medically important fungi. A Guide to Identification. 2th Edition (1987). Elsevier.
6. Ronald M.Atlas. *Handbook of Microbiological Media* (1997) Lawrence C. Parks.
7. Bailey and Scott's. Diagnostic Microbiology. 7th Edition. The C.V. Mosby Company.1986

Bacteriology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +390858930745 Fax +390858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net