MIC Test Strip

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DESCRIPTION

MIC Test Strip is a quantitative assay for determining the Minimum Inhibitory Concentration (MIC) of antimicrobial agents against microorganisms and for detecting the resistance mechanisms.

MIC Test Strip are paper strips with special features* that are impregnated with a predefined concentration gradient of antibiotic, across 15 two-fold dilutions of a conventional MIC method.

On one side of the strip is indicated a MIC scale in $\mu\text{g/mL}$ and a code that identify the antimicrobial agent.

For ESBL (Extended Spectrum Beta Lactamase) and MBL (Metallo Beta Lactamase) detection, the double-sided gradient carries the appropriate diagnostic reagents.

MIC Test Strip are available in a large variety of configurations. Each configuration is available in packages of 10, 30 and 100 tests.

CONTENTS OF THE PACKAGES

The 10-test box contains 10 strips individually packed in desiccant envelopes and an instruction sheet.

The 30-test box contains 30 strips individually packed in desiccant envelopes and an instruction sheet.

The 100-test box contains 10 desiccant envelopes, each containing 10 strips, and an instruction sheet. The 100-test pack also contains a storage tube.

METHOD PRINCIPLE

When the **MIC Test Strip** is applied onto an inoculated agar surface, the preformed exponential gradient of antimicrobial agent is immediately transferred to the agar matrix.

After 18 hours incubation or longer, a symmetrical inhibition ellipse centered along the strip is formed. The MIC is read directly from the scale in terms of $\mu\text{g/mL}$ at the point where the edge of the inhibition ellipse intersects the strip **MIC Test Strip**. Other growth/inhibition patterns may also be seen for resistance detection methods.

COMPOSITION

The strips are made of high-quality paper and each strip is impregnated with a predefined concentration gradient across 15 two-fold dilutions of antibiotic agent.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the evaluation of Minimum Inhibition Concentration (MIC) are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before inoculation.

TEST PROCEDURE

1. Allow unopened envelop to come to room temperature before opening it, for minimising condensation on the strip.
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 mL of a suitable suspension medium. Fastidious microorganisms should be suspended in broth and used within 15 minutes.
3. Compare the turbidity to the appropriate McFarland standard.
4. Dip a sterile swab in the broth culture or in a diluted form thereof and squeeze it on the wall of the test tube to eliminate excess liquid.
5. Drag it along the surface of the medium contained on the plate so as to produce even growth; allow excess moisture to be absorbed and ensure that the surface is completely dry before applying strips.
6. Apply the strip to the agar surface with the MIC scale facing upwards and code of the strip to the outside of the plate, pressing it with a sterile forceps on the surface of the agar and ensure that whole length of the antibiotic gradient is in complete contact with the agar surface. Once applied, do not move the strip.
7. Incubate plates in an inverted position under conditions appropriate for the microorganism.
8. Put the not used strips onto the tube contained in the package.

EVALUATING THE RESULTS

At the end of incubation read the MIC value where the edge of the inhibition ellipse intersects the strip (intersection between two scale segments should be round up to the higher value). MIC break points for defining susceptibility categories as provided by CLSI or EUCAST could be used for interpreting MIC values.

Always round up **MIC Test Strip** half dilution values to the next upper two-fold value before categorisation. An overview of CLSI and EUCAST interpretative criteria is provided in Table no.1.

CLINICAL INTERPRETATION

The test **MIC Test Strip** carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population. The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

QUALITY CONTROL

Each batch of **MIC Test Strip** is subjected to precise and thorough checks in compliance with CLSI standards using the bacterial strains indicated in the table no.1.

PRECAUTIONS

The **MIC Test Strip** cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitisation in sensitive subjects if they come into contact with the skin.

MIC Test Strip are disposable products. **MIC Test Strip** are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

The unopened package of **MIC Test Strip** should be stored at -20°C until the given expiry date.

Leftover **MIC Test Strip** from an opened package must be stored at $2-8^{\circ}\text{C}$ in the airtight tube, containing desiccant, provided in the pack for no more than 7 days.

Do not store near sources of heat and do not expose to excessive temperature variations.

ELIMINATING USED MATERIAL

After use, **MIC Test Strip** and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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- CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 27th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
- CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically; Approved Standard - Tenth Edition. CLSI document M07-A10. Wayne, PA: Clinical and Laboratory Standards Institute; 2015.
- CLSI. Methods for Dilution Antimicrobial Susceptibility Testing of Anaerobic Bacteria. Approved Standard - Seventh Edition. CLSI document M11-A7. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
- CLSI. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard - Third Edition. CLSI document M27-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Fourth Informational Supplement. CLSI document M27-S4. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; Approved Standard - Third Edition. CLSI document M31-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 7.0, 2017. <http://www.eucast.org>
- The European Committee on Antimicrobial Susceptibility Testing. Antifungal Agents. Breakpoint tables for interpretation of MICs. Version 8.0, 2015. <http://www.eucast.org>

- Savini, V. et al. (2016). Daptomycin-resistant *Staphylococcus pettenkoferi* of human origin. *Acta Biochim Pol*; 63(2): 297-301.
- Koeth, L. M. et al. (2016). Evaluation of Ceftolozane-Tazobactam MIC Test Strip Compared to Broth Microdilution MIC for Enterobacteriaceae and *Pseudomonas aeruginosa*. ASM Microbe poster.
- Koeth, L. M. et al. (2016). Multi-Site Evaluation of Dalbavancin and Vancomycin MIC Test Strip Compared To Broth Microdilution MICs. ASM Microbe poster.
- Carretto, E. et al. (2016). Multicentric evaluation of the reliability and the reproducibility of synergy testing using the MIC test strip - synergy application system (MTS-SAS™). ECCMID poster 1009
- Koeth, L. M. et al. (2015) Evaluation of Dalbavancin MIC Test Strip (MTS) Compared to Broth Microdilution MIC for Relevant Gram Positive Isolates. ICAAC poster D-1138.
- Savini, V. et al. (2014). Liofilchem® Chromatic VRE and vancomycin MIC Test Strip detected glycopeptide resistance in a vanB neonatal *Enterococcus faecium* isolate showing alternate vancomycin susceptibility and resistance with bioMérieux Vitek2. *Int J Clin Exp Pathol*; 7(9):6274-6277.
- Canton, R. et al. (2014). Determination of the sensitivity and specificity of the MIC Test Strips (MTS) mechanisms of detection against in-house methods for 644 multi-drug resistant (MDR) strains: a European multi-centre study. ECCMID poster 1009.
- D'Humières, C. et al (2014). Direct determination of antimicrobial susceptibility of Gram Negative Bacilli from respiratory samples on chromogenic agar plates with gradient antibiotic strips. ECCMID poster 348.
- Tyrrell, J. and Walsh, T.R. (2014). Efficacy of Fosmidomycin alone and in combination with colistin, tigecycline and rifampicin against multi-drug resistant and extensively drug resistant Enterobacteriaceae. ECCMID poster 235.
- Walsh, T. et al. (2013). Synergistic effects of rifampicin, nitrofurantoin, fosfomycin and colistin against NDM-1 positive *Klebsiella pneumoniae* and *Escherichia coli* using the new MTS synergy method. ECCMID poster 1553.
- Walsh, T.R. et al. (2013). Evaluation of MIC Test Strips (MTS) gradient-diffusion system for susceptibility testing of NDM-1 positive *Klebsiella pneumoniae* and *Escherichia coli*. ECCMID poster 1554.
- Meletiadiis, J. et al. (2013). Evaluation of gradient concentration strips for in vitro combination testing of antifungal combinations against *Candida* spp. ECCMID poster 976.
- Meletiadiis, J. et al. (2013). Evaluation of strips with concentration gradient of antifungal agents for the susceptibility determination of *Candida* strains ECCMID poster 1583.
- Brocco, F. et al. (2013). Evaluation of a new method for detecting KPC-producing bacteria based on the combination of a screening medium and a gradient diffusion system. ECCMID poster eP690.
- M'Zali, F.H. et al. (2013). A novel, direct susceptibility testing method of sonicated vascular prosthetic graft samples by combination of Liofilchem® MIC Test Strip and Chromatic™ agar plates. ECCMID poster 1562.
- Haldorsen, B. C. et al. (2012). Agreement of the MIC Test Strip vs. Etest in the MIC determination of *Streptococcus pneumoniae*. ECCMID poster 1658.
- Rossolini, G.M. et al. (2011). Evaluation of a new gradient-diffusion system for MIC determination with Gram-negative pathogens. ECCMID, poster 572.
- Stefani, S. et al (2011). A new reliable screening method for the evaluation of VISA and hVISA strains by "Vancomycin-Teicoplanin MIC Test Strip" (VTMTS). ECCMID poster 776.

TABELLA DEI SIMBOLI / TABLE OF SYMBOLS / TABLEAU DES SYMBOLES / SYMBOLE / TABLA DE SIMBOLOS / TABELA DE SÍMBOLOS / SYMBOLFORKLARINGER / TABUĽKA SYMBOLOV / TABULKA SYMBOLŮ / SYMBOLTAVLE / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ / SYMBOLER

LOT **It:** Codice del lotto / **En:** Batch code / **Fr:** Numéro de lot / **De:** Charge / **Es:** Código del lote / **Pt:** Código do lote / **Nn:** Batch-nummer / **Sk:** Číslo šarže / **Cs:** Číslo šarže / **Da:** Lot nr. / **El:** Κωδικός παρτίδας / **Sv:** Lot nummer

It: Fabbricante / **En:** Manufacturer / **Fr:** Fabricant / **De:** Hersteller / **Es:** Fabricante / **Pt:** Fabricante / **Nn:** Produsent / **Sk:** Výrobca / **Cs:** Výrobce / **Da:** Producent / **El:** Κατασκευαστής / **Sv:** Tillverkare

REF **It:** Numero di catalogo / **En:** Catalogue number / **Fr:** Numéro du catalogue / **De:** Artikelnummer / **Es:** Número de catálogo / **Pt:** Número de catálogo / **Nn:** Artikelnummer / **Sk:** Katalógové číslo / **Cs:** Katalogové číslo / **Da:** Produkt nr. / **El:** Αριθμός καταλόγου / **Sv:** Artikel nummer

It: Utilizzare entro / **En:** Use by / **Fr:** Date d'approbation de l'enregistrement / **De:** Siehe Verfallsdatum auf dem Etikett / **Es:** Caducidad / **Pt:** Prazo de validade / **Nn:** Holdbarhetsdato / **Sk:** Použite do / **Cs:** Použitelné do / **Da:** Holdbarhedsdato / **El:** Χρήση έως / **Sv:** Använd senast

IVD **It:** Dispositivo medico diagnostico *in vitro* / **En:** In Vitro Diagnostic Medical Device / **Fr:** Utilisation pour diagnostic *in vitro* / **De:** für In Vitro Diagnostik / **Es:** Dispositivo de diagnóstico *in vitro* / **Pt:** Dispositivo de diagnóstico *in vitro* / **Nn:** *in vitro* medisinsk diagnostikk / **Sk:** Len na diagnostiku In Vitro / **Cs:** Určeno pro *in vitro* diagnostiku / **Da:** Til In Vitro diagnostisk brug / **El:** Ιατρική διαγνωστική συσκευή *in vitro* / **Sv:** För *invitro* diagnostik

It: Contenuto sufficiente per <n> saggi / **En:** Contains sufficient for <n> tests / **Fr:** Contenu suffisant pour <n> tests / **De:** Enthält Material für <n> Tests, siehe Etikett / **Es:** Contenido suficiente para <n> test / **Pt:** Conteúdo suficiente para <n> testes / **Nn:** Inneholdet er tilstrekkelig til <n> prøver / **Sk:** Balenie obsahuje <n> testov / **Cs:** Obsah dostatečný pro <n> testů / **Da:** Antal test indeholdt / **El:** Περιεχόμενο επαρκές για <n> δοκίμια / **Sv:** Innehåller <n> tester

It: Attenzione, vedere le istruzioni per l'uso / **En:** Caution, consult accompanying documents / **Fr:** Attention, consulter les instructions d'utilisation / **De:** Vorsicht, Begleitdokumente beachten / **Es:** Atención, leer las instrucciones de uso / **Pt:** Atenção, ler as instruções de utilização / **Nn:** Forsiktighet, kfr. medfølgende dokumentasjon / **Sk:** Upozornenie, preštudujte priložené informácie / **Cs:** Pozor, nahlédněte do přiložených dokumentů / **Da:** Bemærkning, se produktvejledning / **El:** Προσοχή, δείτε τις οδηγίες χρήσης / **Sv:** Varning, konsultera medföljande dokumentation

It: Limiti di temperatura / **En:** Temperature limitation / **Fr:** Limite de température / **De:** Temperaturbereich / **Es:** Límite temperatura de almacenamiento / **Pt:** Limite temperatura de conservação / **Nn:** Temperaturområde / **Sk:** Teplota skladovania / **Cs:** Rozmezí teplot / **Da:** Opbevaringstemperatur / **El:** Περιορισμοί θερμοκρασίας / **Sv:** Temperatur begränsningar

It: Limite superiore di temperatura / **En:** Upper limit of temperature / **Fr:** Limite supérieure de température / **De:** Temperaturobergrenze / **Es:** Limite superior de temperatura / **Pt:** Limite superior de temperatura / **Nn:** Øvre grense av temperatur / **Sk:** Horný limit teploty / **Cs:** Nejvyšší přípustná teplota / **Da:** Højeste temperatur / **El:** Ανώτερο όριο θερμοκρασίας / **Sv:** Högsta temperaturen

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RUO: For Research Use Only, not CE marked

CE Marking Unless otherwise indicated, all products intended to be used for clinical applications are CE marked according to the European Directive 98/79/EC for *In Vitro* Diagnostic Medical Devices.

MIC Test Strip for ANTIMYCOBACTERIAL susceptibility testing				
Description	µg/mL	code	packaging	ref.
MIC Test Strip ETHAMBUTOL	0.016 - 256	EB	10	921701
			30	92170
			100	921700
MIC Test Strip ETHIONAMIDE	0.016 - 256	ET	10	921721
			30	92172
			100	921720
MIC Test Strip ISONIAZIDE	0.016 - 256	IZ	10	921711
			30	92171
			100	921710

MIC Test Strip for ANTIFUNGAL susceptibility testing				
Description	µg/mL	code	packaging	ref.
MIC Test Strip AMPHOTERICIN B	0.002 - 32	AMB	10	921531
			30	92153
			100	921530
MIC Test Strip ANIDULAFUNGIN	0.002 - 32	AND	10	921551
			30	92155
			100	921550
MIC Test Strip CASPOFUNGIN	0.002 - 32	CAS	10	921541
			30	92154
			100	921540
MIC Test Strip FLUCONAZOLE	0.016 - 256	FLU	10	921471
			30	92147
			100	921470
MIC Test Strip FLUCYTOSINE	0.002 - 32	FC	10	921491
			30	92149
			100	921490
MIC Test Strip ISAVUCONAZOLE	0.002 - 32	IVU	10	921841
			30	92184
			100	921840
MIC Test Strip ITRACONAZOLE	0.002 - 32	ITC	10	921481
			30	92148
			100	921480
MIC Test Strip KETOCONAZOLE	0.002 - 32	KE	10	921511
			30	92151
			100	921510
MIC Test Strip MICAFUNGIN	0.002 - 32	MYC	10	921821
			30	92182
			100	921820
MIC Test Strip POSACONAZOLE	0.002 - 32	POS	10	921521
			30	92152
			100	921520
MIC Test Strip VORICONAZOLE	0.002 - 32	VO	10	921501
			30	92150
			100	921500

MIC Test Strip for Antimicrobial Resistance Testing				
Description	µg/mL	code	packaging	ref.
ESBL: confirmation of Extended Spectrum Beta-Lactamase				
MIC Test Strip CEFEPIME / CEFEPIME + CLAVULANIC ACID (4 µg/mL)	0.25-16 / 0.064-4	FEP/FEL	10 30 100	921611 92161 921610
MIC Test Strip CEFOTAXIME / CEFOTAXIME + CLAVULANIC ACID (4 µg/mL)	0.25-16 / 0.016-1	CTX/CTL	10 30 100	921601 92160 921600
MIC Test Strip CEFTAZIDIME / CEFTAZIDIME + CLAVULANIC ACID (4 µg/mL)	0.5-32 / 0.064-4	CAZ/CAL	10 30 100	921591 92159 921590
MBL: detection of Metallo Beta-Lactamase				
MIC Test Strip IMIPENEM / IMIPENEM + EDTA	4-256 / 1-64	IMI/IMD	10 30 100	921621 92162 921620
MIC Test Strip IMIPENEM / IMIPENEM + EDTA	0.125-8 / 0.032-2	IMI/IMD	10 30 100	921661 92166 921660
MIC Test Strip MEROPENEM / MEROPENEM + EDTA	0.125-8 / 0.032-2	MRP/MRD	10 30 100	921651 92165 921650
KPC detection				
MIC Test Strip MEROPENEM / MEROPENEM + PHENYLBORONIC ACID	0.125-8 / 0.032-2	MRP/MBO	10 30 100	921671 92167 921670
MIC Test Strip ERTAPENEM / ERTAPENEM + PHENYLBORONIC ACID	0.125-8 / 0.032-2	ETP/EBO	10 30 100	921681 92168 921680
AmpC detection				
MIC Test Strip CEFOTETAN / CEFOTETAN + CLOXACILLIN	0.5-32/0.5-32	CTT/CXT	10 30 100	921641 92164 921640
MIC Test Strip ERTAPENEM / ERTAPENEM + CLOXACILLIN	0.125-8 / 0.032-2	ETP/ECX	10 30 100	921691 92169 921690
GRD: Glycopeptide Resistance Detection				
MIC Test Strip VANCOMYCIN / TEICOPLANIN	0.5-32/0.5-32	VA/TEC	10 30 100	921631 92163 921630

Download the currently valid MIC Test Strip specific technical sheets for your applications:

Interpretative Criteria and Quality Control [click here](#)

FDA cleared items [click here](#)

Application Guide [click here](#)

MTS01 - Haemophilus influenzae [click here](#)

MTS02 - Fastidious Gram Positive Bacilli [click here](#)

MTS03 - Fastidious Gram Negative Bacilli [click here](#)

MTS04 - HACEK [click here](#)

MTS05 - Enterococci [click here](#)

MTS06 - Mould [click here](#)

MTS08 - Meningococci [click here](#)

MTS09 - Gram Negative Aerobes [click here](#)

MTS13 - Gonococci [click here](#)

MTS14 - Anaerobes [click here](#)

MTS16 - Campylobacter spp. [click here](#)

MTS19 - Streptococci [click here](#)

MTS20 - Staphylococci [click here](#)

MTS21 - Helicobacter pylori [click here](#)

MTS22- Nontuberculous Mycobacteria [click here](#)

MTS23 - Streptococcus pneumoniae [click here](#)

MTS24 - Yeast [click here](#)

MTS25 - Mycobacterium tuberculosis [click here](#)

MTS26 - ESBL [click here](#)

MTS27 - MBL [click here](#)

MTS28 - GRD [click here](#)

MTS29 - AmpC [click here](#)

MTS30 - Cidal-Static [click here](#)

MTS31 - Synergy Testing [click here](#)

MTS33 - Direct MIC [click here](#)

MTS34 - Enterococci VRE [click here](#)

MTS35 - KPC [click here](#)

MTS37 - Tedizolid [click here](#)

MTS38 - Dalbavancin [click here](#)

MTS39 - Ceftobiprole [click here](#)

MTS40 - Isavuconazole [click here](#)

MTS41 - Vancomycin [click here](#)

MTS42 - Ceftolozane-tazobactam [click here](#)

MTS43 - Ceftazidime-avibactam [click here](#)

MTS44 - Solithromycin [click here](#)

MTS45 - Fosfomycin [click here](#)

MTS47 - Cefepime [click here](#)

 **LIOFILCHEM® s.r.l.**
Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net

 In USA, available for products noted as "FDA Cleared" in the product list.

F00023
Rev.32 / 30.11.2017

liofilchem@liofilchem.net

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Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.
© Copyright LIOFILCHEM 2017