

Evaluation of a new method for detecting KPC-producing bacteria based on the combination of a screening medium and a gradient diffusion system

F. Brocco¹, D. Vitagliano¹, F. Demetrio¹

¹ Liofilchem®, Roseto degli Abruzzi, Italy

INTRODUCTION

Carbapenemes are commonly used to treat infections caused by multidrug-resistant Enterobacteriaceae. An increasing mechanism of carbapenem resistance is the Klebsiella pneumoniae carbapenemase (KPC). The KPC β -lactamase occurs most commonly in *K. pneumoniae*, but it has also been reported sporadically in other species of Enterobacteriaceae (*Klebsiella oxytoca*, *Enterobacter* spp., *Escherichia coli*, *Salmonella* spp., *Citrobacter freundii* and *Serratia* spp.) and *Pseudomonas aeruginosa*. The KPC enzyme confers resistance to all β -lactam agents including penicillins, cephalosporins, monobactams and most β -lactam/ β -lactamase inhibitors. KPC hydrolyzes the carbapenems (ertapenem, imipenem and meropenem) reducing their activity and poses a considerable threat to clinical patient care and public health.

OBJECTIVES

Laboratory identification of KPC-producing clinical isolates will be essential for limiting the spread of strains possessing this resistance mechanism. We evaluated a new method for the detection of potential KPC producers using first Chromatic CRE for the screening and then MTS Ertapenem/Ertapenem+ Boronic acid and Meropenem/Meropenem+ Boronic acid for confirmation.

METHODS

Bacterial isolates: 91 bacterial isolates with known resistance mechanisms (acquired *ampC*, MBL, KPC and ESBL producing Enterobacteriaceae) determined by molecular methods were included in the study.

All the investigated isolates were frozen at -20°C and subcultured on Mueller Hinton Agar plates for two consecutive days prior using to ensure the purity of the culture.

Screening medium

Chromatic CRE (Liofilchem®, Italy) is used as a chromogenic medium designed to detect carbapenem-resistant Enterobacteriaceae.

All plates are incubated aerobically at 36±1°C, read after 18-24 h and interpreted according to manufacturers' instructions.

Gradient diffusion system

MTS Ertapenem/Ertapenem+ Boronic acid (Liofilchem®, Italy) is a double-sided gradient strip of Ertapenem 8-0.125 mg/L and Ertapenem 2-0.032 plus a constant level of Boronic acid.

MTS Meropenem/Meropenem+ Boronic acid (Liofilchem®, Italy) is a double-sided gradient strip of Meropenem 8-0.125 mg/L and Meropenem 2-0.032 plus a constant level of Boronic acid.

All strips are inoculated on Mueller Hinton Agar (Liofilchem®, Italy). Interpretation of results was performed after 24 hours incubation.

All strains were tested on Chromatic CRE medium and with MTS Ertapenem/Ertapenem+ Boronic acid and Meropenem/Meropenem+ Boronic acid on Mueller Hinton Agar. Sensitivity and specificity were calculated.


RESULTS

KPC-producing *K. pneumoniae* grew uninhibited on Chromatic CRE and yielded green or blue colonies (100% sensitivity).

Poor growth/no growth/atypical colony color was evident with KPC-non producing Enterobacteriaceae (77% specificity).


MTS Ertapenem/Ertapenem+ Boronic acid and MTS Meropenem/Meropenem+ Boronic acid were evaluated with KPC positive and negative strains.

Both tests show 100% sensitivity and 89.3% specificity for MTS Ertapenem/Ertapenem+ Boronic acid and 95.6 % specificity for MTS Meropenem/ Meropenem+ Boronic acid, respectively.



growth on Liofilchem® Chromatic CRE

KPC-positive
Klebsiella pneumoniae



Liofilchem® MIC Test Strip
Ertapenem/Ertapenem+ Boronic acid
Meropenem/Meropenem+ Boronic acid

CONCLUSION

The study demonstrates the efficiency of Chromatic CRE intended as a screening step of KPC producing strains, followed by MTS Ertapenem/Ertapenem+ Boronic acid and MTS Meropenem/Meropenem+ Boronic acid as subsequent confirmatory steps; this combination allows the rapid detection of carbapenem-resistant organisms and the immediate implementation of infection control measures avoiding the spread of carbapenem resistance.

REFERENCES

- [K. F. Anderson, D. R. Lonsway, J. K. Rasheed, J. Biddle, B. Jensen, L. K. McDougal, R. B. Carey, A. Thompson, S. Stocker, B. Limbago and J.B.Patel. 2007](#)
Evaluation of Methods to Identify the Klebsiella pneumoniae Carbapenemase in Enterobacteriaceae.
Journal of Clinical Microbiology. Vol.45, N.8: 2723-2725.
- A. Tsakris, I.Kristo, A.Poulou, K. Themeli-Digalaki, A.Ikonomidis, D.Petropoulou, S.Pournaras, D.Sofianou .2009
Evaluation of boronic acid disk tests for differentiating KPC-possessing Klebsiella pneumoniae isolates in the clinical laboratory.
Journal of Clinical Microbiology. Vol.47, N.2: 362-367.
- D.Landman, J.Salamera, M.Singh, J.Quale. 2011.
Accuracy of Carbapenem Nonsusceptibility for Identification of KPC-Possesing Enterobacteriaceae by Use of the Revised CLSI Breakpoints.
Journal of Clinical Microbiology. Vol.49, N.11: 3931-3933.
- Wolter DJ, Kurpiel PM, Woodford N, Palepou MF, Goering RV, Hanson ND.
[Phenotypic and enzymatic comparative analysis of the novel KPC variant KPC-5 and its evolutionary variants, KPC-2 and KPC-4.](#)
Antimicrob Agents Chemother. 2009 Feb;53:557-62.
- Papp-Wallace KM, Bethel CR, Distler AM, Kasuboski C, Taracila M, Bonomo RA.
[Inhibitor resistance in the KPC-2 beta-lactamase, a preeminent property of this class A beta-lactamase.](#)
Antimicrob Agents Chemother. 2010 Feb;54:890-7.