Drigalski Lactose Agar

Selective and differential medium for isolating Enterobacteriaceae and certain non-fermenters from clinical specimens.

Instructions For Use

ENGLISH

DESCRIPTION

Drigalski Lactose Agar is a medium used for the isolation and differentiation of Enterobacteriaceae and non-fermenting gram-negative bacilli (NFGNB), including *Acinetobacter baumannii*, from clinical specimens.

TYPICAL FORMULA*

	(g/litre)
Peptone	15.0
Meat Extract	3.0
Yeast Extract	3.0
Sodium Deoxycholate	1.0
Sodium Thiosulfate	1.0
Lactose	15.0
Crystal Violet	0.005
Bromothymol Blue	0.08
Agar	15.0
Final pH 7.3 ± 0.2 at 25°C	

*Formula may be adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone and meat extract provide of amino acids, nitrogen, carbon, minerals and vitamins for organisms growth. Yeast extract is a source of vitamins of group B. Sodium thiosulfate, desoxycholate and crystal violet suppress the growth of gram-positive bacteria. Lactose is the fermentable carbohydrate. Bromothymol blue is the pH indicator. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 53.1 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the medium by spreading the sample over the agar surface to obtain well isolated colonies. Incubate plates inverted, aerobically at 35 ± 2°C for 18-24 hours.

INTERPRETING RESULTS

Lactose fermenters, such as *E. coli*, *Klebsiella* and *Enterobacter*, produce acid forming yellow colonies surrounded by yellow medium.

Non-fermenters, such as *Proteus*, *Providencia*, *Hafnia*, *Salmonella*, *Serratia*, *Alcaligenes*, *Pseudomonas* and *Acinetobacter*, cultivate with blue-grey to greenish colonies surrounded by blue to blue-green medium.

Additional biochemical or immunological tests may be required for confirmation.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts of plated medium to prevent condensation.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in bottles: 2 years.

Ready-to-use plates: 6 months.

QUALITY CONTROL

Appearance of Dehydrated Medium: Free-flowing, homogeneous, beige.

Appearance of Prepared Medium: Slightly opalescent, greenish-blue.

Expected Cultural Response:

Control strain		Inoculum	Incubation	Criteria	Specification
<i>Escherichia coli</i>	ATCC® 25922	50-100 CFU	18-24 h / 35 ± 2°C	P _R ≥ 0.5	Yellow colonies on yellow medium
<i>Salmonella</i> Typhimurium	ATCC® 14028				Blue-grey colonies on blue medium
<i>Acinetobacter baumannii</i>	ATCC® 19606				Greenish colonies on blue medium
<i>Enterococcus faecalis</i>	ATCC® 29212	10 ⁴ -10 ⁶ CFU		—	Partial to complete inhibition
<i>Staphylococcus aureus</i>	ATCC® 25923			—	Inhibition

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Drigalski Lactose Agar	90 mm Plate	20 plates	10018
Drigalski Lactose Agar	Bottle	6 x 100 ml bottles	402220
Drigalski Lactose Agar	Dehydrated medium	500 g of powder	610016
Drigalski Lactose Agar	Dehydrated medium	100 g of powder	620016

This document is available from the online Support Center:

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