TABLE OF CONTENTS

	English	Example pictures	Italiano	Français	Español	Deutsch	Čeština	Polski
PAGE	1	4	5	8	11	15	18	21

Instructions For Use
ENGLISH

Fastidious Anaerobe Agar w/ Horse Blood

Medium for isolation and susceptibility testing of anaerobes.

INTENDED PURPOSE

Fastidious Anaerobe Agar w/ Horse Blood (FAA-HB) is a plated medium used for the cultivation and antimicrobial susceptibility testing (AST) of anaerobic bacteria.

This medium is recommended by the European Committee on Antimicrobial Susceptibility Testing (EUCAST) for disk diffusion of selected rapidly growing anaerobic bacteria, i.e. *Bacteroides* spp., *Prevotella* spp., *Fusobacterium necrophorum*, *Clostridium perfringens* and *Cutibacterium acnes*.

DESCRIPTION

With increasing antimicrobial resistance, AST of anaerobic bacteria becomes more relevant. FAA-HB permits excellent growth of most anaerobic bacteria and is the medium chosen by EUCAST for the development of a disk diffusion method for anaerobes.

TYPICAL FORMULA* (Per Litre of Purified Water)

Peptone Mix	23.0 g
Sodium Chloride	5.0 g
Soluble Starch	1.0 g
Sodium Bicarbonate	0.4 g
Glucose	1.0 g
Sodium Pyruvate	1.0 g
Cysteine HCl monohydrate	0.5 g
Haemin	0.01 g
Vitamin K	0.001g
L-Arginine	1.0 g
Soluble Pyrophosphate	0.25 g
Sodium Succinate	0.5 g
Agar	12.0 g
Defibrinated Horse Blood	50.0 ml

Final pH 7.2 ± 0.2 at 25°C

*Adjusted and/or supplemented as required to meet performance specifications.

METHOD PRINCIPLE

Peptone mix provides amino acids, nitrogen, minerals, vitamins, and nutrient factors for growth of fastidious organisms. Sodium chloride maintains the osmotic balance of the medium. Starch and sodium bicarbonate act as detoxification agents. The 0.1% of glucose stimulates growth of some anaerobes while preventing the production of high levels of acids and alcohols which would inhibit colonial development. Pyruvate helps neutralize hydrogen peroxide. Vitamin K and sodium succinate provide essential growth factors for some anaerobes. Specific growth promoting agents are cysteine for *Fusobacterium necrophorum*, *Propionibacterium acne* and *Bacteroides fragilis*, arginine for *Eubacterium* spp., soluble pyrophosphate for *Porphyromonas gingivalis* and *Porphyromonas asaccharolyticus*, and haemin for *Porphyromonas melaninogenicus*. Agar is the solidifying agent.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Test tubes, inoculating loops, swabs, incubator, quality control organisms.

SPECIMENS

FAA-HB is used for AST of pure cultures that have been isolated from clinical specimens.

TEST PROCEDURE

Following EUCAST guidelines, FAA-HB plates are inoculated with a bacterial suspension in 0.85% saline from overnight cultures, adhering to the 15-15-15-minute rule (suspension used within 15 minutes, disks applied within 15 minutes of inoculation, and incubation within 15 minutes of applying the disks). The antibiotic disks are applied to the agar surface. Plates are incubated in an inverted position at $36 \pm 1^\circ\text{C}$ for 16-20 hours in an anaerobic atmosphere (80% N₂, 10% H₂, 10% CO₂).

For more details, please refer to the current EUCAST disk diffusion methodology for anaerobic bacteria.

Note:

1. FAA-HB is not intended for the direct culture of clinical specimens. The medium must be inoculated with pure calibrated strains obtained by culture on a Petri plate.
2. Ensure the FAA-HB plates are dried prior to inoculation to avoid excess moisture, which may lead to incorrect results.
3. Do not pre-reduce FAA-HB plates in an anaerobic environment before use.
4. Do not prolong incubation (beyond 20 h) as it will affect zone sizes and invalidate interpretative criteria.

INTERPRETING RESULTS

At the end of the incubation period, the inhibition zone diameters are measured and interpreted according to the current EUCAST breakpoints and QC criteria for disk diffusion of anaerobes.

Refer to the current EUCAST disk diffusion reading guide for anaerobic bacteria.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

SHELF LIFE

2 months.

QUALITY CONTROL

Appearance of Medium: Opaque, ruby red.

The agar depth must be 4.0 ± 0.5 mm.

The performance of the medium can be tested using the following strains:

- *Bacteroides fragilis* ATCC 25285 (NCTC 9343)
- *Clostridium perfringens* ATCC 13124 (NCTC 8237)

Expected Cultural Response:

EUCAST Tentative QC Criteria

		<i>B. fragilis</i> ATCC 25285	<i>C. perfringens</i> ATCC 13124
Antimicrobial Agent	Disk Content	Inhibition Zone Diameter (mm)	
Benzylpenicillin	1 unit	—	23-29
Piperacillin-tazobactam	36 µg	29-35	30-36
Meropenem	10 µg	32-39	34-40
Clindamycin	2 µg	23-29	20-26
Metronidazole	5 µg	29-36	20-26
Vancomycin	5 µg	—	14-20

To monitor the anaerobic atmosphere, a Metronidazole 5 µg disk can be tested with *Clostridium perfringens* DSM 25589. A zone diameter of <25 mm indicates an insufficient anaerobicity, which may affect growth and susceptibility test results of anaerobic bacteria.

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of FAA-HB was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Invalid results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology. The operator should understand the principles of the procedures, including its performance limitations, in advance of operation to avoid potential mistakes.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium.

According to the EUCAST disk diffusion method for anaerobic bacteria, the FAA-HB plates can be used only with selected species of anaerobes.

Various factors have been identified as influencing disk diffusion susceptibility tests. These include the medium, agar depth, disk potency, inoculum concentration, age of inoculum, and pH. It is strongly recommended to follow the EUCAST guidelines to obtain reliable results.

Improper storage of antimicrobial disks may cause a loss of potency and a falsely resistant result. Excessive shrinkage of the medium due to improper storage may lead to falsely susceptible results.

In vitro susceptibility of an organism to a specific antimicrobial agent does not necessarily mean that the agent will be effective *in vivo*. Consult appropriate references for guidance in the interpretation of results.

FAA-HB is intended as an aid in the diagnosis of infectious diseases, requiring further tests to complete the diagnostic results. All identification tests should ideally be performed from non-selective agar.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

See ordering info below. There may be additional product ref. numbers as well. For an updated listing of available products, visit liofilchem.com

Product	Format	Packaging	Ref.
Fastidious Anaerobe Agar w/ Horse Blood (FAA-HB)	Plate 90 mm	20 (2 x 10) plates	10062

Revision History

Revision	Release Date	Change Summary
3	2023-02-06	Product name changed from "Fastidious Anaerobe Agar (FAA)" to "Fastidious Anaerobe Agar w/ Horse Blood (FAA-HB)" Updated references
2	2022-12-16	Updated layout and content in compliance with IVDR 2017/746 Added example photographs

This IFU document and the SDS are available from the online Support Center:

liofilchem.com/ifu-sds

B. fragilis ATCC 25285

C. perfringens ATCC 13124

Fastidious Anaerobe Agar w/ Horse Blood

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1. Bavelaar H et al. Development of a EUCAST disk diffusion method for the susceptibility testing of rapidly growing anaerobic bacteria using Fastidious Anaerobe Agar (FAA): a development study using Bacteroides species. Clin Microbiol Infect. 2021 Nov;27(11):1695.e1-1695.e6. DOI: [10.1016/j.cmi.2021.03.028](https://doi.org/10.1016/j.cmi.2021.03.028)
2. EUCAST disk diffusion of anaerobic bacteria. Version 2.0, January 2023.
3. EUCAST disk diffusion reading guide for selected rapidly growing anaerobic bacteria on Fastidious Anaerobe Agar with 5% horse blood (FAA-HB). Version 2.0, January 2023.
4. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 13.0, 2023.
5. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 13.0, 2023.

Above documents are available from the EUCAST website at <http://www.eucast.org>.

Table of Symbols / Tabella dei Simboli / Tableau des Symboles / Tabla de Símbolos / Tabelle der Symbole / Tabulka Symbolů / Tabela Symboli

	Batch code / Codice del lotto / Numéro de lot / Código de Lote / Batch-Code / Kód šarže / Kod partii
	Catalogue number / Numero di catalogo / Numéro de catalogue / Numero di catalogo/número de catálogo / Katalognummer / Katalogové číslo / Numer katalogowy
	<i>In Vitro</i> Diagnostic Medical Device / Dispositivo Medico Diagnostico <i>in vitro</i> / Dispositif médical de diagnostic <i>in vitro</i> / Dispositivo Medico Diagnostico <i>in vitro</i> / Dispositivo médico de diagnóstico <i>in vitro</i> / <i>In-vitro-Diagnostikum</i> / Diagnostický lékařský pomůcka <i>in vitro</i> / Wyrób medyczny do diagnostyki <i>in vitro</i>
	Manufacturer / Fabbicante / Fabricant / Fabricante / Hersteller / Výrobce / Producent
	Use by / Utilizzare entro / Utiliser par / Usar por / Verwenden von / Použití do / Używany przez
	Fragile, handle with care / Fragile, maneggiare con cura / Fragile, manipuler avec soin / Frágil, manejar con cuidado / Zerbrechlich - Vorsichtig behandeln / Křehké, zacházejte opatrně / Delikatne, obchodź się ostrożnie
	Temperature limitation / Limiti di temperatura / Limite de température / Límite de temperatura / Temperaturbegrenzung / Teplotní limity / Ograniczenie temperatury
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contient des éléments suffisants pour <n> tests / Contiene suficiente para <n> pruebas / Enthält ausreichend für <n> Tests / Obsahuje dostatek pro <n> testů / Zawiera ilość wystarczającą na <n> testów
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter la notice d'utilisation / Consulte instrucciones para su uso / Gebrauchsanweisung beachten / Přečtěte si návod k použití / Zapoznaj się z instrukcją użytkowania
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilice / Nicht wiederverwenden / Nepoužívejte opakovaně / Nie używać ponownie

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