## Chromatic™ Candida

Chromogenic selective medium for the isolation and differentiation of *Candida* spp directly from clinical and nonclinical specimens.

### DESCRIPTION

Chromatic™ Candida is a chromogenic selective medium used for the isolation and differentiation of *Candida* species directly from clinical and nonclinical specimens permitting to distinguish among *C.albicans*, *C.tropicalis*, *C.krusei*, *C.glabrata*, *C.dubliniensis* and *C.parapsilosis*.

Although *Candida albicans* remains the most common cause of human Candidiasis, the frequency of infection attributed to other members of the genus is also increasing. Effective treatment requires both early diagnosis and prompt initiation of therapy against fungal infection.

### TYPICAL FORMULA (g/l)

Peptone	10.0
Chloramphenicol	0.5
Chromogenic Mix	25.2
Agar	15.0

Final pH 6.1 ± 0.2 at 25°C

### METHOD PRINCIPLE

Peptone provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Chloramphenicol is the selective agent inhibiting most of the bacteria. Chromogenic mix allows to identify the *Candida* genus on the basis of the color and morphology of the colonies. Agar is the solidifying agent.

### PREPARATION

Dehydrated medium Suspend 50.7 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE.

Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

### TEST PROCEDURE

Inoculate the medium by direct streaking, spread plating or membrane filtration method. Incubate aerobically at 30-37°C for 24-48 hours.

### INTERPRETING RESULTS

After incubation observe the color and the morphology of the colonies and interpret the results as indicated in the ID table.

#### ID Table.

Microorganism	Typical colony color
<i>Candida albicans</i>	Green
<i>Candida dubliniensis</i>	Yellow-green
<i>Candida glabrata</i>	Beige
<i>Candida krusei</i>	Pink, pale edges
<i>Candida parapsilosis</i>	Pale pink-white
<i>Candida tropicalis</i>	Blue

See pictures in Appendix I.

### APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.

Prepared medium: slightly opalescent, very light beige.

**STORAGE**

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

**SHELF LIFE**

Dehydrated medium: 4 years.  
Medium in bottles: 1 year.  
Ready-to-use plates: 6 months.

**QUALITY CONTROL**

Plates are inoculated with the microbial strains indicated in the QC table.  
Inoculum for productivity: 50-100 CFU  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>6</sup> CFU.  
Incubation conditions: aerobically at 35 ± 2°C for 24-48 hours.

**QC Table.**

Microorganism		Growth	Specification
<i>Candida albicans</i>	ATCC® 10231	Good	Green colonies
<i>Candida krusei</i>	ATCC® 14243	Good	Pink colonies
<i>Candida parapsilosis</i>	ATCC® 22019	Good	Pale pink-white colonies
<i>Candida tropicalis</i>	ATCC® 750	Good	Blue colonies
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---
<i>Staphylococcus aureus</i>	ATCC® 25923	Inhibited	---

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**BIBLIOGRAPHY**

- Odds, F.C. And Bernaerts. 1994. CHROMagar Candida, a new differential medium for presumptive identification of clinically important *Candida* species. J. Clin. Microbiol. 32: 1923-1929.
- Wingard, JR. Importance of *Candida* species other than *C. albicans* as pathogens in oncology patients. Clin Infect Dis. 1995; 20: 115-25.
- Pfaller, Huston and Coffman. 1996. J. Clin. Microbiol. 32: 1923-1929.
- Maertens JA. History of the development of azole derivatives. J Clin Microbiol Infect. 2004; 10: 1-10.

**PRESENTATION**

		Contents	Ref.
Chromatic™ Candida	90 mm ready-to-use plates	20 plates	11612
Chromatic™ Candida	60 mm ready-to-use plates	20 plates	163692
Chromatic™ Candida	Bottles	6 x 100 ml bottles	481110
Chromatic™ Candida	Dehydrated medium	500 g of powder	610613
Chromatic™ Candida	Dehydrated medium	100 g of powder	620613

**TABLE OF SYMBOLS**

<b>LOT</b> Batch code	<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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**BIBLIOGRAFIA**

- Odds, F.C. And Bernaerts. 1994. CHROMagar Candida, a new differential medium for presumptive identification of clinically important Candida species. J. Clin. Microbiol. 32: 1923-1929.
- Wingard, JR. Importance of Candida species other than C. albicans as pathogens in oncology patients. Clin Infect Dis. 1995; 20: 115-25.
- Pfaller, Huston and Coffman. 1996. J. Clin. Microbiol. 32: 1923-1929.
- Maertens JA. History of the development of azole derivatives. J Clin Microbiol Infect. 2004; 10: 1-10.

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- Odds, F.C. And Bernaerts. 1994. CHROMagar Candida, a new differential medium for presumptive identification of clinically important Candida species. J. Clin. Microbiol. 32: 1923-1929.
- Wingard, JR. Importance of Candida species other than C. albicans as pathogens in oncology patients. Clin Infect Dis. 1995; 20: 115–25.
- Pfaller, Huston and Coffman. 1996. J. Clin. Microbiol. 32: 1923-1929.
- Maertens JA. History of the development of azole derivatives. J Clin Microbiol Infect. 2004; 10: 1–10.

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		Contenido	Ref.
Chromatic™ Candida	Placas listas para uso de 90 mm	20 placas	11612
Chromatic™ Candida	Placas listas para uso de 60 mm	20 placas	163692
Chromatic™ Candida	Botellas	6 x 100 ml botellas	481110
Chromatic™ Candida	Medio deshidratado	500 g de polvo deshidratado	610613
Chromatic™ Candida	Medio deshidratado	100 g de polvo deshidratado	620613

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